

## Product datasheet for **RG210055**

### EXOSC10 (NM\_001001998) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	EXOSC10 (NM_001001998) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	EXOSC10
Synonyms:	p2; p3; p4; PM-Scl; PM/Scl-100; PMSCL; PMSCL2; RRP6; Rrp6p
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide  
Sequence:**

>RG210055 representing NM\_001001998  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCGCCACCCAGTACCCGGGAGCCAGGGTCTGTGCGCGACCAGCGCAACCAATCCGACGGAGAGA  
 TGGTGCTGCCAGGCTTCCCGACGCCGACAGCTTTGTGAAGTTTCTCTGGGTCCGTGGTGGCAGTCAC  
 CAAGGCATCTGGGGCCCTACCACAGTTTGGCGATGAGTATGATTTTTACCGAAGTTTTCTGGCTTCCAA  
 GCATTTTGCAGAACACAGGGAGACAGGTTGCTTCAGTGCATGAGCAGAGTAATGCAGTACCATGGGTGTC  
 GCAGCAACATTAAGGATCGAAGTAAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGT  
 TGATGTAATCTGGAGAGAGTGGGTATTTACTGGATGAAGCCTCAGGTGTAAACAAGAATCAACAGCCT  
 GTCCTCCCTGCCGGCTTGCAGGTCCCAAAACGGTAGTGTCCAGCTGGAACCGTAAGGCAGCAGAATATG  
 GCAAAAAAGCAAATCTGAAACTTCCGGCTGCTTCATGCAAAAAATATCATCCGACCTCAGCTCAAGTT  
 TCGAGAGAAGATTGACAATTCACACACACCATTTCTTCTAAAACTTTCATCAAACCCAATGCTCAGAAA  
 CCTCTCCCTCAAGCTCTCTAAGGAAAGCGGGAAACGCCACAGGATCGTCTGAGGACTTGGAGCTCC  
 CCCCTGCACTGGCTGATTTTCATCCATCAGCAGAGAACCCAGCAGGTTGAGCAAGACATGTTTGCACATCC  
 TTATCAATATGAACTAAATCACTTTACCCAGCAGATGCAGTGCCTTCAAAGCCACAACCCAGTTATAC  
 AGACCTATAGAAGAGACACCATGCCATTTTCATATCCTCCCTGGATGAACTCGTGGAACTCAACGAAAAGC  
 TCTTGAATTGTCAGGAATTTGCAGTTGACTTGGAGCACCCTTACAGGAGCTTCTGGGACTGACCTG  
 CCTGATGCAAATTTCTACTCGGACGGAAGACTTCATCATTGACACCCTCGAGCTTCAAGTGACATGTAC  
 ATTTCTCAATGAGAGCCTCACAGACCCAGCCATCGTTAAGGTCTTTCATGGTGTGATTGACATAGAAT  
 GGCTCAGAAAAGACTTTGGTTGTATGTAGTAAACATGTTTGAATCATCAGGCACGCCTTTTAA  
 CCTGGGCAGGCACTCACTCGATCATCTCCTGAAACTCTACTGCAACGTGGACTCAAACAAGCAATATCAG  
 CTGGCTGATTGGAGAATACGCCCTCTGCCCGAGGAGATGCTCAGCTACGCCCGGGATGACACCCATTACC  
 TGCTATATATCTATGACAAAATGAGGCTGGAGATGTGGGAGCGCGGCAACGGGCAGCCCGTGCAGCTGCA  
 GGTGGTGTGGCAACGGAGCAGGGACATCTGCCTCAAGAAATTCATCAAACCTATCTTACGGATGAGTCC  
 TACCTTGAATCTATAGGAAGCAGAAGAAGCACCTTAACACACAGCAGTTGACAGCCTTTCAGCTGCTGT  
 TTGCCTGGAGGGATAAAACAGCTCGCAGGGAAGATGAAAGTTACGGATATGACTGCCAAACCACATGAT  
 GCTGAAAATAGCTGAAGAACTGCCTAAGGAACCTCAGGGCATCATAGCTTGTGCAACCCAGTACCGCC  
 CTTGTGCGGCAGCAGATCAACGAAATGCACCTTTAATCCAGCAGGCCCGAGAGATGCCCTGCTCAAGT  
 CTGAAGTTGCAGCCGGAGTGAAGAAGAGCGGACCCTGCCAGTCTGAGAGATTGGAGAATGTTCTCTT  
 TGGACCTCACGACTGCTCCCATGCCCTCCGGATGGCTATCCAATCATCCAACCCAGTGGATCTGTGCCA  
 GTTCAGAAGCAGGCGAGCCTCTCCCTGATGAAAAAGAAGATAACTTGTGGGTACCACATGCCTGATTG  
 CCACAGCTGTATCACGTTATTTAATGAACCTAGTGTGAAGACAGTAAAAAGGGTCCATTGACAGTTGC  
 ACAGAAAAAAGCCAGAACATCATGGAGTCTTTGAAAATCCATTTAGGATGTTTCTGCCCTCACTGGGA  
 CACCGTGTCCCGTCTCTCAGGCAGCGAAGTTCGATCCATCAACCAAAATCTATGAAATCAGCAACCGTT  
 GGAAGCTGGCCAGGTACAAGTACAAAAGACTCTAAAGAAGCTGTCAAGAAGAAGGCAGCTGAGCAAAC  
 AGCTGCCCGGGAACAGGCAAAGGAGGCGTGCAAGCTGCAGCAGAACAGGCCATCTCCGTCGCAGCAG  
 GTCGTGCTAGAAAATGCTGCAAAGAAGAGAGAGCGAGCAACAAGCGACCCAAGGACCACAGAACAGAAAC  
 AAGAGAAGAAACGACTCAAAATTTCCAAGAAGCCAAAGGACCCAGAGCCACCAGAAAAAGAGTTTACGCC  
 TTACGACTACAGCCAGTCAAGTCAAGGCTTTTGTGAAAACAGCAAATCCAAAGTTTCTTCTCAGTTT  
 GATCCAAATAAACAGACCCCGTCTGGCAAGAAATGCATTGCGACCAAAAAAATTAACAGTCGGTGGGAA  
 ACAAAGCATGCTCTTCCAAGTGGAAAGTCAAGCAGAGGCTTCAGGTACAACCTGGCCACAGAGA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG210055 representing NM\_001001998  
 Red=Cloning site Green=Tags(s)


MAPPSTREPRVLSATSATKSDGEMVLPGFADDSFVKFALGSVVAVTKASGGLPQFGDEYDFYRSFPGFQ  
 AFCETQGDRLQCMRSVMQYHGCRSNIKDRSKVTELEDFKDLLVDANDVILERVGILLDEASGVNKNQQP  
 VLPAGLQVPKTVVSSWNRKAAEYGKKAKSETFRLLHAKNIIRPQLKFKREKIDNSNTPFLPKIFIKPNAQK  
 PLPQALSKERRERPQDRPEDLDVPPALADFIHQRTQQVEQDMFAHPYQYELNHFTPADAVLQKPQPLY  
 RPIEETPCHFISLDELVELNEKLLNCQEFADVLEHHSYRSFGLTCLMQISTRTEDFIIDTLELRSDMY  
 ILNESLTDPAIVKVFHGADSDIEWLQKDFGLYVVMFDMTHQAARLLNLGRHSLDHLKLYCNVDSNKQYQ  
 LADWRIRPLPEEMLSYARDDTHYLLIYDKMRLEMWERNQGPVQLQVVWQSRDICKKFIKPIFTDES  
 YLELYRKQKKHLNTQQLTAFQLLFAWRDKTARREDESYGYVLPNHMLKIAEELPKPEQGIACCNPVPP  
 LVRQQINEMHLLIQQAREMPLLKSEVAAGVKKSGPLPSAERLENVLFPHDCSHAPPDGYPIIPTSGSVP  
 VQKQASLFPDEKEDNLLGTTCLIAATAVITLNFPSAEDSKKGPLTVAQKKAQNIMESFENPFMFLPSLG  
 HRAPVSQAAKFDPSTKIYEISNRWKLAVQVQVQKDSKEAVKKKAAEQTAAREQAKEACKAAAEQAI SVRQQ  
 VVLENAAKKRERATSDPRTEQKQEKRLKISKPKDPEPPEKEFTPYDYSQSDFKAFAGNSKSKVSSQF  
 DPNKQTPSGKKCIAAKKIKQSVGNKSMFPTGKSDRGFRYNWPQR

TRTRPLE - GFP Tag - V

**Restriction Sites:**

SgfI-MluI


**Cloning Scheme:**



<b>ACCN:</b>	NM_001001998
<b>ORF Size:</b>	2655 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001001998.1</a> , <a href="#">NP_001001998.1</a>
<b>RefSeq Size:</b>	2834 bp
<b>RefSeq ORF:</b>	2658 bp
<b>Locus ID:</b>	5394
<b>UniProt ID:</b>	<a href="#">Q01780</a>
<b>Cytogenetics:</b>	1p36.22
<b>Protein Pathways:</b>	RNA degradation

**Gene Summary:**

Putative catalytic component of the RNA exosome complex which has 3'->5' exoribonuclease activity and participates in a multitude of cellular RNA processing and degradation events. In the nucleus, the RNA exosome complex is involved in proper maturation of stable RNA species such as rRNA, snRNA and snoRNA, in the elimination of RNA processing by-products and non-coding 'pervasive' transcripts, such as antisense RNA species and promoter-upstream transcripts (PROMPTs), and of mRNAs with processing defects, thereby limiting or excluding their export to the cytoplasm. The RNA exosome may be involved in Ig class switch recombination (CSR) and/or Ig variable region somatic hypermutation (SHM) by targeting AICDA deamination activity to transcribed dsDNA substrates. In the cytoplasm, the RNA exosome complex is involved in general mRNA turnover and specifically degrades inherently unstable mRNAs containing AU-rich elements (AREs) within their 3' untranslated regions, and in RNA surveillance pathways, preventing translation of aberrant mRNAs. It seems to be involved in degradation of histone mRNA. EXOSC10 has 3'-5' exonuclease activity (By similarity). EXOSC10 is required for nucleolar localization of C1D and probably mediates the association of MTREX, C1D and MPP6 with the RNA exosome involved in the maturation of 5.8S rRNA.[UniProtKB/Swiss-Prot Function]

**Product images:**


Circular map for RG210055