

Product datasheet for **RG210057**

ZNF282 (NM_003575) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF282 (NM_003575) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ZNF282
Synonyms:	HUB1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG210057 representing NM_003575
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGCAGTTTGTGTCAACACGGCCGAGCCTCAGCAGCTGGGCATCCAGGGCCTGGGGCTGGACAGCGGA
 GCTGGAGCTGGGCCAGGCTCTGCCCCGAGGAGGTCTGCCACCAGGAGCCGGCCTGCGCGGGAAAT
 GGCCGAGGGAATGCCGCCATGCAGGCTCAAGAATGGGACATGGACGCCCGCGGCCAATGCCTTTTCAG
 TTCCACCCCTTTCCAGATAGGGCACCTGTCTTCCCCGACCGCATGATGCGAGAGCCCCAGTTGCCACAG
 CAGAGATCTCACTCTGGACTGTGGTGGCTGCCATTAGGCTGTGGAGAGGAAGGTGGATGCCAGGCCAG
 CCAGCTGCTGAACCTGGAGGGGCGACGGGGACAGCCGAGAAGAAGCTGGCCGACTGTGAAAGACGGCC
 GTGGAATTTGGGAACCACATGGAGAGCAAGTGGCCGTGCTGGGACCCTGCTGCAGGAGTACGGGCTGC
 TGCAGAGCGGCTGGAGAACTTGGAGAACTTGTGCGCAACAGGAACCTTCTGGGTCTGCGGCTGCCCC
 GGGCAGCAAGGGGAGGCCCCCAAGGTTCCAGTGACTTTTGTGACATTGCTGTGTACTTCTCCGAAGAC
 GAGTGGAGAAGTGGACGAATGGCAGAAGGAGCTTTATAACAACCTTGTTAAGGAGAATACAAAACCC
 TCATGTCCCTGGACGCGGAGGGCTCAGTCCCCAAGCCAGATGCTCCAGTCCAGGCTGAGCCCAGGGAAGA
 ACCTTGTGTGTGGGAGCAGCGCCACCCGAAGAGAGAGAAATCCCAATGGATCCCGAAGCAGGAGCAGAG
 CCCCTGGTGCCTGCGCAGGATGCGTCTCCAGGTGAAGCGTGAGGACACCCTGTGTGTCGGGGTCAAG
 GGGGCTGGAGGAAAGAGCCATCCCTACGGAATCCATTACCGACTCCCCAATTTCTGCCAGGACCTCTT
 GTCCCGGATTAACAGGAGGAGCATCAGTGCCTGTGGGATCAGCAGGATTTGGCAGACAGAGATATCC
 ACGGATCCCAATTCAGAGTCTCTCATCTCAGCACATGACATTTTGTGATGGATCAAGCAGGAGGAGCAGC
 CATACCATGGGACACGCGACTCAATGGACGGAGAGCTTGGATTAGACTCTGGCCCTAGTGACAGCCT
 GCTGATGGTGAAGAACCACCCCGGCCCGCCACAGCCCCAGCCCCAGCCCCAGCCACCCGAGCCGCGAG
 CTGCAGTCGACAGCCCCAGCCCCAGAGCCTGCCCCCATCGCGGTGGCCGAGAACCAGGGCGGCCCGGA
 GCCGAGGGCTGCTGGACGACGGTTTCCAGGTGCTGCCCGGGAGCGTGGCTCCGGCGAGGCGCGCGGG
 TGGGGACCGCAGCACCGGGGCGCGGGGGCGATGGGGGCGGTGGGGGCGCGCGCGGAGGCGGGGACG
 GGGGACGGCGCGGCTGTGGCAGCTGCTGCCCTGGCGGGCTGCGCGGAGCCTCCTCTGCACGGCGCCC
 GCAGCAAGCCCTACTCGTGCCCGAGTGCAGCAAGAGCTTCGGCGTGCGCAAGAGCCTCATCATCCACCA
 CCGCAGCCACACCAAGGAGCGGCCCTACGAGTGCCTGAGTGCAGAGAAGAGCTTCAACTGCCACTCGGGC
 CTCATCCGCCACCAGATGACGCACCGCGCGGAGCGGCCCTACAAGTGTCTGGAGTGCAGAGAGCTACA
 GCCGTAAGGAGCACCTGCAGAACCACAGCGGCTGCACACGGGCGAGCGGCTTTCCAATGTGACTGTG
 CGGCAAGAGCTTCATCCGCAAGCAGAACCTGCTCAAGCACCAGCGCATCCACACGGGCGAGCGCCCTAC
 ACGTGGGCGAGTGGGCAAGAGCTTCCGCTACAAGGAGTTCGCTCAAGGACCACCTGCGCGTGCACAGCG
 GCGGCCGGGCCCCGCGGCCACGGCAGCTCCCGCCGCTCCTGAGCGAGAC

ACGCGTACGCGGGCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG210057 representing NM_003575
Red=Cloning site Green=Tags(s)


MQFVSTRPQPQLGIQGLDGSWSWAQALPPEEVCHQEPALRGEMAEGMPPMQAQEWMDARRPMPFQ
FPPFPDRAPVFPDRMRPQLPTAEISLWTVVAAIQAVKVDASQLLNLEGRGTAEKKLADCEKTA
VEFGNHMESKWAVLGTLLQEYGLLQRRLENLENLLRNRNFWLRLPPGSKGEAPKVPVTFVDIAVYFSED
EWKNLDEWQKELYNL VKENYKTLMSLDAEGSVKPDAPVQAEPREPCVWEQRHPEEREIPMDPEAGAE
PLVPAQDASSQVKREDTLCVRGQRGLEERAIPTESITDSPISAQDLLSRIKQEEHQCVWDQDLADRDIP
TDPNSELISAHILSWIKQEEQYPWPGRDSMDGELGLDSGPSDSLMMVKNPPAPPQPQPQPQPQ
LQSQPQPQSLPPIAVAENPGPPSRGLLDDGFQVLPGERGSGEAPPGGDRSTGGGGGDDGGGGGAEAGT
GAGGGCGSCCPGGLRRSLLLHGARSKPYSCPECGKSFGVRKSLIIHRSHTKERPYECAECEKSFNCHSG
LIRHQMTHERGERPYKCSECKTYSRKEHLQNHQRLHTGERPFQCALCGKSFIRKQNLKQHRIHTGERPY
TCGECGKSFYKESLKDHLRVHSGGPGGAPRQLPPPPERD

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_003575

ORF Size: 2013 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_003575.4](#)

RefSeq Size: 3736 bp

RefSeq ORF: 2016 bp

Locus ID: 8427


UniProt ID: [Q9UDV7](#)

Cytogenetics: 7q36.1

Protein Families: Transcription Factors

Gene Summary: Binds to the U5 repressive element (U5RE) of the human T cell leukemia virus type I long terminal repeat. It recognizes the 5'-TCCACCCC-3' sequence as a core motif and exerts a strong repressive effect on HTLV-I LTR-mediated expression.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG210057