

Product datasheet for **RG210149**

OR3A1 (NM_002550) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OR3A1 (NM_002550) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	OR3A1
Synonyms:	OLFRA03; OR17-40; OR17-82; OR40
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG210149 representing NM_002550 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAGCCAGAATCTGGGGCCAATGGAACAGTCATTGCTGAGTTCATCCTGCTGGGCTTGCTGGAGGCGC
CAGGGCTGCAGCCAGTTGTCTTTGTGCTCTTCTCTTTGCCTACCTGGTCACGGTCAGGGGCAACCTCAG
CATCCTGGCAGCTGTCTTGGTGGAGCCAACTCCACCCCCATGACTTCTTCTGGGGAACCTATCA
GTGCTGGATGTTGGGTGCATCAGCGTCACTGTTCCATCAATGTTGAGTCGTCTCCTGTCCGCAAGCGTG
CAGTTCCTGTGGGGCCTGCCTTACCCAGCTCTTCTTCTCCATCTGTTGTTGGAGTGGACTGCTTCTCT
GCTGACCGCCATGGCCTATGACCGATTCTGGCCATCTGCCGGCCCTCACCTACAGCACCCGCATGAGT
CAGACAGTCCAGAGGATGTTGGTGGCTGCGTCTGGGCTTGCTTTACCAACGCACTGACCCACACTG
TGGCCATGTCCACGCTCAACTTCTGTGGCCCAATGTGATCAATCACTTCTACTGTGACCTCCCACAGCT
CTTCCAGCTCTCCTGCTCCAGCACCCAATCAATGAGCTGCTGCTTTTGTGCTGGGTTTTATAATGGCA
GGTACCCCATGGCTCTCATTGTCTCCTATATCCAGTGGCAGCTGCAGTCTGCGAATTCGCTCTG
TAGAGGGCAGGAAGAAAGCCTTCTCCACATGTGGCTCCACCTCACTGTGGTTGCCATATTCTATGGTTC
AGGTATCTTTAACTATATGCGACTGGGTTCAACCAAGCTTTCAGACAAGGATAAAGCTGTTGGAATTTTC
AACACTGTCATCAATCCCATGCTGAACCAATCATCTACAGCTTCAGAAACCCTGATGTGCAGAGTGCCA
TCTGGAGGATGCTCACAGGGAGGCGGTCACTGGCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG210149 representing NM_002550
 Red=Cloning site Green=Tags(s)

MQPESGANGTVIAEFILLGLLEAPGLQPVVFLFLFAYLVTVRGNLSILA AVLVEPKLHTPMYFFLGNLS
 VLDVGCISVTPSMLSRLLSRKRAVPCGACL TQLFFFHLFVGVDCFLLTAMAYDRFLAICRPLTYSTRMS
 QTVQRMLVAASWACAFTNALTHTVAMSTLNFCGPNVINHFYCDLPQLFQLSCSSTQLNELLLFAVGFIMA
 GTPMALIVISYIHVAAAVLRIRSVVEGRKAFSTCGSHLTVVAIFYGSIFNYMRLGSTKLSDKDKAVGIF
 NTVINPMLNPIIYSFRNPDVQSAIWRMLTGRRSLA

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_002550

ORF Size: 945 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_002550.2 , NP_002541.2
RefSeq Size:	948 bp
RefSeq ORF:	948 bp
Locus ID:	4994
UniProt ID:	P47881
Cytogenetics:	17p13.3
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Olfactory transduction
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG210149