

Product datasheet for **RG210351**

CCDC157 (NM_001017437) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC157 (NM_001017437) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CCDC157
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG210351 representing NM_001017437
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGTGGCCCTGTGAGCACTATGACCATGTTCCGGGTGACCCCGAGTTCACGCAGCTGTCCCATGCCG
 TGCTGCTGGAGCTGGTCATCGACAGGCTCCTGCTGCTGCTTCAAAGCTGTATGAGCTACTTGGAGAACCT
 TGGCTCAGAGCAGATGATGCCCCCTGCACAGGCTGCGGGGCCCTGCATGTCCGTGGGGCTCACGGTGCGG
 CGCTTCTGGGACAGCCTGCTGAGGCTGGGCACGCTCCACCAGCAGCCACTCCCCAGAAAGGGGCAAACC
 AAAGGGAGACTCCACCTCCAAGCCACCACCAAGGGCGAGCCAGCCAGGAGCCCTGAATATCTGACTAC
 CAAGTTAATCAAGCCCTCCTCCCAGTGTAGGCTTGCCCGAGCCTGCCAAGAGCCAGAGAGCATCCCT
 GTCAGAGCCTCCCTGCAGTTCACAGCCACGACCTTCAAGAACCAGGAGTGTCCACTCCAGACCATTG
 AGACGGCCCTGGTGCCCTGTGACGCATGCGCCAGCGTCCAGGGAAGCCTGCAGAAGTGGCAAGGTGGT
 CATCAGCCTGTGTAGAGCCAGAACCTGCCCTCGTCCTTAGGCCAGTTCAGCAACTGGTGCAGGACAGC
 ATGGGGCTCAGGCCACTGCCGGCTGCCACCGTGGGCCGCTGGGCAGCAGAGCAGAGGAAAGACCTGACGC
 GCCTCAGTAAGCATGTGGAGGCCCTCAGGGCCAGCTGGAGGAGGCTGAAGGGCAGAAGGATGGCCTGAG
 GAAGCAGGCGGGCAAGCTGGAGCAGGCGCTGAAACAGGAGCAGGGGGCACGGCGGCACAGGCGGAGGAG
 GATGAGCAGTGCCTGTCTGAGTGGGAGCAGGACAAACAGCAGCTGCTCACAGAAACAAGTACCTAAAGA
 CAAAGATGGCCACCCTGGAGAGAGAAGTAAACAGCAGCGGGAGTCCACACAGGCTGTGGAGGCAAAGGC
 CCAGCAGCTGCAGGAGGAAGGTGAGCGCAGGGCGCAGCGGAGAGGCAGGTGCAGCAGCTGGAGGAGCAG
 GTGCAGCAGTTGGAGGCGCAGGTGCAGCTGTTGGTGGTGGCTGGAGGGCGCTGGCCAGCAGGTCTGCT
 GGGCCAGCAGGAGCTGGATAAGGAGAAGGCCCTGTCGACAGCATGGTCCGCCACCAGGAGTCTCTGCA
 GGCCAAGCAGCGAGCCCTGCTAAAGCAGCTGGACAGCCTGGACCAGGAACGTGAGGAGCTGCGGGGCAGC
 CTGGACGAGGCTGAGGCCAGCGGGCCCGCTGGAGGAGCAGCTGCAGAGCGAGCGGGAGCAGGGCAAT
 GCCAGCTCAGGGCCAGCAGGAGCTGCTGCAGAGCCTGCAGAGGAGAAGCAAGGCTGGAGCAGGCGAC
 TACGGACTGCGGCTAACCATCCTGGAGCTAGAACGGGAGCTGGAGGAGCTGAAGGAGCGGGAGCGGCTG
 CTGGTGGCCTTCCCAGACCTGCACAGGCCACCGAGCCAGATCCATGGAGGCAGATCCAGCAGCGTGG
 AATCCCAGATAACATGTCCACAGACTCTGGCAACGTACAGATCACATGGAGAGCAAGTGCAGGCCAA
 CGACATCCGCATCCGGTCTACAGGAGGAGAACGGGCGGCTCCAATCAATGCTGTCCAAAATCCGGGAA
 GTGGCCAGCAGGTTGGCCTCAAGCTGATCCCGCAGGACCGGCTCTGGTCCCTTCCAGCAAGGGAACCC
 AGGGAGCAACACCACAGTCCAGGCCAAGAGCACATCCCCAGGGCCTCTGGGCAGACAGCACCTTCTAG
 CAGCAGGACGGGTAGGACCCCTGCTGGGCCAGCCCTGCACATCCCCACCTCGGCAGCCCTGCACATCCCCA
 CCTCGGCAGCCCTGCACATCCCCACCTCGGCAGCCCTGCACATCCCCATCTCGGCAGCCCTGCAGCCAGC
 CCAGCAAGTCTTGTGGAGGTTGTGACCCACTTGACACCTGCACCCAGAACCCTCAAGGTCTTGGT
 CAGGCTGAGAAAGAGACTGTCACCTGGCCGGGGACAGGCCAGCTCTGCACACCAGCCCCAGGAGCGGCC
 ATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG210351 representing NM_001017437
 Red=Cloning site Green=Tags(s)


MVALLEHYDHPGDPEFTQLSHAVLLELVIDRLLLLLQSCMSYLENLGSEQMPPPAQAAGPCMSVGLTVR
 RFWDSLRLGLHQQPLPQKGANQRETPTSKPTTKGEPARSPEYLTTKL IKPSSPVLGLPQTCQEPESIP
 VRASLQFPATTFKNTRS VHSQT IETALVPCDACASVQGS LQKVGKVVISLCQS QNLPSSLGQFQQLVQDS
 MGLRPLPAATVGRWAAEQRKDLTRL SKHVEALRAQLEEAEGQKDGLRKQAGKLEQALKQEQGARRRQAE
 DEQCLSEWEHDKQQLLTETSDLKTKMATLERELKQRESTQAVEAKAQQLQEEGERRAAAERQVQLEEQ
 VQQLEAQVQLLVGRLEGAGQVCWASTELDKKARVDSMVRHQESLQAKQRALLKQLDSDLQEREELRGS
 LDEAEAQRARVEEQQSEREQGCQLRAQQELLQSLQREKQGLEQATTDLRLTILELERELEELKERERL
 LVAFPDLHRPTETQIHGGRSSSVESQITCPTDSGNVTDHMERQVQANDIRIRVLQEENGLRQSMISKIRE
 VAQQGGLKLIPQDRLWSPSSKGTQGATPPVQAKSTSPGPLGRQHLPSRTGRTL LGQPCTSPRPQCTSP
 PRQPCTSPRPQCTSPSRQPCSQPSKSLLEGVTHLDTCTQNP I KVLVRLRKRLSPGRGQASSAHQPQER
 M

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



ACCN: NM_001017437

ORF Size: 2103 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001017437.1](#), [NP_001017437.1](#)

RefSeq Size: 2763 bp


RefSeq ORF: 2259 bp

Locus ID: 550631

UniProt ID: [Q569K6](#)

Cytogenetics: 22q12.2

Product images:



Circular map for RG210351