

Product datasheet for **RG210415**

RHCE (NM_020485) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RHCE (NM_020485) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	RHCE
Synonyms:	CD240CE; RH; Rh4; RH30A; RHC; RHCE(152N); RHE; RhIVb(J); RHIXB; RHNA; RHPI; RhVI; RhVIII
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG210415 representing NM_020485 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCTCTAAGTACCCGCGTCTGTCCGGCGCTGCCTGCCCTCTGCGCCCTAACACTGGAAGCAGCTC
TCATTCTCTCTTCTATTTTTTACCCACTATGACGCTTCCTTAGAGGATCAAAGGGGCTCGTGGCATC
CTATCAAGTTGGCCAAGATCTGACCGTATGGCGGCCATTGGCTTGGGCTTCTCACCTCGAGTTTCCGG
AGACACAGCTGGAGCAGTGTGGCCTTCAACCTTTCATGCTGGCGCTTGGTGTGCAGTGGCAATCTGCG
TGGACGGCTTCTGAGCCAGTTCCTTCTGGGAAGGTGGTCATCACACTGTTCAAGTATTCGGCTGGCCAC
CATGAGTGCTATGTCGGTGTGATCTCAGCGGGTGTCTTGGGAAGGTCAACTTGGCGCAGTTGGTG
GTGATGGTGTGGTGGAGGTGACAGCTTATGGCACCTGAGGATGGTCATCAGTAATATCTTCAACACAG
ACTACCACATGAACCTGAGGCCTTCTACGTGTTGCGAGCCTATTTGGGCTGACTGTGGCTGGTGCCT
GCCAAAGCCTCTACCAAGGGAACGGAGGATAATGATCAGAGAGCAACGATACCCAGTTTGTCTGCCATG
CTGGGCGCCCTCTTCTTGTGGATGTTCTGGCCAAGTGTCAACTCTGCTCTGCTGAGAAGTCCAATCCAAA
GGAAGAATGCCATGTTCAACACCTACTATGCTCTAGCAGTCAGTGGTGACAGCCATCTCAGGGTCATC
CTTGGCTCACCCCAAAGGAAGATCAGCATGACTTATGTGCACAGTGGTGTGGCAGGAGGCGTGGCT
GTGGGTACCTCGTGTACCTGATCCCTTCTCCGTGGCTTGGCATGGTGTGGTCTTGTGGCTGGGCTGA
TCTCCATCGGGGAGCCAAGTGCCTGCCGGTGTGTTGTAACCGAGTGTGGGATTACCCACATCTCCGT
CATGCACTCCATCTTCACTTGTGGTGTGCTTGGAGAGATCACCTACATTGTGCTGCTGGTGTCTCAT
ACTGTCTGGAACGGCAATGGCATGATTGGCTTCCAGGTCCTCCTCAGCATTGGGAACTCAGCTTGGCCA
TCGTGATAGCTCTCACGTCTGGTCTCCTGACAGGTTTGTCTCTAAATCTCAAAATATGGAAGCACCTCA
TGTGGCTAAATATTTGATGACCAAGTTTTCTGGAAGTTTCTCATTGGCTGTTGGATT

ACGGTACGGCGCCGCTCGAG - GFP Tag - GTTTAA



[View online >](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_020485.3](#), [NP_065231.2](#)

RefSeq Size: 1635 bp

RefSeq ORF: 1254 bp

Locus ID: 6006


UniProt ID: [P18577](#)

Cytogenetics: 1p36.11

Protein Families: Transmembrane

Gene Summary: The Rh blood group system is the second most clinically significant of the blood groups, second only to ABO. It is also the most polymorphic of the blood groups, with variations due to deletions, gene conversions, and missense mutations. The Rh blood group includes this gene which encodes both the RhC and RhE antigens on a single polypeptide and a second gene which encodes the RhD protein. The classification of Rh-positive and Rh-negative individuals is determined by the presence or absence of the highly immunogenic RhD protein on the surface of erythrocytes. A mutation in this gene results in amorph-type Rh-null disease. Alternative splicing of this gene results in multiple transcript variants encoding several different isoforms. [provided by RefSeq, Aug 2016]

Product images:



Circular map for RG210415