

Product datasheet for **RG210494**

SLC26A7 (NM_052832) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC26A7 (NM_052832) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SLC26A7
Synonyms:	SUT2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG210494 representing NM_052832
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGACAGGAGCAAAGAGGAAAAAGAAAAGCATGCTTTGGAGCAAGATGCATACCCCCAGTGTGAAGACA
TTATACAGTGGTGTAGAAGGCGACTGCCATTTTGGATTGGGCACCACATTACAATCTGAAAAGAACTT
GCTCCAGACACTGTGTCTGGGATAATGTTGGCAGTTCACAGGTGACCCAAGGATTGGCCTTTGCTGTT
CTCTCATCTGTGCACCCAGTGTGTTGTTTATATGGGTCTCTGTTTCTGCCATAATTTATGCCATATTTG
GAATGGGACATCATGTTGCCACAGGCACCTTTCCTTGACATCCTTAATATCAGCCAACGCCGTGGAACG
GATTGTCCCTCAGAACATGCAGAATCTCACCACACAGAGTAACACAAGCGTGTGGGCTTATCCGACTTT
GAAATGCAAAGGATCCACGTTGCTGCAGCAGTTTCTTCTGGGAGGTGTGATTCAGGTGGCCATGTTTG
TGCTGCAACTGGCAGTCCACATTTGTGGTACAGAGCCTGTGATCAGCGCAATGACAACGGGGCTGC
CACCCATGTGGTACTTACAAGTCAAATATCTCTGGGAATGAAAATGCCATATATATCCGGACCACTT
GGATTCTTTTATATTTATGCATATGTTTTGAAAACATCAAGTCTGTGCGACTGGAAGCATTGCTTTTAT
CCTTGCTGAGCATTGTGGTCTTGTCTTGTAAAGAGCTGAATGAACAGTTTAAAGGAAAATTAAGT
TGTTCTTCTGTAGATTTAGTTTTGATTATTGCTGCATCATTTGCTTGTATTGCACCAATATGAAAAAC
ACATATGGATTAGAAGTAGTTGGTCATATCCACAAGGAATCCCTCACCTAGAGCTCCCCGATGAACA
TCCTCTCTGCGGTGATCACTGAAGCTTTCGGAGTGGCAGTGTAGGCTATGTGGCCTCACTGGCTTTGC
TCAAGGATCTGCCAAAAATCAAATATCAATTGATGACAACCAGGAATTTTGGCCATGGCCTCAGC
AATATAGTTTCTCATTTTTCTTCTGCATACCAAGTGTGCTGCCATGGGAAGGACGGCTGGCCTGTACA
GCACAGGAGCGAAGACACAGGTGGCTTGTCTAATATCTTGCATTTTCGTCCTTATAGTCATCTATGCAAT
AGGACCTTTGCTTTACTGGCTGCCATGTGTCTTGCAGCATTATTGTTGTGGGACTGAAGGGAATG
CTAATACAGTTCGAGATTTAAAAAATATTGGAATGTGGATAAAATCGATTGGGAATATGGGTACGTA
CATATGATTTTACAATATGCTTTGCTGCCAATGTGGGACTGCTGTTTGGTGTGTTGTACCATAGCTAT
AGTGATAGGACGCTTCCCAAGAGCAATGACTGTAAGTATAAAAAATGAAAAGAAATGGAATTTAAAGT
AAGACAGAAATGGACAGTGAACCCTGCAGCAGGTGAAAATTATCTCAATAAACAACCCGCTTGTTTTCC
TGAATGCAAAAAATTTTATACTGATTTAATGAACATGATCCAAAAGGAAAATGCCTGTAATCAGCCACT
TGATGATATCAGCAAGTGTGAACAAAACACATTGCTTAATCCCTATCCAATGGCAACTGCAATGAAGAA
GCTTACAGTCTGCCCTAATGAGAAGTGTATTTAATCCTGGATTGCAGTGGATTTACCTTTTTGACT
ATTCTGGAGTCTCCATGCTTGTGAGGTTTACATGGACTGTAAAGGCAGGAGTGGATGTATTGTTAGC
CCATTGTACAGCTTCTTGATAAAAGCAATGACGTATTATGAAAACCTAGACTCAGAGAAAACCAATTTT
TTTGAATCGGTATCTGCTGCAATAAGTCATATCCATTCAAATAAGAATTTGAGCAAACCTCAGTGACCACA
GTGAAGTC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG210494 representing NM_052832
 Red=Cloning site Green=Tags(s)


MTGAKRKKKSM LWSKMHTPQCEDI IQWCRRLPILDWAPHYNLKENLLPDTVSGIMLAVQQVTQGLAFAV
 LSSVHPVFGLYGLFP AIIY AIFGMGHHVATGTFALTSLISANAVERIVPQNMQNLTTQSNTSVLGLSDF
 EMQRIHVAAAVSFLGGVIQVAMFVLQLGSATFVVTPEVISAMTTGAATHVVT SQVKYLLGMKMPYISGPL
 GFFYIYAYVFENIKSVRLEALLSLLSIVVLVLVKELNEQFKRKIKVVL PVDLVLIIAASFACYCTNMEN
 TYGLEVVGHIPQGIPSPRAPPMNILSAVITEAFGVALVGYVASLALAQGS AKKFKYSIDDNQEFLAHGLS
 NIVSSFFFCIPSAAMGR TAGLYSTGAKTQVACLISCFVLIVIY AIGPLLYWLP MCVLASIIVVGLKGM
 LIQFRDLKKYWNVDKIDWGIWVSTYVFTICFAANVGLLFGVVCTIAIVIGRFPRAMTVS IKNKMEFEKV
 KTEMDSETLQQVKIISINPLVFLNAKKFYDLMNMIQKENACNQLDDISKCEQNTLLNSL SNGNCNEE
 ASQSCPNEKCYLILDCSGFTFFDYSGVSM LVEVYMDCKGRSVDVLLAHCTASLIKAMTY YGNLDSEKPIF
 FESVSAATSHIHSNKNL SKLSDHSEV

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_052832

ORF Size: 1968 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_052832.4](#)

RefSeq Size: 5280 bp

RefSeq ORF: 1971 bp

Locus ID: 115111

UniProt ID: [Q8TE54](#)

Cytogenetics: 8q21.3

Protein Families: Transmembrane

Gene Summary: This gene is one member of a family of sulfate/anion transporter genes. Family members are well conserved in gene structure and protein length yet have markedly different tissue expression patterns. This gene has abundant and specific expression in the kidney. Alternatively spliced transcript variants that encode different isoforms have been described. [provided by RefSeq, Aug 2013]