

Product datasheet for **RG210628**

DDX4 (NM_024415) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DDX4 (NM_024415) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	DDX4
Synonyms:	VASA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG210628 representing NM_024415
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGGAGATGAAGATTGGGAAGCAGAAATCAACCCTCATATGTCTTCTATGTTCCCATATTTGAGAAGG
 ATAGGTATTCTGGAGAAAATGGAGACAATTTAACAGGACTCCAGCTTCATCATCAGAAATGGATGATGG
 ACCTTCTCGAAGAGATCATTTCATGAAAAGTGGATTTGCCTCTGGGCGGAATTTTGGAAACAGAGATGCT
 GGTGAGTGAATAAGCGAGATAATACATCCACAATGGGTGGTTTTGGAGTTGAAAGAGTTTTGGAAACA
 GAGGTTTTTCAAACAGCAGGTTTGAAGATGGTGATAGCTCTGGTTTTCTGGAGAGAGTCTAGTAATGACTG
 CGAAGATAATCCAACACGGAACAGAGGGTTTTCCAAGAGAGGCGGCTATCGAGATGGAAATAATTCAGAA
 GCTTCAGGGCCATACAGAAGAGGTGGAAGAGGTAGTTCCGAGGTTGCCGTGGAGGATTTGGTCTAGGAA
 GTCCAAAATAATGACTTAGACCCAGACGAATGTATGCAGCGCACTGGTGGCCTTTTTGGTTCTAGAAGACC
 AGTATTAAGTGGCACAGGTAATGGTGATACTTCTCAAAGCAGAAGTGGCAGTGGAAAGTGAACGAGGTGGT
 TACAAAGGTTTAAATGAAGAAGTAATAACAGGCTCTGGAAAGAATTCCTGGAAAGTGAAGAAGGAGGAG
 GAGAAAGTAGTGATACTCAAGGACCAAAAGTGACCTACATACCCCTCCTCCACCTGAGGATGAGGACTC
 CATCTTTGCACATTATCAGACAGGCATAAACTTCGACAAATACGACACTATTCTGTGGAAGTGTCTGGA
 CATGATGCACCACCAGCAATTCTGACTTTTGAAGAAGCTAATCTCTGTCAGACACTGAATAACAACATTG
 CTAAGCTGGTTATACTAAGCTTACTCCTGTGCAAAAATACAGTATTCTATCATACTTGCGAGCAGAGA
 TTTGATGGCTTGGCTCAAACAGGGTCTGGGAAGACTGCGGCTTTTCTCTACCAATTTTGGCTCATATG
 ATGCATGATGGAATAACTGCCAGTCGTTTTAAAGAGTTGCAGGAACCAGAGTGTATTATTGTAGCACCAA
 CTCGAGAATTGGTCAACCAGATTTATTTGGAAGCCAGAAAATTTCTTTTGGGACTTGTGAAGAGCTGT
 TGTATATATGGGGGAACCCAGCTGGGACATTCATTTCGACAAATAGTACAAGGCTGTAATATATTATGT
 GCTACTCTGGAAGACTGATGGATATCATAGGCAAAAGAAAAGATTGGTCTCAAACAGATCAAATACTTAG
 TTTTGGATGAAGCTGATCGCATGTTGGATATGGGTTTTGGTCCAGAAATGAAGAAGTTAATTTCTTGCCC
 AGGAATGCCATCAAAGGAACAGCGCCAAACCCTTATGTTCAAGTCAACTTTTCCAGAGGAAATTCAAAGG
 TTGGCTGCAGAGTTTTTAAAGTCAAATTATCTGTTTTGTTGCTGTTGGACAAGTGGGTGGAGCATGTAGAG
 ATGTTCCAGCAGACCGTTCTCAAGTTGGCCAGTCTCAAAAAGAGAAAAGCTCGTTGAAATTTCTGCGAAA
 CATAGGGGATGAAAGAACTATGGTCTTTGTTGAAACTAAGAAAAAGCAGATTTTATTGCAACTTTTCTT
 TGTCAGAAAAAATATCAACTACAAGTATTCATGGTGATCGGGAACAGAGAGAGCGGGAGCAAGCTCTTG
 GAGATTTTCGCTTTGGAAGTGCCAGTCTTGTGTTGCTACTTCAGTAGCTGCCAGAGGGCTGGATATTGA
 AAATGTGCAACATGTTATCAATTTTGTCTTCTTCTACCATTGATGAATATGTTTCATCGAATTTGGCGT
 ACTGGTCGTTGTGGGAATACTGGCAGAGCAATTTCTTTTTTGTCTTGAATCGGATAACCATTTAGCAC
 AGCCTCTAGTAAAAGTATTGACAGATGCTCAACAGGATGTTCTGATGTTGGAAGAAATTCCTTTAG
 TACATACATTCTGGCTTCAAGTGTAGTACAAGAGGAAACGTGTTTGCATCAGTTGATACCAGAAAGGGC
 AAGAGCACTTTGAACACAGCTGGGTTTTCTTCTTCAAGCTCCCAATCCAGTAGATGATGAGTCATGGG
 AT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG210628 representing NM_024415
 Red=Cloning site Green=Tags(s)

```


MGDEDWEAEINPHMSSYVPIFEKDRYSGENDNFNRTPASSEMDDGPSRRDHFMKSGFASGRNFGNRDA
GECNKRDNTSTMGGFVVGKSFGNRGSNSRFEDGDSSGFWRSSNDCEPNPTRNRGFSKRGGYRDGNNSE
ASGPYRRGGRGSFRGCRGGFGLGSPNNDLPDECMQRTGGLFGSRRPVLSGTGNGDTSQSRSGSGSERGG
YKGLNEEVITGSGKNSWKSEAEGGESSDTQGPKVTYIPPPPEDEDSIFAHYQTGINFDKYDTILVEVSG
HDAPPAILTFEEANLCQTLNNNIAKAGYTKLTPVQKYSIPIILAGRDLMACAQTGSGKTA AFLLPILAHM
MHDGITASRFKELQEPECIIVAPTRELVNQIYLEARKFSFGTCVRAVVIYGGTQLGHSIRQIVQGCNILC
ATPGRMLDIIGKEKIGLKQIKYLVLDEADRMLDMGFGPEMKKLISCPGMPKSKEQRQTLMF SATFPEEIQR
LAAEFLKSNYLFVAVGQVGGACRDVQQT VLVQVQF SKREKLV EILRNIGDERTMVFVETKKKADFIATFL
CQEKISTTSIHGDREQREREQALGDFRFGKCPVLVATSV AARGLDIENVQHVINFDLPSTIDEYVHRIGR
TGRCGNTGRAISFFDLESDNHLAQPLVKVLTDAQQDVPALWEEIAFSTYIPGFSGSTRGNV FASVDTRKG
KSTLNTAGFSSSQAPNPVDESWD
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:




ACCN: NM_024415

ORF Size: 2172 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	<p>This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.</p>
Components:	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_024415.3
RefSeq Size:	2471 bp
RefSeq ORF:	2175 bp
Locus ID:	54514
UniProt ID:	Q9NQI0
Cytogenetics:	5q11.2
Gene Summary:	<p>DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This gene encodes a DEAD box protein, which is a homolog of VASA proteins in Drosophila and several other species. The gene is specifically expressed in the germ cell lineage in both sexes and functions in germ cell development. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009]</p>

Product images:



Circular map for RG210628