

Product datasheet for **RG210642**

NDUFS5 (NM_004552) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: NDUFS5 (NM_004552) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: NDUFS5
Synonyms: CI-15k; CI15K
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG210642 representing NM_004552
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCTTTCTTGGACATCCAGAAAAGTTTCGGCCTTAACATAGATCGATGGTTGACAATCCAGAGTGGTG
AACAGCCCTACAAGATGGCTGGTCGATGCCATGCTTTTAAAAAGAATGGATAGAATGTGCACATGGAAT
CGGTTACTCGGCAGAGAAAGAGTGAAGATAGAATATGATGATTTTCGTAGAGTGTTCCTCGGCAG
AAAACGATGAGACGTGCAGGTACCATCAGGAAGCAGCGGGATAAGCTGATAAAGGAAGGAAAGTACACCC
CTCCACCTCACCACATTGGCAAGGGGAGCCTCGGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG210642 representing NM_004552
Red=Cloning site Green=Tags(s)
MPFLDIQKRFGLNIDRWLTIQSGEQPYKMAGRCHAFEKEWIECAHGIGYTRAEKECKIEYDDFVECLLRQ
KTMRRAGTIRKQRDKLIKEGKYTPPHHIGKGEPRP

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_004552

ORF Size: 318 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_004552.3](#)

RefSeq Size: 540 bp

RefSeq ORF: 321 bp

Locus ID: 4725
UniProt ID: [O43920](#)
Cytogenetics: 1p34.3
Protein Pathways: Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease
Gene Summary: This gene is a member of the NADH dehydrogenase (ubiquinone) iron-sulfur protein family. The encoded protein is a subunit of the NADH:ubiquinone oxidoreductase (complex I), the first enzyme complex in the electron transport chain located in the inner mitochondrial membrane. Alternative splicing results in multiple transcript variants and pseudogenes have been identified on chromosomes 1, 4 and 17. [provided by RefSeq, May 2010]

Product images:



Circular map for RG210642