

Product datasheet for RG210750

galectin 9 (LGALS9) (NM 002308) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: galectin 9 (LGALS9) (NM_002308) Human Tagged ORF Clone

Tag: TurboGFP Symbol: galectin 9

Synonyms: HUAT; LGALS9A

Mammalian Cell

Selection:

Neomycin

Vector: pCMV6-AC-GFP (PS100010)

E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide >RG210750 representing NM_002308

Sequence: Red=Cloning site Blue=ORF Green=Tags(s)

 ${\tt TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC}$

GCCGCGATCGCC

ATGGCCTTCAGCGGTTCCCAGGCTCCCTACCTGAGTCCAGCTGTCCCCTTTTCTGGGACTATTCAAGGAG
GTCTCCAGGACGGACTTCAGATCACTGTCAATGGGACCGTTCTCAGCTCCAGTGGAACCAGGTTTGCTGT
GAACTTTCAGACTGGCTTCAGTGGAAATGACATTGCCTTCCACTTCAACCCTCGGTTTGAAGATGGAGGG
TACGTGGTGTGCAACACGAGGCAGAACGGAAGCTGGGGGCCCGAGGAGAGAAGACACACATGCCTTTCC
AGAAGGGGATGCCCTTTGACCTCTGCTTCCTGGTGCAGAGCTCAGATTTCAAGGTGATGGTGAACGGGAT
CCTCTTCGTGCAGTACTTCCACCGCGTGCCCTTCCACCGTGTGGACACCCATCTCCGTCAATGGCTCTGTG
CAGCTGTCCTACATCAGCTTCCAGCCTCCCGGCGTTGGCCTGCCAACCCGGCTCCCATTACCCAGACAG
TCATCCACACAGTGCAGAGCGCCCCTGGACAGATGTTCTCTACTCCCGCCATCCCACCTATGATGTACCC
CCACCCCGCCTATCCGATGCCTTTCATCACCACCATTCTGGGAGGGCTGTACCCATCCAAGTCCATCCTC
CTGTCAGGCACTGTCCTGCCCAGTGCTCAGAGGTTCCACATCAACCTGTGCTCTGGGAACCACACTCGCCT
TCCACCTGAACCCCCGTTTTGATGAGAATGCTGTGGTCCGCAACACCCAGATCGACAACTCCTGGGGGTC
TGAGGAGCGAAGTCTGCCCCGAAAAAATGCCCTTCGTCCGTGGCCAGAGCTTCTCAGTGTGGATCTTTGTGT
GAAGCTCACTGCCTCAAGGTGGCCGTGGATGGTCAGCACCCTTTTGAATACTACCATCGCCTGAGGAACC
TGCCCACCATCAACAGACTGGAAGTGGGGGGGGCGACATCCAGCTTTTGAATACTACCATCGCCTGAGGAACC
TGCCCACCATCAACAGACTGGAAGTGGGGGGGGCGACATCCAGCTTGTGCCCATGTGCAGACA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Protein Sequence:

>RG210750 representing NM_002308 Red=Cloning site Green=Tags(s)

MAFSGSQAPYLSPAVPFSGTIQGGLQDGLQITVNGTVLSSSGTRFAVNFQTGFSGNDIAFHFNPRFEDGG YVVCNTRQNGSWGPEERKTHMPFQKGMPFDLCFLVQSSDFKVMVNGILFVQYFHRVPFHRVDTISVNGSV QLSYISFQPPGVWPANPAPITQTVIHTVQSAPGQMFSTPAIPPMMYPHPAYPMPFITTILGGLYPSKSIL LSGTVLPSAQRFHINLCSGNHIAFHLNPRFDENAVVRNTQIDNSWGSEERSLPRKMPFVRGQSFSVWILC

EAHCLKVAVDGQHLFEYYHRLRNLPTINRLEVGGDIQLTHVQT


TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:





ACCN: NM_002308

ORF Size: 969 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts

of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by

calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.



Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

17q11.2

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 002308.4</u>

 RefSeq Size:
 1739 bp

 RefSeq ORF:
 972 bp

 Locus ID:
 3965

 UniProt ID:
 000182

Domains: Gal-bind lectin


Gene Summary: The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-

cell and cell-matrix interactions. The protein encoded by this gene is an S-type lectin. It is overexpressed in Hodgkin's disease tissue and might participate in the interaction between the H&RS cells with their surrounding cells and might thus play a role in the pathogenesis of this disease and/or its associated immunodeficiency. Multiple alternatively spliced transcript

variants have been found for this gene. [provided by RefSeq, Jul 2008]

Product images:

Cytogenetics:



Circular map for RG210750