

Product datasheet for **RG210801**

KIAA0495 (NM_207306) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: KIAA0495 (NM_207306) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: KIAA0495
Synonyms: MGC138189
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG210801 representing NM_207306
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTGTCTTTTGTCCAGCTCAGCCGCTCGGATCTCGCTGCCACCAGCCTTACTGCACACCCGTGCCACC
CGCCCTTGCCCGTCAGCCTCAGCCTCAGCTCATTCTCTCAGGCAGTCCCAGCATTGGCACGGTACCTCC
TCCCGTGTGGCCACACGTCTCTGCTCCCTGTCAACCTCCTGCCATCAGCACCACCACCGACTTG
TCTGCCGGGAGGATGCAACACCATCTGCCTCCACCGCCACCTTTCAGTGTTCCTGCTTCCAAGTAA
AGATACCAGCAGTGCCCTCAGAGCAGACCAGCCAGAGTTTTCTGAAGCCTCCACAGGGTCTCCCCGG
AGGAGGCCGAGATGCTCTCGTGACTTTGGAGCCGGGGTGGCTGGCCACCTTGGCCTGGGATCTTTGGG
GTGGGTTCTGGGAGCCCGCACTGCTGCAGAGTGTGGTGGATGAGAACAGCTGCTGCTGTACGTGG
TGGAGGACCAGCTGTGTGATGTGGAGCAAGCCTTCAGAGCTGAGCATTGGGCCAACACCAGAGTGGTTC
GGGAACAGGATGCAGACATTGTCCTCAACGACCAGCGGCACTT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG210801 representing NM_207306
Red=Cloning site Green=Tags(s)
MCLSSSAASDLAATSLTAHPCHPPLPRQPQQLILSGSPSIGTVPPPAVGHTSLLPVNPPAISTTTSDL
SAREDATPSASTGHLVFPFAFQVKIPAVPSEQTSQSFSEASHRVLPGGGPRCSRDFGAGVAGHLGLGIFG
VFGSPALLQSVVDENSCCLLYVVEDQLCDVEQAFRAEHLGQHQS SSGTGRHCPQRPAAAL

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_207306

ORF Size: 603 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_207306.2](#), [NP_997189.1](#)


RefSeq Size: 4042 bp

RefSeq ORF: 605 bp

Locus ID: 57212

Cytogenetics: 1p36.32

Product images:



Circular map for RG210801