

Product datasheet for **RG210929**

PHACTR1 (NM_030948) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PHACTR1 (NM_030948) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PHACTR1
Synonyms:	DEE70; dJ257A7.2; EIEE70; RPEL; RPEL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG210929 representing NM_030948
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGGATTATCCCAAATGGATTATTTCTGGATGTAGAGTCTGCTCACAGACTCTTGGATGTTGAGTCAG
 CTCAAAGATTCTTCTACAGTCAAGGAGCTCAAGCTCGCCGGGCGACCCTGCTCCTGCCTCCACATTAAT
 GGCGGCATCCTCGGAGGATGATATAGACCGGCGGCCATCCGGAGAGTGCCTCCAAGAGCGACACGCCG
 TACCTCGCAGAGGCCAGGATCTCCTTTAACCTGGGGCAGCTGAGGAAGTGGAGAGGCTGGCGGCGATGC
 GTTCTGACTCCCTCGTCCAGGCACCCACACCCACCCATCCGCAGGAGAAGTAAGTTTGCCAACTGGG
 AAGGATTTTCAAGCCTTGAAATGGAGGAAGAAGAAAAGCGAAAAGTTCAAACACACGTACGACGCCCTG
 GAAAGGAAATATCTATGAGGCAAAGCAGAGAAGAGCTGATAAAGCGAGGAGTCTGAAGGAAATCTATG
 ATAAAGATGGGGAAGTCTCTATATCCAATGAAGAGGACTCCCTAGAAAATGGGCAGTCCCTGAGCTCCAG
 CCAGCTGTCTCTGCTGCCCTGTCCGAAATGGAGCCAGTCCCAATGCCAGGGATCCCTGCTCATATGAG
 GTGCTCCAACCGTCAGACATCATGGATGGGCCAGATCCTGGCGCCCTGTGAAATTGCCTTGTCTGCCAG
 TGAAACTGTGCGCTCCGTACCTCCAAAGAAAGTCATGATCTGTATGCCCGTGGGGGGGCCAGACCTCTC
 ACTGGTGTCTACACAGCCCAGAAGAGTGGCCAGCAGGGTGTGGCCAGCACCACCACACTGTCTGCC
 TCCCAGATCCAGCACCAGCTGCAGTACGGCAGCCACGGCCAGCACCTCCCTCCACCACCGGCTCCCTCC
 CCATGCACCCCTCGGGCTGCAGAATGATAGACGAGCTCAACAAAACGCTGGCCATGACCATGCAGAGGCT
 GGAAAGCTCTGAGCAGCGGGTCCCCTGTTCCACTTCTTACCACAGCTCTGGGTGCACTCGGGTGATGGG
 GTCACCAAAGCAGGACCTATGGGCCTTCCAGAAATAAGACAAGTGCCAACTGTTGTGATTGAATGTGATG
 ACAATAAAGAAAATGTGCCTCATGAGTCAGACTACGAAGACTCTTCTGCTGTATACAAGAGAAGAGGA
 GGAAGAGGAGGAGGACGAAGACGACGACAGCTCATTATACACCAGCTCCCTGGCCATGAAGGTCTGCAGG
 AAGGACTCCTTAGCCATCAAACCTCAGCAACAGGCCCTCCAAGCGAGAGCTGGAAGAAAAGAACATCCTTC
 CCAGGCAGACGGATGAGGAGCGGCTGGAGCTGAGGCAACAGATTGGCACCAAGCTCACCAGCGGCTGAG
 CCAGAGGCCAACTGCAGAGGAACTGGAACAGAGGAACATTTTGAACCTCGGAATGAACAAGAGGAACAG
 GAGGAGAAGAGAGAGATCAAGAGGAGGCTAACCCGAAAGCTCAGTCAAAGGCCACGGTGAAGAGCTTC
 GGGAAAGAAAGATCCTCATCCGCTTCAGTACTACGTGGAGGTGGCTGACGCTCAGGACTATGACCGCAG
 GGCAGATAAGCCGTGGACCCGCTCACCCTGCAGACAAAGCTGCCATCCGAAAGGAGCTCAATGAATTC
 AAAAGCACTGAGATGGAAGTTCATGAATTGAGTAGACACTTAACAAGTTTCACCGACCT

ACGCGTACGCGGCGGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:


>RG210929 representing NM_030948
 Red=Cloning site Green=Tags(s)

MDYPKMDYFLDVESAHRLLDVESAQRFYFYSQGAQARRATLLLPTLMAASSEDIDRRPIRRVRKSDTP
 YLAEARISFNLGAAEEVERLAAMRSDSLVPGHTPPIRRRSKFANLGRIFPKWKWRKKSEKFKHTSAAL
 ERKISMRSREELIKRGVLKEIYDKDGELSI SNEEDSLENGQSLSSQLSLPALSEMEPVMPRPDPCSYE
 VLQPSDIMDGPDPGAPVKLPCLPVKLSPLPPKVMICMPVGGPDL SLVSYTAQKSGQQGVAQHHTVLP
 SQIQHQLQYGSQGHL PSTTGLPMHPSGCRMIDELNKTAMTMQRLESSEQRVPCSTSYHSSGLHSGDG
 VTKAGPMGLPEIRQVPTVVEICDDNKENVPHESDYEDSSCLYTREEEEEEDEDDSSLYTSSLAMKVCR
 KDSLAIKLSNRPSKRELEENILPRQTDEERLELRQIGTKLTRRLSQRPTAEELQRNLIKPRNEQEEQ
 EEKREIKRRLTRKLSQRPTVEELRERKILIRFSDYVEVADAQDYDRRADKPWTRLTAADKAAIRKELNEF
 KSTEMEVHEL SRHLTRFHRP

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_030948

ORF Size: 1740 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_030948.6](#)

RefSeq Size: 2150 bp

RefSeq ORF: 1743 bp


Locus ID: 221692

UniProt ID: Q9C0D0

Cytogenetics: 6p24.1

Gene Summary: The protein encoded by this gene is a member of the phosphatase and actin regulator family of proteins. This family member can bind actin and regulate the reorganization of the actin cytoskeleton. It plays a role in tubule formation and in endothelial cell survival. Polymorphisms in this gene are associated with susceptibility to myocardial infarction, coronary artery disease and cervical artery dissection. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Apr 2016]

Product images:



Circular map for RG210929