

## Product datasheet for **RG210981**

### OR2H3 (OR2H2) (NM\_007160) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	OR2H3 (OR2H2) (NM_007160) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	OR2H3
Synonyms:	dJ271M21.2; FAT11; hs6M1-12; OLF2R2; OLF2R42B; OR2H3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG210981 representing NM_007160 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGTTAACCAAAGCTCCACACCGGGCTTCCTCCTTCTGGGCTTCTCTGAACACCCAGGGCTGAAAGGA  
CTCTCTTCGTGGTTGTCCTCACTTCCTACCTCCTAACCCCTAGTGGGCAACACACTCATCATCCTGCTGTC  
TGGCTGGACCCCAAGCTCCACTCTCCAATGTACTTTTCTCTCCAACCTCTCCTTCTGGACCTCTGT  
TTCACCACGAGTTGTGTTCCCAAATGCTGGTCAACCTCTGGGGCCAAAGAAGACCATCAGCTTCTCG  
ACTGCTCTGTCCAGATCTTCATCTTCTGTCCCTGGGACAAGTGCATCCTCTTGACAGTGATGGC  
TTTTGATCGCTACGTGGCTGTCTGCCAGCCCTCCACTATGCCACCATCATCCACCCCGCTGTGCTGG  
CAGCTGGCATCTGTGGCCTGGGTCAATTGGGCTAGTGGAGTCAGTGGTCCAGACACCATCCACCTGCACC  
TGCCCTTCTGCCCGATCGGCAGGTGGATGATTTTGTCTGTGAGGTCCCAGCTCTAATTCGACTCTCCTG  
TGAAGACACCTCTACAATGAGATCCAGGTGGCTGTTGCCAGTGTCTTCATCTTGTTGTGCCTCTCAGC  
CTCATCCTTGTCTTTACGGAGCCATTACCTGGCAGTGTGAGGATTAAGTCTGCAAAAGGGCGGAGGA  
AAGCTTTGGGACCTGCTCCTCCCATCTCACTGTGGTCAACCTCTTCTACAGCTCAGTCATTGCTGTCTA  
CCTCCAGCCCAAAAATCCCTATGCCAAGAGAGGGGCAAGTCTTTGGTCTTCTATGCAGTGGGCACT  
CCTTCACTTAACCCTCTCATATACACCCTGAGGAACAAGGAGGTAACCGAGGCATTACAGGAGATTGCTGG  
GGAAGGAAATGGGGCTCACACAAAGC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG210981 representing NM\_007160  
 Red=Cloning site Green=Tags(s)

MVNQSSTPGFLLLGSEHPGLERTLFVVVLTSYLLTLVGNTLIILL SALDPKLSHPMYFFLSNLSFLDLC  
 FTTSCVPQMLVNLWGPKKTISFLDCSVQIFIFLSLGTTECILLTVMAFDRYVAVCQPLHYATIIHPRLCW  
 QLASVAWVIGLVESVVQTPSTLHLPFCPDRQVDDFVCEVPALIRLSCEDTSYNEIQVAVASVFILVVPLS  
 LILVSYGAITWAVLRINSAGRRKAFGTCSHLLTVVTLFYSSVIAVYLPKPNPYAQERKFFGLFYAVGT  
 PSLNPLIYTLRNKEVTRAFRRLLGKEMGLTQS

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_007160

**ORF Size:** 936 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_007160.3](#), [NP\\_009091.3](#)

**RefSeq Size:** 1063 bp

**RefSeq ORF:** 939 bp

**Locus ID:** 7932

**UniProt ID:** [O95918](#)


**Cytogenetics:** 6p22.1

**Protein Families:** Druggable Genome, Transmembrane

**Protein Pathways:** Olfactory transduction

**Gene Summary:** Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

### Product images:



Circular map for RG210981