

Product datasheet for **RG211064**

ZNF182 (NM_001007088) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF182 (NM_001007088) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ZNF182
Synonyms:	HHZ150; KOX14; Zfp182; ZNF21
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG211064 representing NM_001007088
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCAAACCCAGGGCTAGTGACATTTGAAGATGTAGCTGTGGATTTACCCAGGAGGAGTGGCAGT
 ACCTGAACCCACACAGAGGACCCTGTACAGAGACGTGATGCTGGAGACCTACAGCAACTTGGTCTTTGT
 GGGCAGCAGGTTACAAACCAAACCTCATCTCAAGTTGGAGGTAGAAGAATGCCCGCAGAAGGAAAA
 ATCCCATTTTGGAACTTTCCAGAAGTCTGTCAAGTTGATGAACAGATTGAGAGGCAACATCAGGATGACC
 AAGATAAATGTCTGCTGATGCAAGTTGGATTTTCTGACAAGAAAACAATTATACCAAGAGTGCTCGTGA
 CTGTCAATGAGTTTGGAAACATACTTCATCTGAGTACAAACCTTGTGCTTCAATACAAAGACCCGATAAA
 CACGAATCATTTGGAAAATAATGGTAGATAATTTAGACTATTTAGTAGAAGTTCTGCAGAAAAATAAT
 ATGATAATGGATGTGCAAAATTTGTTCTCCATACTGAGTATGAGAAAAACAATCTTGGAAATGAAGCCCTA
 TGGCTATAAAGAGTGTGGAAAGGTCTTAGGCGAAAGAAAGGCCTTAGTCTACATCAGAGAATTAATAAT
 GGAGAGAAACCCTTTGAATGTACTGCATGTAGGAAAACCTTCAGCAAGAAGTACACCTCATTGTACATT
 GGAGAACTCATACGGGAGAGAAACCTTTTGGATGTACCGAATGTGAAAAAGCTTTTAGCCAAAAATCTCA
 GCTCATTATACACTTGGCAACTCATACAGGAGAGAGACCCTTTGGAGTGTCTGAATGTGGAAAAGCCTTC
 AGAGAAAAGTCAACTGTCATCATACATTACAGGACTCACACAGGAGAAAAACCTTATGAATGTAATGAAT
 GTGGAAAAGCCTTCACTCAGAAGTCCAACCTCATTGTCCATCAGAAAACCCACACTGGAGAGAAAACCTA
 TGAATGCACTAAATGTGGAGAATCTTTCATACAGAAGCTTGATCTAATTATACATCATAGTACCCATACA
 GGAAAGAAACCCCATGAATGTAATGAGTGAAGAAAACCTTTCAGTGATAAGTCAACTCTCATTATACACC
 AGAGAACTCATACGGGAGAGAAACCTCATAAATGTACTGAATGTGGGAAGTCTTTCAATGAGAAGTCAAC
 CCTCATTGTGCATCAAAGAACTCATACAGGAGAGAAACCTATGAATGTGATGTGTGGAAAAACCTTC
 ACGCAAAAGTCAAACCTTGGTGTACATCAGAGAACTCATTACAGGAGAGAAACCCTTTGAATGTAATGAAT
 GTGAGAAAAGCCTTCTCTCAGAAGTCTACCTCATGCTGCATCAGAGAGGTATACAGGAGAGAAACCTTA
 CGAGTGAATGAATGTGAAAAAGCATTTTACAGAAATCATATCTCATTATACATCAAAGAACTCATACA
 GAAGAAAACCTTATAAATGTAATGAATGTGGCAAAGCCTTCAGAGAAAAGTCAAAGCTCATTATACATC
 AGAGAATTCATACAGGAGAGAAACCTTATGAATGTCTGTGTGTTGAAAAGCTTTTAGCCAGAAGTCA
 GCTCATAATACATCAGAGAACGCACACGGGAGAGAAACCTTATGCATGCCTGAGTGTGGCAAAGCCTTC
 CGAGAAAAGTCAACATTCCTGTACATCAAAGAACTCATACTGGAGAGAAACCTTATAAATGTACAGAAT
 GTGGAAAAGCCTTACCCAAAAATCAAACCTTATTGTACATCAGCGAACCCATGCAGGAAAAGAAAGCCCA
 TGGAAAGAGCCCACTCGGAAGTCAAAGTTCATGGCACAT

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG211064 representing NM_001007088
 Red=Cloning site Green=Tags(s)

MAKPQGLVTFEDVAVDFTQEWEQYLNPPQRTLRYDVMLETYSNLVFGQVQVTKPNLILKLEVEECPAEGK
 IPFWNFPEVCQVDEQIERQHQQDDQDKLLMQVGFSDKKTIIITKSARDCHEFGNHLHLSTNLVASIQRPDK
 HESFGNNMVDNLDLFSRSSAENKYDNGCAKLFFHTEYEKTNPGMKPYGYKECGKGLRRKGLSLHQRIKN
 GEKPFECTACRKTFSKKSHLIVHWRHTGKPFGEKTECGKAFSQQSQLIIHLRHTHTGERPFECPECGKAF
 REKSTVIIHYRHTGKPYECNECGKAFQKSNLIVHQKTHTEKTYECTKCGESFIQKLDLIIHHSHTHT
 GKKPHECNECKTFSDKSTLIIHQRTHTGKPHKCTECGKSFNEKSTLIVHQRTHTGKPYECDVCGKTF
 TQKSNLGVHQRTHSGKPFECNECEKAFSQQSYLMLHQRGHTGKPYECNECEKAFSQQSYLIIHQRTHT
 EEKPYKNECGKAFREKSKLIIHQRIHTGKPYECPVCKAFSQQSQLIIHQRTHTGKPYACTECGKAF
 REKSTFTVHQRHTGKPYKCTECGKAFTQKSNLIVHQRTHTAGKKAHGRGHTRKSKFMAH

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001007088

ORF Size: 1860 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001007088.2](#)

RefSeq Size: 3504 bp

RefSeq ORF: 1863 bp


Locus ID: 7569

UniProt ID: [P17025](#)

Cytogenetics: Xp11.23

Gene Summary: Zinc-finger proteins bind nucleic acids and play important roles in various cellular functions, including cell proliferation, differentiation, and apoptosis. This gene encodes a zinc finger protein, and belongs to the krueppel C2H2-type zinc-finger protein family. It may be involved in transcriptional regulation. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, May 2010]

Product images:



Circular map for RG211064