

## Product datasheet for **RG211070**

### **TMEM195 (AGMO) (NM\_001004320) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	TMEM195 (AGMO) (NM_001004320) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TMEM195
Synonyms:	TMEM195
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG211070 representing NM\_001004320  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGAAGAACCAGAAGCCAGCAGGATGTTTCAGTTTCCAGGGATTTCCGATGTTGTTTTACACGATGA  
 AACCCAGTGAAACTTCATTCCAAACATTAGAAGAGGTGCCTGATTATGTAAAAAGGCAACTCCATTTTT  
 CATTTCCTTTGATGCTGCTTGAACCTGTTGTCTGAGCTGGATTCTCAAAGGAAAGCCACCAGGTCGCCTGGAT  
 GATGCTTTAACGTCAATCTCAGCTGGTGTCTGTCTCGACTTCCAAGTCTATTTTTCAGGAGCATTGAAC  
 TGACCAGTTATATTTATCTGGGAGAACTACAGGCTGTTCAATTTGCCTGGGATTCTCCATGGACTTG  
 GTATTCAGCCTTCTTAGGAGTTGACTTTGGCTACTACTGGTTCATCGTATGGCTCATGAAGTTAATATT  
 ATGTGGGCGGACACCAACACATCATAGTTCTGAAGACTATAACTTATCCACAGCACTGAGACAGTCTG  
 TCCTCCAGATATACTTCTGGATTTTCTACTCTCCCCTGGCCCTTTCATACCCCTTCAGTATATGC  
 TGTTTCATCTTCAATCAATCTTCTTTACCAATTTTGGATCCATACAGAGGTCATCAATAACCTTGGTCCT  
 TTGGAAGTGAATCTTAATACTCCTAGCCATCATAGGGTTCATCATGGCAGAAATCGTTATTGCATAGACA  
 AAAATTATGCTGGTGTCTTATTATTTGGGATAAAATTTTGGGACATTTGAAGCAGAAAATGAAAAAGT  
 TGTATATGGCTTAACACATCCCATTAATACATTTGAACCAATCAAAGTGCAGTTCATCACTTATTTTCC  
 ATATGGACTACATTCTGGGCCACACCTGGATTCTTCAATAAGTTTTCTGTCAATTTAAGGGACCGGGAT  
 GGGGTCCAGGTAACCAAGACTTGGTCTCAGTGAAGAAATTCAGAGGTCACCGGCAAAGAAGTTCCTT  
 CTCATCATCTTTCAGCTATTAAGATATATACAGTTGTACAGTTGCTCTGATGTTGGCATTATAT  
 GAAGAGACCTTTGAGATACAGCTGCACTGTGCAAGTACTCTCCTTCTGAGGGTTTGCCTCATTATCC  
 TGACCTTGACTTCCATTGGATTTCTTCTGGATCAAAGACCAAGGCAGCTATTATGGAACCTCCGTTG  
 CTTGATGTTCTTAATGCTGTACCGATTTGGTACCTGAAGCCTCTTGTCCCTTATTGTATCTGCTTTT  
 GAGATTGTTTTTCCATTGCTTCTGGGAGTTAGAAGCATGAAACAACTCACCTCTCACCTT  
 GGAAA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:**

>RG211070 representing NM\_001004320  
 Red=Cloning site Green=Tags(s)

MKNPEAQQDVSVSQGFRMLFYTMKPSSETSFTLEEVDPYVKKATPFFISLMLELVVSWILKGKPPGRLD  
 DALTSISAGVLSRLPSLFFRSIELTSYIYIWIENYRLFNLPWDSPWTWYSAFLGVDFGYWYFHRMAHEVNI  
 MWAGHQTHHSSSEDYNLSTALRQSVLQIYTSWIFYSPLALFIPPSVYAVHLQFNLLYQFWIHTEVINLGP  
 LELILNTPSHHRVHHGRNRYCIDKNYAGVLIWIDKIFGTFEAENEKVYGLTHPINTFEPKIKVQFHHLFS  
 IWTTFWATPGFFNKFVIFKPGWGPGRKRLGLSEEIPEVTGKEVPFSSSSQLLKITYTVVQFALMLAFY  
 EETFADTAALSQVTLRLRVCFIILTLSIGFLLDQRPKAAIMETLRCLMFLMLYRFGHLKPLVPSLSSAF  
 EIVFSICIAFWGVRSMKQLTSHPWK

**TRTRPLE** - GFP Tag - V

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**


**ACCN:** NM\_001004320

**ORF Size:** 1335 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001004320.2](#)

RefSeq Size: 2475 bp  
 RefSeq ORF: 1338 bp  
 Locus ID: 392636  
 UniProt ID: [Q6ZNB7](#)  
 Cytogenetics: 7p21.2  
 Protein Families: Transmembrane  
 Gene Summary:

The protein encoded by this gene is a tetrahydrobiopterin- and iron-dependent enzyme that cleaves the ether bond of alkylglycerols. Sequence comparisons distinguish this protein as forming a third, distinct class of tetrahydrobiopterin-dependent enzymes. Variations in this gene have been associated with decreased glucose-stimulated insulin response, type 2 diabetes, and susceptibility to intracranial aneurysms. [provided by RefSeq, Aug 2012]

**Product images:**



Circular map for RG211070