

## Product datasheet for **RG211117**

### **C18orf21 (NM\_031446) Human Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	C18orf21 (NM_031446) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	C18orf21
Synonyms:	HsT3108; PNAS-124; PNAS-131; XTP13
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG211117 representing NM_031446 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGACAGAAGCACTACCTTGAGGCTGCAGCGCGGGACTGCACGACAGCTGCCCGGGCCAAGCCCGCT  
ACCTCCTCTGGGCTACACTTCGTCGCACGATGATAAGAGCACTTTTGAAGAAACGTGTCCATACTGTTT  
CCAGCTGTTGGTCTGGATAACTCTCGAGTGCCTCAAACCCAAAGCCAGGTTGACACCCAAAATACAG  
AACTTCTTAATCGAGAAGCGAGAACTATACACTCAGTTTTAAAGAAGCAAAAATGGTAAAAAGTTCA  
AAGACTCCAAAAGTGTATTGTTGATCACTTGTAACATGCAACAGAACAGTGAAACATCATGGTAAAAG  
TAGAAGCTTTGTGCAACATTGAAGAGCAATCCTGCCACTCCTACAAGTAACTCAGCCTGAAGACCA  
GAGAGAAGGACTGCAAACCCAAATCATGACATGTCTGGCTCGAAAAGCAAGAGCCAGCATCGGTTTTCA  
GAACACCTACATCTGGACAGTCAGTATCTACTTGCTCCTCAAAGAACCAGCAAAAACAAAGAACTT  
CTCTCAACTAAAAATGTTACTTAGTCAGAATGAATCCCAAAGATTCCAAAGGTGGACTTCAGAAATTC  
TTATCTTCTGAAGGTGGACTTTAAAA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG211117 representing NM\_031446  
 Red=Cloning site Green=Tags(s)

MRQKHYLEAAARGLHDSCPGQARYLLWAYTSSHDDKSTFEETCPYCFQLLVLDNSRVRLKPKARLTPKIQ  
 KLLNREARNYTL SFKEAKMVKKFKDSKSVLLITCKTCNRTVKHHGKSRSFVSTLKSNPATPTSKLSLKT  
 ERRRTANPNHDMGSGKSPASVFRTPTSQQSVSTCSSKNTSKTKKHFSQLKMLLSQNESQKIPKVDFRNF  
 LSSLKGGLLK

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_031446

**ORF Size:** 660 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_031446.5](#)

**RefSeq Size:** 1061 bp


**RefSeq ORF:** 663 bp

**Locus ID:** 83608

**UniProt ID:** [Q32NC0](#)

**Cytogenetics:** 18q12.2

**Product images:**



Circular map for RG211117