

Product datasheet for **RG211132**

GLUD1 (NM_005271) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GLUD1 (NM_005271) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	GLUD1
Synonyms:	GDH; GDH1; GLUD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG211132 representing NM_005271
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTACCCTACCTGGGCGAAGCGCTGTTGCTGTCCCGGGCCGGCCCGCTGCCCTGGGCTCGGCGTCCG
 CCGACTCGGCCGGTGTCTGGGCTGGGCCCGGGACAGCCCGCCGCCCGCCCGCAGCCGGGGCTGGCATT
 GGCCGCCCGCGCCACTACAGCGAGGCGGTGGCCGACCGCGAGGACGACCCCAACTTCTTCAAGATGGTG
 GAGGGCTTCTCGATCGCGGCGCCAGCATCGTGGAGGACAAGCTGGTGGAGGACCTGAGGACCCGGGAGA
 GCGAGGAGCAGAAGCGGAACCGGTGCGCGGCATCTGCGGATCATCAAGCCCTGCAACCATGTGCTGAG
 TCTCTCTTCCCATCCGGCGGACGACGGCTCTGGGAGGTATCGAAGGCTACCGGGCCAGCACAGC
 CAGCACCAGCCCTGCAAGGGAGGTATCCGTTACAGCACTGATGTGAGTGTAGATGAAGTAAAAGCTT
 TGGCTTCTGATGACATAACAAGTGTGCAGTGGTTGATGTGCCGTTTGGGGTGTAAAGCTGGTGTAA
 GATCAATCCCAAGAACTATACTGATAATGAATTGAAAAGATCACAAGGAGGTTACCATGGAGCTAGCA
 AAAAAGGGCTTTATTGGTCTGGCATTGATGTGCCTGCTCCAGACATGAGCACAGGTGAGCGGGAGATGT
 CCTGGATCGCTGATACCTATGCCAGCACCATAGGGCACTATGATATTAATGCACACGCCTGTGTTACTGG
 TAAACCCATCAGCCAAGGGGAATCCATGGACGCATCTCTGCTACTGGCCGTGGTGTCTTCCATGGGATT
 GAAAATTCATCAATGAAGCTTCTTACATGAGCATTTTAGGAATGACACCAGGGTTTGGAGATAAAACAT
 TTGTTGTTCAAGGATTTGGTAATGTGGGCTACACTCTATGAGATATTTACATCGTTTTGGTGCTAAATG
 TATTGCTGTTGGTGAAGTCTGATGGGAGTATATGGAATCCAGATGGTATTGACCCAAAGAACTGGAAGAC
 TTCAAATTGCAACATGGGTCCATTCTGGGCTTCCCAAGGCAAAGCCCTATGAAGGAAGCATCTTGGAGG
 CCGCATGTGACATACTGATCCCAGCTGCCAGTGAGAAGCAGTTGACCAAATCCAACGCACCCAGGATCAA
 AGCCAAGATCATTGCTGAAGGTCCAATGGGCCAACAACCTCAGAAAGCTGACAAGATCTTCTGGAGAGA
 AACATTATGTTATTCCAGATCTCTACTTGAATGCTGGAGGAGTGACAGTATCTTACTTTGAGTGGCTGA
 AGAATCTAAATCATGTCAGCTATGGCCGTTGACCTTCAAATATGAAAGGGATTCTAACTACCACTTGCT
 CATGTCTGTTCAAGAGAGTTTAGAAAGAAAATTTGAAAGCATGGTGGAACTATTTCCATTGTACCCACG
 GCAGAGTTCCAAGACAGGATATCGGGTGCATCTGAGAAAGACATCGTGCCTCTGGCTTGGCATAACAA
 TGGAGCGTTCTGCCAGGCAAATATGCGCACAGCCATGAAGTATAACCTGGGATTGGACCTGAGAACAGC
 TGCCTATGTTAATGCCATTGAGAAAGTCTTCAAAGTGTACAATGAAGCTGGTGTGACCTTCA

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG211132 representing NM_005271
 Red=Cloning site Green=Tags(s)

MYRYLGEALLSRAGPAALGSASADSAALLGWARGQPAAPQPGLALAARRHYSEAVADREDDPNFFKMV
 EGFDRGASIVEDKLVEDLRTRESEEQKRNVRGILRIKPCNHVLSLSPFIRDDGSWEVIEGYRAQHS
 QHRTPCGGIRYSTDVSVDEVKALASLMTYKCAVVDVPPGGAKAGVKINPKNYTDNELEKITRRFTMELA
 KKGFIGPGIDVPAPDMSTGEREMSWIADTYASTIGHYDINAHACVTGKPISSQGGIHRISATGRGVFHGI
 ENFINEASYMSILGMPGFGDKTFVVQGFVNLHSMRYLHRFGAKCIAVGESDGIWNPDPGIDPKELED
 FKLQHGSI LGFPKAKPYEGSILEADCDILIPAASEKQLTKSNAPRVKAKIIAEGANGPTTPEADKIFLER
 NIMVIPDLYLNAGGVTVSYFEWLKLNHVSYGRLTFKYERDSNYHLLMSVQESLERKFGKHGGTIPVPT
 AEFQDRISGASEKDIVHSLAYTMERSARQIMRTAMKYNLGLDLRTAAYVNAIEKVFVYNEAGVTFT

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_005271

ORF Size: 1674 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	NM_005271.5
RefSeq Size:	3051 bp
RefSeq ORF:	1677 bp
Locus ID:	2746
UniProt ID:	P00367
Cytogenetics:	10q23.2
Domains:	GLFV_dehydrog, GLFV_dehydrog_N
Protein Families:	Druggable Genome
Protein Pathways:	Alanine, aspartate and glutamate metabolism, Arginine and proline metabolism, D-Glutamine and D-glutamate metabolism, Metabolic pathways, Nitrogen metabolism
Gene Summary:	This gene encodes glutamate dehydrogenase, which is a mitochondrial matrix enzyme that catalyzes the oxidative deamination of glutamate to alpha-ketoglutarate and ammonia. This enzyme has an important role in regulating amino acid-induced insulin secretion. It is allosterically activated by ADP and inhibited by GTP and ATP. Activating mutations in this gene are a common cause of congenital hyperinsulinism. Alternative splicing of this gene results in multiple transcript variants. The related glutamate dehydrogenase 2 gene on the human X-chromosome originated from this gene via retrotransposition and encodes a soluble form of glutamate dehydrogenase. Related pseudogenes have been identified on chromosomes 10, 18 and X. [provided by RefSeq, Jan 2016]

Product images:



Circular map for RG211132