

Product datasheet for **RG211443**

HNF 4 alpha (HNF4A) (NM_001030003) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HNF 4 alpha (HNF4A) (NM_001030003) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	HNF 4 alpha
Synonyms:	FRTS4; HNF4; HNF4a7; HNF4a8; HNF4a9; HNF4alpha; MODY; MODY1; NR2A1; NR2A21; TCF; TCF-14; TCF14
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG211443 representing NM_001030003
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGTCAGCGTGAACCGCCCCCTCGGGCTCCAGTGGAGAGTTCTTACGACACGTCCCCATCAGAAGGCA
 CCAACCTCAACGCGCCCAACAGCCTGGGTGTCAGCGCCCTGTGTGCCATCTGCGGGGACCGGCCACGGG
 CAAACTACTCGGTGCCTCGAGCTGTGACGGCTGCAAGGGCTTCTCCGGAGGAGCGTGCGGAAGAACCAC
 ATGTACTCCTGCAGATTTAGCCGGCAGTGCCTGGTGGACAAGACAAGAGGAACCAGTGCCGCTACTGCA
 GGCTCAAGAAATGCTTCCGGGTGGCATGAAGAAGGAAGCCGTCAGAAATGAGCGGGACCGGATCAGCAC
 TCGAAGGTCAAGCTATGAGGACAGCAGCCTGCCCTCCATCAATGCGCTCCTGCAGCGGAGGTCTGTCC
 CGACAGATCACCTCCCCGTCTCCGGGATCAACGGCGACATTCCGGCGAAGAAGATTGCCAGCATCGCAG
 ATGTGTGTGAGTCCATGAAGGAGCAGCTGCTGGTTCTCGTTGAGTGGGCCAAGTACATCCAGCTTTCTG
 CGAGCTCCCCCTGGACGACCAGGTGGCCCTGCTCAGAGCCATGCTGGCGAGCACCTGCTGCTCGGAGCC
 ACCAAGAGATCCATGGTGTTCAGGACGTGCTGCTCCTAGGCAATGACTACATTGCTCCCTCGGCACTGCC
 CCGAGCTGGCGGAGATGAGCCGGGTGCCATACGCATCCTTGACGAGCTGGTGTGCCCTTCCAGGAGCT
 GCAGATCGATGACAATGAGTATGCCTACCTCAAAGCCATCATCTTCTTTGACCCAGATGCCAAGGGGCTG
 AGCGATCCAGGGAAGATCAAGCGGCTGCGTTCCAGGTGCAGGTGAGCTTGAGGACTACATCAACGACC
 GCCAGTATGACTCGCGTGGCCGCTTTGGAGAGCTGCTGCTGCTGCCACCTTGACAGCATCACCTG
 GCAGATGATCGAGCAGATCCAGTTCATCAAGCTCTTCGGCATGGCCAAGATTGACAACCTGTTGCAGGAG
 ATGCTGCTGGGAGGGTCCCCAGCGATGCACCCATGCCACCACCCCTGCACCCTCACCTGATGCAGG
 AACATATGGGAACCAACGTCATCGTTGCCAACACAATGCCCACTCACCTCAGAACGAGATGTCCAC
 CCCTGAGACCCACAGCCCTCACCGCCAGGTGGCTCAGGGTCTGAGCCCTATAAGCTCCTGCCGGGAGCC
 GTGCCACAATCGTCAAGCCCTCTCTGCCATCCCCAGCCGACCATCACCAAGCAGGAAGTTATC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG211443 representing NM_001030003
 Red=Cloning site Green=Tags(s)

MVSVNAPLGAPVESSYDTSPSEGTNLNAPNSLGVSAICAIGDRATGKHYGASSCDGCKGFFRRSVRKNH
 MYSCRFSRQCVDKDKRNQCRYCRLKKCFRAGMKKEAVQNERDRISTRSSYEDSSLPSINALLQAEVLS
 RQITSPVSGINGDIRAKKIASIADVCEMKEQLLVLVEWAKYIPAFCELPDQVALLRAHAGEHLLLGA
 TKRSMVFKDVLVLLGNDYIVPRHCPELAEMSRVSRILDELVLPFQELQIDDNEYAYLKAIIFDPDAKGL
 SDPGKIKRLRSQVQSLEDYINDRQYDSRGRFGELLLLLPTLQSIWQMIQIQIFIKLFGMAKIDNLLQE
 MLLGGSPSDAPHAHPLHPLMQEHMGTNVIVANTMPHTLSNGQMSTPETPQSPPPGGSGSEPYKLLPGA
 VATIVKPLSAIPQPTITKQEVI

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001030003

ORF Size: 1326 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001030003.3](#)

RefSeq Size: 1339 bp

RefSeq ORF: 1329 bp

Locus ID: 3172

UniProt ID: [P41235](#)


Cytogenetics: 20q13.12

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Nuclear Hormone Receptor, Transcription Factors

Protein Pathways: Maturity onset diabetes of the young

Gene Summary: The protein encoded by this gene is a nuclear transcription factor which binds DNA as a homodimer. The encoded protein controls the expression of several genes, including hepatocyte nuclear factor 1 alpha, a transcription factor which regulates the expression of several hepatic genes. This gene may play a role in development of the liver, kidney, and intestines. Mutations in this gene have been associated with monogenic autosomal dominant non-insulin-dependent diabetes mellitus type I. Alternative splicing of this gene results in multiple transcript variants encoding several different isoforms. [provided by RefSeq, Apr 2012]

Product images:



Circular map for RG211443