

## Product datasheet for **RG211488**

### **KCNIP4 (NM\_025221) Human Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	KCNIP4 (NM_025221) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	KCNIP4
Synonyms:	CALP; KCHIP4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG211488 representing NM_025221 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAATGTGAGGAGGGTGGAAAGCATTTTCGGCTCAGCTGGAGGAGGCCAGCTCTACAGGCGGTTTCCTGT  
ACGCTCAGAACAGCACCAAGCGCAGCATTAAAGAGCGGCTCATGAAGCTTTGCCCTGCTCAGCTGCCAA  
AACGTCGTCCTGCTATTCAAACAGCGTGGAAGATGAACTGGAGATGGCCACCGTCAGGCATCGGCCT  
GAAGCCCTTGAGCTTCTGGAAGCCAGAGCAAATTTACCAAGAAAGAGCTTCAGATCCTTACAGAGGAT  
TTAAGAATGAATGCCCCAGTGGTGTGTTAATGAAGAAACCTTCAAAGAGATTTACTCGCAGTTCTTTCC  
ACAGGGAGACTCTACAACATATGCACATTTTCTGTTCAATGCATTTGATACAGACCAACAATGGAGCTGTG  
AGTTTCGAGGATTTATCAAAGGTCTTTCCATTTTGTCTCCGGGGACAGTACAAGAAAACTCAATTGGG  
CATTTAATCTGTATGACATAAATAAAGATGGCTACATCACTAAAGAGGAAATGCTTGATATAATGAAAGC  
AATATACGATATGATGGGTAAATGTACATATCCTGTCTCAAAGAAGATGCTCCCAGACAACACGTTGAA  
ACATTTTTTCAGAAAATGGACAAAAATAAAGATGGGGTTGTTACCATAGATGAGTTTCATTGAAAGCTGCC  
AAAAAGATGAAAACATAATGCGCTCCATGCAGCTCTTTGAAAATGTGATT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG211488 representing NM\_025221  
Red=Cloning site Green=Tags(s)


MNVRVESISAQLEEASSTGGFLYAQNSTKRSIKERLMKLLPCSAAKTSSPAIQNSVEDELEMATVRHRP  
 EALELLEAQSKFTKKELQILYRFGKNECPSGVVNEETFKEIYSQFFPQGDSTTYAHFLFNAFDTDHNGAV  
 SFEDFIKGLSILLRGTVQEKLNWAFNLYDINKDGYITKEEMLDIMKAIYDMMGKCTYPVLKEDAPRQHV  
 TFFQKMDKNKDGVVTTIDEFIESCQKDENIMRSMQLFENI

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**



**ACCN:** NM\_025221

**ORF Size:** 750 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_025221.6</a>
<b>RefSeq Size:</b>	1556 bp
<b>RefSeq ORF:</b>	753 bp
<b>Locus ID:</b>	80333
<b>UniProt ID:</b>	<a href="#">Q6PIL6</a>
<b>Cytogenetics:</b>	4p15.31-p15.2
<b>Domains:</b>	EFh
<b>Protein Families:</b>	Druggable Genome, Ion Channels: Other
<b>Gene Summary:</b>	This gene encodes a member of the family of voltage-gated potassium (Kv) channel-interacting proteins (KCNIPs), which belong to the recoverin branch of the EF-hand superfamily. Members of the KCNIP family are small calcium binding proteins. They all have EF-hand-like domains, and differ from each other in the N-terminus. They are integral subunit components of native Kv4 channel complexes. They may regulate A-type currents, and hence neuronal excitability, in response to changes in intracellular calcium. This protein member also interacts with presenilin. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]