

Product datasheet for **RG211690**

Kv beta 2 (KCNAB2) (NM_003636) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Kv beta 2 (KCNAB2) (NM_003636) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Kv beta 2
Synonyms:	AKR6A5; HKvbeta2; HKvbeta2.1; HKvbeta2.2; KCNA2B; KV-BETA-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG211690 representing NM_003636 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTATCCAGAATCAACGACGGGCTCCCCGGCTCGGCTCTCGCTGCGGCAGACGGGCTCCCCGGGATGA
TCTACAGTACTCGGTATGGGAGTCCCAAAGACAGCTCCAGTTTTACAGGAACCTGGGCAAGTCTGGCCT
GCGGGTCTCCTGCCTGGGACTTGAACATGGGTGACCTTCGGAGGCCAGATCACCGATGAGATGGCAGAG
CAGCTCATGACCTTGGCCTATGATAATGGCATCAACCTCTCGATACAGCAGAAGTCTACGCAGCCGGCA
AGGCTGAAGTGGTACTGGGAAACATCATTAAAGAAGAAAGGATGGAGGCGTCCAGCCTCGTCATCACCAC
CAAGATCTTCTGGGGCGAAAGGCGGAGACGGAGCGGGGCTGTCCAGGAAGCACATAATCGAAGGTCTG
AAAGCTTCCCTGGAGCGACTGCAGCTGGAGTACGTGGATGTGGTGTGGCAACCGCCGGACCCCAACA
CCCCGATGGAAGAGACCGTCCGCGCCATGACCCACGTCATCAACCAGGGGATGGCCATGTACTGGGGCAC
GTCACGCTGGAGCTCCATGGAGATCATGGAGGCCACTCCGTGGCCCGGCAGTTCAACCTGACCCCGCCC
ATCTGCGAGCAGGCTGAGTACCACATGTTCCAGCGTGAGAAAGTGGAGGTGCAGCTGCCGGAGCTGTTCC
ACAAGATAGGAGTGGGCGCCATGACCTGGTCCCCTCTGGCCTGTGGCATTGTTTCTGGCAAGTACGACAG
TGGCATCCCACCCTACTCAAGAGCCTCCTTGAAGGGCTACCAGTGGCTGAAGGACAAGATCCTCAGTGAG
GAGGGCCGCGCCAGCAAGCCAAGCTGAAGGAGCTGCAGGCCATCGCCGAGCGCCTGGGTCACCCCTGC
CCAGCTGGCCATAGCCTGGTGCCTGAGGAATGAGGGAGTCAGCTCCGTGCTCCTGGGGGCTCCAATGC
GGACCAGCTCATGGAGAACATTGGGCAATACAGTCCCTTCCGAAACTGTCGTCTTCCATTATCCACGAG
ATTGATAGTATTTGGGCAATAAACCTACAGCAAAAAGGACTACAGATCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG211690 representing NM_003636
 Red=Cloning site Green=Tags(s)

MYPESTTGSPARLSLRQTGSPGMIYSTRYGSPKRQLQFYRNLGKSGLRVSLGLGTWVTFGGQITDEMAE
 QLMTLAYDNGINLFDTAEVYAAGKAEVVLGNIKIKKGWRRSSLVITTKIFWGGKAETERLSRKHIIEGL
 KASLERLQLEYVDVVFANRPDPNTPMEETVRAMTHVINQGMAMYWGTSRWSSMEIMEAYSVARQFNLTTP
 ICEQAEYHMFQREKVEVQLPELFHKIGVGAMTWSPLACGIVSGKYDSGIPPYSRASLKGQWLKDKILSE
 EGRRQQAKLKELQAI AERLGCTLPQLAI AWCLRNEGVSSVLLGASNADQLMENIGAIQVLPKLSSSIIE
 IDSILGNKPYSKKDYRS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_003636

ORF Size: 1101 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_003636.1
RefSeq Size:	3171 bp
RefSeq ORF:	1104 bp
Locus ID:	8514
UniProt ID:	Q13303
Cytogenetics:	1p36.31
Domains:	aldo_ket_red
Protein Families:	Druggable Genome, Ion Channels: Other
Gene Summary:	<p>Voltage-gated potassium (Kv) channels represent the most complex class of voltage-gated ion channels from both functional and structural standpoints. Their diverse functions include regulating neurotransmitter release, heart rate, insulin secretion, neuronal excitability, epithelial electrolyte transport, smooth muscle contraction, and cell volume. Four sequence-related potassium channel genes - shaker, shaw, shab, and shal - have been identified in Drosophila, and each has been shown to have human homolog(s). This gene encodes a member of the potassium channel, voltage-gated, shaker-related subfamily. This member is one of the beta subunits, which are auxiliary proteins associating with functional Kv-alpha subunits. This member alters functional properties of the KCNA4 gene product. Alternative splicing of this gene results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Dec 2010]</p>

Product images:



Circular map for RG211690