

Product datasheet for **RG211713**

ORC4L (ORC4) (NM_181742) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ORC4L (ORC4) (NM_181742) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ORC4L
Synonyms:	ORC4L; ORC4P
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG211713 representing NM_181742 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGCAGTCGTAATCAAAGAGTAACAGCTTAATTCACACAGAGTGCCTTTCACAGGTACAAAGAATTT
 TACGTGAAAGATTTTGTCTGAGAGTCCACATAGTAACCTATTTGGAGTGCAGTACAATACAAACACTT
 AAGTGAGCTGCTGAAAAGAACTGCTCTCCATGGAGAGAGTAACCTGTCTTATTATCGGACCCCGAGGA
 TCAGGAAAACTATGTTAATAAATCATGCTTTGAAAGAACTCATGAAATAGAAGAAGTGAGTGAAAATG
 TATTACAAGTTCACTTAAATGGACTGCTGCAGATCAATGACAAAATCGCCCTAAAGGAAATCACAAGGCA
 GTTAAATCTGGAAATGTAGTTGGAGATAAAGTTTTTGAAGCTTTGCTGAAAACCTTTCATTTCTCTG
 GAAGCTTTAAAAAAGGTGACCGAACTAGCAGTTGCCAGTGATCTTCATATTAGATGAATTTGATCTTT
 TTGCTCATCATAAAAAACCAACACTTCTCTATAATCTTTTTGACATTTCTCAGTCTGCACAGACCCCAAT
 AGCAGTTATTGGTCTTACATGTAGATTGGATATTTTGAAGCTCTAGAAAAAGAGTGAAGTCAAGATTT
 TCTCACCAGGAGATACACTTAATGAATTCATTTGGTTTTCCACAGTATGTTAAAAATTTAAAGAACAGT
 TATCTCTACCTGCAGAGTTTCCAGACAAGGTTTTTGTGAGAAGTGAATGAAAATGTTTCAGTATCTCTC
 AGAAGATAGAAGTGTGCAAGAAGTACTACAGAAGCATTTCAATATCAGCAAAAACCTGCGGTCAATACAC
 ATGCTATTGATGCTTGCTTTAAATCGAGTAACAGCATCGCACCCATTTATGACTGCCGTAGATCTAATGG
 AAGCAAGCCAACCTGTGTAGCATGGACTCGAAAGCAAATATTGTACATGGTCTATCAGTCTTGAAAACTCG
 TCTTATAATAGCAATGAACATTTAAATGACATCTATGAGGAAGAGCCATTTAATTTTCAAATGGTCTAT
 AATGAGTTTCAGAAGTTTGTTCAAAGGAAAGCACATTCGGTTTATAATTTTGAAGAACTGTTGTCATGA
 AGGCTTTTGAACACTTGACGCAATTAGAATTAATAAGCCCATGGAAAGAACTTCAGGAAATTCACAGAG
 AGAGTACCAGCTGATGAACTGCTTTTGGATAATACTCAAATTATGAATGCTCTGCAGAAATATCCCAAC
 TGTCTACAGATGTGAGGCAGTGGGCAACATCCTCACTAAGCTGGTTA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG211713 representing NM_181742
 Red=Cloning site Green=Tags(s)


MSSRKSKSNSLIHTECLSQVQRILRERFCRQSPHSNLFQVQVQYKHLSELLKRTALHGESNSVLIIGPRG
 SGKTMLINHAKELMEIEEVSENVLQVHLNGLLQINDKIALKEITRQLNLENNVGDVFGSFAENLSFLL
 EALKKGDRTS SCPVIFILDEFDLFAHHKNQTLLYNLFDISQSAQTPIAVIGLTCRLDILELLEKRVKSRF
 SHRQIHLMNSFGFPQYVKIFKEQLSLPAEFPDKVFAEKWNENVQYLSEDRSVQEVLPQKHFNISKNLRLSH
 MLLMLALNRVTASHPFMTAVDLMEASQLCSMDSKANIVHGLSVLEICLIAMKHLNDIYEEEPFNFMVY
 NEFQKFVQRKAHSVYNFEKPVVMKAFAEHLQQLELIKPMERTSGNSQREYQLMKLLLDNTQIMNALQKYPN
 CPTDVRQWATSSLSWL

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_181742

ORF Size: 1308 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_181742.3](#)

RefSeq Size: 2795 bp

RefSeq ORF: 1311 bp

Locus ID: 5000


UniProt ID: [O43929](#)

Cytogenetics: 2q23.1

Protein Pathways: Cell cycle

Gene Summary: The origin recognition complex (ORC) is a highly conserved six subunit protein complex essential for the initiation of the DNA replication in eukaryotic cells. Studies in yeast demonstrated that ORC binds specifically to origins of replication and serves as a platform for the assembly of additional initiation factors such as Cdc6 and Mcm proteins. This gene encodes a subunit of the ORC complex. Several alternatively spliced transcript variants, some of which encode the same protein, have been reported for this gene. [provided by RefSeq, Oct 2010]

Product images:



Circular map for RG211713