

Product datasheet for **RG211782**

GPATCH4 (NM_015590) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	GPATCH4
Synonyms:	GPATC4
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG211782 representing NM_015590
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTTCCGCGGAGCGGACCCAGCAATGACGAGGTCAAGAGTCGTGGGATGAAGTTTGCTGAGGAGCAGC
TGCTAAAGCATGGATGGACTCAAGGCAAAGGCCTCGGCCGGAAGGAGAATGGTATCACTCAGGCTCTCAG
GGTGACACTGAAGCAAGACACTCATGGGTAGGACATGACCTGCCAAGGAGTTCACAAACCACTGGTGG
AATGAGCTCTTCAACAAGACTGCGGCCAACTTGGTAGTGAAACTGGGCAGGATGGAGTACAGATAAGGA
GCCTTTCTAAGGAGACCACCGTTATAATCATCCCAAGCCCACTTGTGTATCAGAAGTTTGTGAAGAT
GGCTACATTGACTTCAGGTGGAGAGAAGCCAAACAAGACTTGAGAGCTGCAGTGATGACGACAACCAAG
GGGTCCAAGTCCCCAAAGATTCTGACTGATGAGATGCTGCTCCAAGCCTGTGAGGGGCGAACAGCACACA
AGGCTGCCCGTCTTGGGATCACAATGAAGGCCAAGCTTGTCTCGCTAGAGGCCAGGAGCAGGCCTTCTCT
GGCTCGTCTCAAAGGCCAGGACCCTGGGGCCCCTCAACTGCAGTCAGAGAGCAAGCCCCCAAAAAAAG
AAAAAGAAAAGGAGGCAGAAAAGAGGAGGAAGAAGCTACAGCATCTGAAAGGAATGATGCAGATGAGAAGC
ACCCAGAACATGCTGAGCAGAACATCAGAAAAAGCAAGAAGAAAGGAGGACATCAAGAAGGAAAGGT
CTCAGATGAAAGAGAGGGTACAATAAGGGAATGAGAAGGAGGACGCTGCAGGAACAAGTGGGCTTGGG
GAATTGAATAGCAGAGAGCAAAACCAATCAGTCCCTCAGGAAAGGGAAGAAAAAGAGAGGTGGCACCATG
AAGAGGAGAAGATGGGGTCTTGGAGGAAGGAGGAAAGGCAAGGAGGCTGCAGGCAGTGTCAGGACAGA
GGAGGTAGAGAGCAGGGCATATGCTGACCCATGCAGCCGAAGAAAGAGAGGCAGCAACAGGAGGAGGAG
GACTTGAACCTAGAAGATAGAGGTGAGGAACTTTAGGTGGTGAACAGGGAAGCAGAGAGCAGAGCA
TGCAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG211782 representing NM_015590
Red=Cloning site Green=Tags(s)

MFPRSGPSNDEVKSRGMKFAEEQLLKHWGTQGKGLGRKENGITQALRVTLKQDTHGVGHDPAGEFTNHHW
NELFNKTAANLVVETGQDGVQIRSLKETTRYNHPKPNLLYQKFVKMATLTSGGEKPNKDLESCSDDDNQ
GSKSPKILTDEMLLQACEGRTAHKAARLGITMKAKLARLEAQEAFLARLKGQDPGAPQLQSESKPPKKK
KKKRRQKEEEEATASERNDADKHEHAEQNIKRSKKKKRRHQEGKVSDEREGTTKGNEKEDAAGTSLG
ELNSREQTNQSLRKGGKKRWHHEEKMGVLEEGKGKEAAGSVRTEEVESRAYADPCSRKKRQQQEEE
DLNLEDRGEEFRWWNQSRQSMQ

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_015590

ORF Size: 1125 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_015590.3](#), [NP_056405.2](#)


RefSeq Size: 1906 bp

RefSeq ORF: 1128 bp

Locus ID: 54865

Cytogenetics: 1q22-q23.1

Product images:



Circular map for RG211782