

Product datasheet for **RG211789**

OR5L2 (NM_001004739) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: OR5L2 (NM_001004739) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: OR5L2
Synonyms: HSHTPCR16; HTPCR16; OR11-153
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG211789 representing NM_001004739
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCAAGGAAAACGCACCACTGTGGCTGAGTTCATTCTCCTTGGACTATCAGATGTCCTGAGTTGA
GAGTCTGCCTCTCCTGCTGTTCTTCTCATCTATGGAGTCACGTTGTTAGCCAATCTGGGCATGACTGC
ACTGATTCAGGTCAGCTCTCGGCTCCACACCCCGTACTTTTTCTCAGCCACTTGCCTTTGTAGAT
TTCTGCTACTCCTCAATAATTGTGCCAAAGATGTTGGCTAATATCTTTAACAAGGACAAAGCCATCTCT
TCCTAGGGTGCATGGTGAATTCTACTTGTGTTTGCACATGTGGAGTCACTGAGGCTTCTCTGCTGGCCGT
GATGGCCTATGACCGCTTGTGGCCATCTGTAACCCCTGCTGTACATGGTGACCATGTCTCAGAAGCTG
CGTGTGGAGCTGACCTCTTGCTGCTACTTCTGTGGACGGTGTGTTCTCTGATCACTCGTCCTTAGCTC
TTAGGATCCTCTTCTATAGATCTAATGTGATTAACCACTTCTTCTGTGATCTACCCCTCTCCTAAGTCT
TGCTTGCTCTGATGTCACCTGTGAATGAGACACTGCTGTTCTGCTGGCCACTTTGAATGAGAGTGTACC
ATCATGATCATCCTCACCTCCTACCTGCTAATTCTCACCCTATCCTGAAGATACACTCTGCAGAGAGCA
GGCACAAGCTTTCTCCACCTGTGCTCCACCTCACAGCCATCACTGTCTCCATGGAACAATCCTTTA
CATTTATTGCAGGCCGAGTTCAGGCAACAGTGGAGATGTTGACAAAGTGGCCACCGTGTCTACACAGTT
GTGATCCCATGCTGAACCCCTGATCTACAGCCTGAGAAATAAGGATGTGAACAAGCTCTCAGAAAAG
TGATGGGCTCCAAAATCACTCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG211789 representing NM_001004739
Red=Cloning site Green=Tags(s)

MGKENCTTVAEFILLGLSDVPEL RVCLFLLFLLIYGVTL LANLGMTALIQVSSRLHTPVYFFLSHLSFVD
 FCYSSIIVPKMLANIFNKDKAISFLGCMVQFYLFCTCGVTEVFLAVMAYDRFVAICNPLLYMVTMSQKL
 RVELTSCCYFCGTVCSLIHSSLALRILFYRSNVINHHFFCDLP LLSLACSDVTVNETLLFLVATL NESVT
 IMIILTSYLLILTTILKIHSAESRHKAFSTCASHLTAITVSHGTILYICRPSSGNSGDV DKKVATV FYTV
 VIPMLNPLIYSLRNKDVNKALRKVMGSKIHS

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_001004739

ORF Size: 933 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001004739.1 , NP_001004739.1
RefSeq Size:	936 bp
RefSeq ORF:	936 bp
Locus ID:	26338
UniProt ID:	Q8NGL0
Cytogenetics:	11q12.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Olfactory transduction
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]