

Product datasheet for **RG211818**

CCNYL1 (NM_152523) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	CCNYL1
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PSI00010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG211818 representing NM_152523
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGCCCGAAGATTTAGCTTTGGAGTCAAACCTTCTGACCATCCAAGGGCAAGCACAAATTTCTGAGCA
AATCTCAAACGGATGTGCGAGAAAAGAGGAAGACCAATTTGAACCATGTATCTCCAGGGCAGCTTAC
TAAAAAGTATAGCTCATGCTCAACAATATTTCTAGATGACAGCACAGTCAGCCAGCCTAATCTTAGAACC
ACAGTAAATGTGTGACCTTAGCAATATATTACCACATAAAGAACAGAGATGCAATAGATCCCTGGATA
TTTTTGATGAGAGATCACATCCACTTACACGAGAAAAGTTCCAGAGGAATACTTTAAGCATGATCCTGA
GCACAAATTTATTACAGATTTGTTCTGACTCTTTTAGTGCTGCACAGCTAACAGCTGAATGTGCAATA
GTAACCTTGGTTTACTTAGAAAGGCTTTTAACTTATGCTGAAATCGACATTTGTCCCACTGGAATA
GGATTGTTCTGGGAGCCATTCTTCTGCCTCCAAGTTTGGGACGATCAGGCTGTATGGAATGTGGACTA
CTGCCAGATCCTCAAGGACATTACAGTTGAGGACATGAATGAAATGGAAGGCATTTCTGGAGCTTCTT
CAGTTTAATATTAATGTTCTGCCAGTGTTATGCCAATACTACTTTGACCTTCGCTCCTTAGCAGATG
ACAACAACCTGAATTTCTATTTGCTCCTTAGCAAAGAAAGAGCACAGAACCTAGAGGCTATTTCTAG
ATTGTGTGAAGACAAAGACTTGTGTAGAGCCGCTATGAGAAGGTCTTTCAGTGCTGATAACTTCATTGGT
ATTCAGCGCTCTAAAGCCATCCTCTCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence:

>RG211818 representing NM_152523

Red=Cloning site Green=Tags(s)

MPEDLALESNPSDHPRASTIFLSKSQTDVREKRKSNHLNHVSPGQLTKKYSSCSTIFLDDSTVSQPNLRT
TVKCVTLAIYYHIKNRDANRSLDIFDERSHPLTREKVPEEYFKHDPPEHKFIYRFVRTLSAAQLTAECAI
VTLVYLERLLTYAEIDICPTNWKRIVLGAILLASKVWDDQAVWNVDCQILKDITVEDMNEMERHFLELL
QFNINVPASVYAKYYFDLRLADDNNLNFAPLSKERAQNLEAISRLCEDKDL CRAAMRRSF SADNFIG
IQRSKAILS



TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



ACCN:

NM_152523

ORF Size:

867 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:


This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_152523.2</u>
RefSeq Size:	1953 bp
RefSeq ORF:	870 bp
Locus ID:	151195
UniProt ID:	<u>Q8N7R7</u>
Cytogenetics:	2q33.3
Domains:	CYCLIN, cyclin

Product images:



Circular map for RG211818