

Product datasheet for **RG211883**

ACSL6 (NM_001009185) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ACSL6 (NM_001009185) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ACSL6
Synonyms:	ACS2; FACL6; LACS2; LACS5; LACS 6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG211883 representing NM_001009185
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCTGACCTTCTCCTCGTGTGCGGGGGCTCCCTCTGGCTATTGCTAGAGTTTGTCTCTCACTTCTGG
AGAAGATGCAGACACAGGAGATCCTGAGGATACTGCGACTGCCTGAGCTAGGTGACTTGGGACAGTTTTT
CCGCAGCCTCTCGGCCACCACCCTCGTGTAGTATGGGTGCCCTGGCTGCCATCCTTGCCTACTGGTCACT
CACCGGCCAAAGGCCCTTGCAGCCGCCATGCAACCTCCTGATGCAGTCAGAAGAAGTAGAGGACAGTGGCG
GGGCACGGCGATCTGTGATTGGGTCTGGCCCTCAGCTACTTACCCACTACTATGATGATGCCCGGACCAT
GTACCAGGTGTTCCGCCGTGGGCTTAGCATCTCAGGGAATGGGCCCTGTCTTGGTTTCAGGAAGCCTAAG
CAGCCTTACCAGTGGCTGTCTACCAGGAGGTGGCCGACAGGGCTGAATTTCTGGGGTCCGGACTTCTCC
AGCACAATTGTAAGCATGCACTGATCAGTTTATTGGTGTTCACAAAATCGGCCAGAGTGGATCAT
TGTGGAGCTGGCCTGCTACACATATTCCATGGTGGTGGTCCCCTCTATGACACCCTGGGCCCTGGGGCT
ATCCGCTACATCATCAATACAGCGGACATCAGCACCGTATTGTGGACAAACCTCAGAAGGCTGTGCTTC
TGCTAGAGCATGTGGAGAGGAAGGAGACTCCAGGCCTCAAGCTGATCATCCTCATGGACCCATTCGAAGA
AGCCCTGAAAGAGAGAGGGCAGAAGTGCAGGGTGGTCAATTAAGTCCATGCAGGCCGTGGAGGACTGTGGC
CAAGAGAAATCACCAGGCTCCTGTGCCCCCGCAGCCTGATGACCTTCCATTGTGTGTTTCACAAGCGGCA
CGACAGGGAACCCAAAAGGTGCGATGCTCACCCATGGGAACGTGGTGGCTGATTTCTCAGGCTTTCTGAA
AGTGACAGAGAAAGTGATCTTCCGAGACAGGACGATGTGCTCATCTCCTTCCCTGGCTCTGGCTCACATG
TTTGAGAGAGTAATCCAGTCTGTGCTATTGCCACGGAGGGCGTGTGGCTTCTCCAGGGAGATATCC
GCCTTCTCTCAGATGACATGAAGGCTCTATGCCCCACCCTTCCCTGTGGTCCCACGACTGCTGAACCG
GATGTACGACAAGATCTTCAGCCAGGCAACACACCATTAAAGCGCTGGCTCCTGGAGTTTGCAGCAAAG
CGTAAGCAAGCCGAGGTCCGGAGTGAATCATCAGGAATGATAGTATCTGGGATGAACTCTTCTTAATA
AGATTCAGGCCAGTCTTGGTGGGTGTGTGCGGATGATTGTTACTGGAGCAGCCCCAGCATACCAACAGT
TCTGGGATTTCTCCGGCAGCTTAGGGTGCCAGGTTTATGAAGTTATGGCCAACTGAGTGCACAGCT
GGATGTACCTTACCACCTCCTGGCAGTGGACCTCAGGGCACGTAGGGCGCCACTTCCCTGCAATCATA
TCAAGCTCGTTGATGTTGAGGAAGTGAAGTACTGGGCTGCAAAGGAGAGGGAGAGATATGTGTGAGAGG
ACCAAATGTGTTCAAAGGCTACTTGAAGATCCAGACAGGACGAAGGAGGCCCTGGACAGCGATGGCTGG
CTTCACACTGGAGACATCGGAAAATGGCTGCCCGCAGGAACTCTTAAAATTATTGATCGGAAAAAGCATA
TATTTAAACTTGCTCAGGGAGAATATGTTGCACCCGAGAAGATTGAGAACATCTACATCCGGAGCCAACC
TGTGGCGCAAATCTATGTCCATGGGACAGCTTAAAGGCCTTTTTGGTAGGCATTGTTGTGCTGACCCCT
GAAGTTATGCCCTCCTGGGCCAGAAAGAGAGGAATTGAAGGAACATATGCAGATCTCTGCACAAATAAGG
ATCTGAAGAAAGCCATTTTGAAGATATGGTGAGGTTAGGAAAAGAAAGTGGACTCCATTCTTTGAGCA
GGTTAAAGCCATTCACATCCATTCTGACATGTTCTCAGTTCAAAATGGCTTGCTGACACCAACACTAAAA
GCTAAGAGACCTGAGCTGAGAGAGTACTTCAAAAAACAAATAGAAGAGCTTTACTCAATCTCCATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG211883 representing NM_001009185
 Red=Cloning site Green=Tags(s)

```

MLTFFLVSGGSLWLFVEFVLSLLEKMQTQEILRILRRLPELGDLGQFFRSLSATTLVSMGALAAILAYWFT
HRPKALQPPCNLLMQSEEVEDSGGARRSVIGSGPQLLTHYYDDARTMYQVFRRLSISGNGPCLGFRKPK
QPYQWLSYQEVADRAEFLGSGLLQHNCKACTDQFIGVFAQNRPEWIIVELACYTYSMVVPLYDTLGPGA
IRYIINTADISTVIVDKPQKAVLLLEHVERKETPGLKLIILMDPFEEALKERGQKCGVVIKSMQAVEDCG
QENHQAPVPPQDDLSIVCFTSGTTGNPKGAMLTHGNVVADFSGFLKVKTEKVIFFRQDDVLI SFLPLAHM
FERVIQSVVYCHGGRVGFQGDIRLLSDDMKALCPTIFPVVPRLLNRMYDKIFSQANTPLKRWLLEFAAK
RKQAEVRSIIRNDIWDDELFFNKIQASLGGCVRMIVTGAAPASPTVLGFLRAALGCQVYEGYQTECTA
GCTFTTPGDWTSGHVGAPLPCNHIKLVDVEELNYWACKGEGEICVRGNVFKGYLKDPDRTEALDSDGW
LHTGDIGKWLPA GTLKIIDRKKHIFKLAQGEYVAPEKIENIYIRSQPVAQIYVHGDLSKAFLVGIVVDPD
EVMPSWAQKRGIEGTYADLCTNKDLKKAILEDMVR LGKESGLHSFEQVKAIHIHSDMFSVQNGLLTPTLK
AKRPELREYFKQIEELYSISM
  
```

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_001009185

ORF Size: 2166 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001009185.3](#)

RefSeq Size: 3047 bp

RefSeq ORF: 2169 bp

Locus ID:	23305
UniProt ID:	Q9UKU0
Cytogenetics:	5q31.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Adipocytokine signaling pathway, Fatty acid metabolism, Metabolic pathways, PPAR signaling pathway
Gene Summary:	<p>The protein encoded by this gene catalyzes the formation of acyl-CoA from fatty acids, ATP, and CoA, using magnesium as a cofactor. The encoded protein plays a major role in fatty acid metabolism in the brain. Translocations with the ETV6 gene are causes of myelodysplastic syndrome with basophilia, acute myelogenous leukemia with eosinophilia, and acute eosinophilic leukemia. Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Apr 2011]</p>