

## Product datasheet for **RG211982**

### CD32B (FCGR2B) (NM\_004001) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CD32B (FCGR2B) (NM_004001) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CD32B
Synonyms:	CD32; CD32B; FCG2; FCGR2; FCGR2C; FcRII-c; IGFR2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG211982 representing NM_004001 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGAATCCTGTCATTCTTACCTGTCCTTGCCACTGAGAGTGACTGGGCTGACTGCAAGTCCCCCAGC  
CTTGGGGTCATATGCTTCTGTGGACAGCTGTGCTATTCTGGCTCCTGTTGCTGGGACACCTGCAGCTCC  
CCCAAAGGCTGTGCTGAAACTCGAGCCCCAGTGGATCAACGTGCTCCAGGAGGACTCTGTGACTTGACA  
TGCCGGGGGACTCACAGCCCTGAGAGCGACTCCATTCAAGTGGTCCACAATGGGAATCTCATTCCCACCC  
ACACGCAGCCAGCTACAGGTTCAAGGCCAACAAATGACAGCGGGGAGTACACGTGCCAGACTGGCCA  
GACCAGCCTCAGCGACCCTGTGCATCTGACTGTGCTTTCTGAGTGGCTGGTGTCCAGACCCTCACCTG  
GAGTTCAGGAGGGAGAAACCATCGTGCTGAGGTGCCACAGCTGGAAGGACAAGCCTCTGGTCAAGGTCA  
CATTCTCCAGAAATGGAAAATCCAAGAAATTTTCCCGTTCGGATCCCAACTTCTCCATCCCAAGCAAA  
CCACAGTCACAGTGGTGATTACCACTGCACAGGAAACATAGGCTACACGCTGACTCATCCAAGCCTGTG  
ACCATCACTGTCCAAGCTCCCAGCTCTTACCAGTGGGGATCATTGTGGCTGTGGTCACTGGGATTGCTG  
TAGCGGCCATTGTTGCTGCTGTAGTGGCCTTGATCTACTGCAGGAAAAAGCGGATTTCACTCTCCCAGG  
ATACCTGAGTGCAGGAAATGGGAGAGACCCTCCCTGAGAAACCAGCCAATCCCCTAATCCTGATGAG  
GCTGACAAAGTTGGGGCTGAGAACACAATCACCTATTCCTTCTCATGCACCCGGATGCTCTGGAAGAGC  
CTGATGACCAGAACCCTATT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG211982 representing NM\_004001  
Red=Cloning site Green=Tags(s)

MGILSFLPVLATESDWADCKSPQPWGHMLLWTAFLAPVAGTPAAPPKAVLKLEPQWINVLQEDSVTLT  
 CRGTHSPESDSIQWFHNGNLIPHTHTQPSYRFKANNNDSGEYTCQTGQTSLSDPVHLTVLSEWLVLQTPHL  
 EFQEGETIVLRCHSWKDKPLVKVTFVFQNGKSKKFSRSDPNFSIPQANHSHSGDYHCTGNIGYTLYSSKPV  
 TITVQAPSSSPMGIIVAVVTGIAVAIAVAVALIYCRKKRISALPGYPECREMGETLPEKPANPTNPDE  
 ADKVGAENTITYSLLMHPDALEEPDDQNRI

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_004001

**ORF Size:** 930 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_004001.4](#)

**RefSeq Size:** 1633 bp

**RefSeq ORF:** 933 bp

**Locus ID:** 2213

**UniProt ID:** [P31994](#)

**Cytogenetics:** 1q23.3


**Domains:** ig, IGc2, IG

**Protein Families:** ES Cell Differentiation/IPS, Transmembrane

**Protein Pathways:** B cell receptor signaling pathway, Fc gamma R-mediated phagocytosis, Systemic lupus erythematosus

**Gene Summary:** The protein encoded by this gene is a low affinity receptor for the Fc region of immunoglobulin gamma complexes. The encoded protein is involved in the phagocytosis of immune complexes and in the regulation of antibody production by B-cells. Variations in this gene may increase susceptibility to systemic lupus erythematosus (SLE). Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2010]

### Product images:



Circular map for RG211982