

Product datasheet for **RG211984**

SLC26A10 (NM_133489) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	SLC26A10
Synonyms:	FLJ32854
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >RG211984 representing NM_133489
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAGGCTTGATTTAGCATCCTTGATGTCAGCTCCTAAGAGTCTGGGAAGTGCATTTAAGTCTGGAGGT
 TGGACAAGGCCCTCCACAGCACACCTTTCCATCCACTTCTATCCAGGCATGGCTTTTGCTCTCCT
 GGCCTCCGTGCCCGGTGTTTGGACTCTACACTTCTTTCTTCCCGCTCTCATCTACAGCTTGCTAGGT
 ACTGGGAGACACCTGTCCACAGGAACCTTTCGCATACCTCAGCCTCATGACAGGCTCGGCCGTCGAGCGGC
 TGGTGCCGAACCCCTCGTGGGAATCTGAGCGGAATCGAGAAGGAGCAGCTGGACGCTCAACGGGTTGG
 GGTAGCCGCGGCCGCTGGCCTTCGGGAGCGGGCGTTGATGCTGGGATGTTCTGCTGCAGCTCGGCGTC
 TTGTCCACCTTTTGTCCGAGCCTGTGGTCAAGGCGCTGACCAGCGGGCGCGCTGCACGTGCTCTTGT
 CCCAGCTGCCGAGCCTTTGGGTTGTCCCTCCGCGCCAGATCGGCTGCTTCTCTCTTCAAGACGCT
 GGCCTCCTTGCTGACTGCGCTGCCTCGGAGCAGTCCGGCCGAACCTGACCATCTCCGCGCTCAGCTGGCG
 CTGCTCGTGCCGGTCAAGGAATTGAACGTGAGATTCCGAGACCGGTACCCACGCCGATCCCGGGGAAG
 TCGTCTTGGTGCTTCTGGCCTCCGTGCTGCTTACCTCTTCTGTGGACACAAGATACCAAGTCCAGAT
 AGTGGGGCTGTTGCCTGGAGGATTTCCCAACCCCTCCTCCCAACCTGGCTGAGCTGCCAGGATTCTG
 GCTGACTCGCTGCCATTGCACTGGTTAGTTTTCGGGTGCTGCCTCCCTGGCCTCCATCCATGCAGACA
 AGTATAGCTACACTATTGACTCCAACAGGAGTTCTGGCAGATGGTGCCTCCAACCTCATCTCCTCCCT
 CTTCTCTTGTTCCTCAACTCGGCTACGCTGGCCACCACCAATCTACTGGTGGATGCTGGTGGGAAAACA
 CAGCTGGCAGGCTCTTCTCCTGCACAGTGGTCTGTCGGTGTGCTGTGGCTGGGGCCCTCTTTTACT
 ATCTGCCCAAGGCTGTCTGGCTTGATCAACATCTCCAGCATGCGCCAGGTGTTCTGCCAGATGCAGGA
 ACTTCCACAACATGAGCAGATCAGCCGAGTGGACTTTCTCCTCCAGGTCCCGGGGCTCTGCATCCTGAGC
 TATCCAAACCACTGTACTTTGGGACCCGTGGGCGAGTTTCGCTGCAACCTGGAGTGGCACCTGGGGCTCG
 GAGAAGGAGAAAAGGAGACTTCAAAGCCAGATGGCCAATGGTTCAGTTGCTGAGCCTGTGAGGTGGT
 GGTCTAGACTTCAGTGGTGTACCTTTGCAGATGCTGCTGGGGCCAGAGAAGTGGTGCAGGTGAGGGAG
 AGGCTGGCCAGCCGATGTCGAGATGCTAGGATCCGCCTCCTCTGGCTCAGTGAATGCCTTGGTGACAG
 GGACACTGACCCGGGTAGGACTCCTGGACAGGGTGAATCCAGATCAGCTGTTTGTGAGTGTGACGATGC
 AGCTGCTTATGCCCTGGGAGCCTGTTAAGGGGAGTAGCACCAGGAGCGGGAGCCAGGAGGCACTGGGC
 TCGCGCAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG211984 representing NM_133489
 Red=Cloning site Green=Tags(s)

MRLDLASLMSAPKSLGSAFKSWRLDKAPSPQHTFPSTIPGMAFALLASVPPVFGLYTSFFPVLIYSLLG
 TGRHLSTGTFAILSLMTGSAVERLVPEPLVGNLSGIEKEQLDAQRVGVAAAVAFGSGALMLGMFVLQLGV
 LSTFLSEPVVKALTSAGALHVLLSQLPSLLGLSLPRQIGCFSLFKTLASLLTALPRSSPAELTISALSLA
 LLVPVKELNVRFRDRLPTPIPEVVLLASVLCFTSSVDTRYQVQIVGLLPGGFPQPLLPNLAEPRIL
 ADSLPIALVSFAVSASLASIHADKYSYTTIDSNQEFLAHGASNLISLFSFPPNSATLATTNLLVDAGGKT
 QLAGLFSCTVVLVLLWLPFFYYLPKAVLACINISSMRQVFCQMQLPQLWHISRVDLFLQVPLCILS
 YPTPLYFGRGQFRCNLEWHLGLGEGEKETSKPDGPMVAEPPVRVVLDLDFSGVTFADAAGAREVVQVRE
 RLASRCRDARIRLLLAQCNALVQGTLTRVGLLDRTVPDQLFVSVQDAAYALGSLLRGSSTRSGSQEALG
 CGK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_133489

ORF Size: 1689 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_133489.2](#), [NP_597996.2](#)

RefSeq Size: 2446 bp

RefSeq ORF:	1692 bp
Locus ID:	65012
UniProt ID:	Q8NG04
Cytogenetics:	12q13.3
Protein Families:	Transmembrane
Gene Summary:	Chloride/bicarbonate exchanger.[UniProtKB/Swiss-Prot Function]

Product images:


Circular map for RG211984