

Product datasheet for **RG212188**

ZNF98 (NM_001098626) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ZNF98
Synonyms:	F7175; ZNF739
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >RG212188 representing NM_001098626
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGCCAGGACCCCTTGAAGCCTAGAAATGGGAGTGTTGACATTTAGGGATGTGGCCTTAGAATTCTCTC
 TGGAGGAGTGGCAATGCCTGGACACCGCACAGAGAATTTATATAGGAATGTGATGTTAGAGAACTACAG
 AAACCTGGTCTTTGTGGGTATTGCTGCCTCTAAGCCAGACCTGATCACCTGTCTGGAGCAAGGAAAAGAA
 CCTTGGAAATGTGAAGAGACATGAGATGGTAACTGAACCCCACTGTATATTCTTATTTTGGCCAAAGACC
 TTTGGCCAAAGCAGGGCAAAAAAATTTTCCAAAAAGTGATACTGAGAACATATAAAAAATGTGGACG
 TGAAAAATTTACAGTTAAGAAAAATCTGTAAAAGCATGGATGAGTGTAAAGTGACAAAAGAAATGTTACAAT
 GGACTTAACCAAGTGTGACAACTACCCAGAACAAAAATTTCAATATGACAAATATGTGAAAGTCTTTC
 ATAAATTTTCAAATTCAAACAGACATAAGATAGGACATACTGGAAAGAAATCTTCAAGTGTAAAGAATG
 TGAAAAGTCATTTTGCATGCTTTCACACTAGCTCAACATAAAAGAATTCATAGTGGAGAGAAACCTTAC
 AAATGTAAAGAATGTGGGAAAGCCTATAATGAGGCCTCAAACCTTTCTACACATAAAAGAATTCATACTG
 GAAAGAAACCTACAAATGCGAAGAGTGTGAAAAGCCTTTAACCGGCTCTCACACCTTACTACACATAA
 GATAATTCATACTGGAAAGAAACCTACAAATGTGAGGAGTGTGGCAAAGCTTTTAACCAATCTGCAAAC
 CTTACTACACATAAGAGAATTCATACTGGAGAGAAACCTACAAATGTGAAGAATGTGGCAGAGCTTTTA
 GCCAGTCCTCAACCTTACTGCACATAAGATAATTCATGCTGGAGAGAAACCTTACAAATGTGAAGAATG
 TGGCAAAGCTTTTAGCCAGTCTCAACCTTACTACACATAAGATAATTCATACTGGAGAGAAATCTTAC
 AAATGTGAAGAATGTGTAAGGCCTTTAGCCGGTTATCCACCTTACTACACATAAGAGAATTCATTCTG
 GAGAGAAACCTACAAGTGTGAAGAATGTGGCAAAGCTTTTAACCAATCTCAACCTTACTACACATAA
 GAGAATTCATGCTGGAGAGAAATTCACAAATGTGAAGTATGTAGCAAAGCCTTTAGCCGGTTCTCACAC
 CTTACTACACATAAGAGAATTCATACTGGAGAGAAAGCCTACAAATGTGAAGAATGTGGCAAAGCTTTTA
 ACCTATCATCACACTTACTACACATAAGATAATTCATACTGGAGAGAAACCTACAAATGTGAAGAATG
 TGGCAAAGCTTTTAACCAAGTCTCACTCTTTTAAACATAAGGTAATTCATACTGGAGAGAAAGCCTAC
 AAATGTGAAGAATGTGGCAAAGCTTTTAACCAAGTCTCACTACACCTTACTACACATAAGATGATTCTATG
 GAGAGAAACCTACAAGTGTGAAGAATGCGGCAAAGCCTTTAACAACCTCTATTCTTAACAGACATAA
 GATGATTCTACTGGAGAGAACTCTACAAACCTGAAAGTTGTAAACATGCTTGTGACAACATTGCAAAG
 ATTTCCAAATATAAAAGAAATGTGCTGGTGAGAAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG212188 representing NM_001098626
 Red=Cloning site Green=Tags(s)

MPGPLGSEMGVLTFRDVALEFSLEEWQCLDTAQQLYRNVMLENYRNLVFGIAASKPDLITCLEQKKE
 PWNVVRHEMVTEPPVVSYFAQDLWPKQGKKNYFQKVLRTYKKCGRENQLRKYCKSMDECKVHKECYN
 GLNQCLTTTQNKIFQYDKYVKVFHKFSNSNRHKIGHTGKKSFKCKECEKSFCMLSHLAQHKRIHSGEKPY
 KCKEKGKAYNEASNLSTHKRIHTGKKPYKCEECGKAFNRLSHLTTHKIIHTGKKPYKCEECGKAFNQSAN
 LTTHKRIHTGEKPYKCEECGRAFSQSSTLTAHKIIHAGEKPYKCEECGKAFSSTSSTLTHKIIHTGEKPY
 KCEECGKAFSRLSHLTTHKRIHSGEKPYKCEECGKAFQSSTLTHKRIHAGEKPYKCEVCSKAFSRFSH
 LTTHKRIHTGEKPYKCEECGKAFNLSSTLTHKIIHTGEKPYKCEECGKAFNQSSTLSKHKVIHTGEKPY
 KCEECGKAFNQSSTLTHKMIHTGEKPYKCEECGKAFNNSILNRHKMIHTGEKLYKPESCNACDNIAK
 ISKYKRNCAKEK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001098626

ORF Size: 1716 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001098626.2](#)

RefSeq Size: 2260 bp

RefSeq ORF:	1719 bp
Locus ID:	148198
UniProt ID:	A6NK75
Cytogenetics:	19p12
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG212188