

Product datasheet for **RG212403**

PRR3 (NM_025263) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PRR3 (NM_025263) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: PRR3
Synonyms: CAT56
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG212403 representing NM_025263
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCGAAACGAAAGAAGCAGAATCATCACCAGCCACCGACACAGCAGCAGCCCCGCTGCCGAGCGGG
 AAGAGACTGGAGATGAGGAGGATGGGAGTCCCATCGGACCACCCAGCCTTCTGGGCCCTCCCCCATGGC
 CAATGGAAAACCTGGCGACCCTAAGTCAGCTCTTCACAGAGGTCCTCCAGGATCAAGGGGACCCTGATT
 CCACCACTGCTGAGTCTCCACCTCCTCCTTGGGGTAGAGGCCAATTCGGAGAGGGCTTGCCCCAGGT
 CTAGCCCATATGGTCGTGGTTGGTGGGGAGTCAATGCAGAACCTCCTTTCCGGGGCCAGGCCATGGGG
 TCCACCAGGGGAAGCTTTCACAAGGAACAGAGAAACCTCGAAGGCTCAAAGCTGGTCTCTTATCAAG
 AATACCTGCCCGCCAAGGATGACCCCCAGGTTATGGAAGACAAATCCGACCGCCCTGTCTGCCGACATT
 TTGCCAAAAGGGCCACTGTCGATATGAGGACCTCTGTGCCTTCTACCACCCAGGCGTCAATGGACCTCC
 TCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG212403 representing NM_025263
 Red=Cloning site Green=Tags(s)

MPKRKKQNHQPPTQQQPPLPEREETGDEEDGSPIGPPSLLGPPPMANGKPGDPKSALHRGPPGSRGPLI
 PPLL~~SLPPP~~WGRGPIRRGLGPRSSPYGRGWGVNAEPPFP~~GP~~GHGGPTRGSFHKEQRNPRRLKSWSLIK
 NTCPPKDDPQVMEDKSDRPVCRHF~~AKK~~GHCYEDLCAFYHPGVNGPPL

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



Cloning Scheme:



ACCN: NM_025263

ORF Size: 564 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_025263.4](#)

RefSeq Size: 2261 bp


RefSeq ORF: 567 bp

Locus ID: 80742

UniProt ID: [P79522](#)

Cytogenetics: 6p21.33

Product images:



Circular map for RG212403