

Product datasheet for **RG212570**

MSK2 / RSK-B (RPS6KA4) (NM_003942) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MSK2 / RSK-B (RPS6KA4) (NM_003942) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MSK2 / RSK-B
Synonyms:	MSK2; RSK-B; S6K-alpha-4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG212570 representing NM_003942
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGGGACGAGGACGACGATGAGAGCTGCGCCGTGGAGCTGCGGATCACAGAAGCCAACCTGACCGGGC
 ACGAGGAGAAGGTGAGCGTGGAGAACTTCGAGCTGCTCAAGGTGCTGGGCACGGGAGCCTACGGCAAGGT
 GTTCTGGTGCAGGAGGCGGGCGGGCACGACGCGGGGAAGCTGTACGCCATGAAGGTGCTGCGCAAGGCC
 GCGCTGGTGCAGCGGCCAAGACGCAAGAGCACACGCGCACCGAGCGCTCGGTGCTGGAGCTGGTGCGCC
 AGGCGCCCTTCTGGTACGCTGCACTACGCTTCCAGACGGATGCCAAGCTGCACCTCATCCTGGACTA
 TGTGAGCGGGGGAGATGTTACCCACCTCTACCAGCGCAGTACTTCAAGGAGGCTGAGGTGCGCGTG
 TATGGGGTGAGATCGTGTGCCCTGGAACACCTGCACAAGCTCGGCATCATTTACCGAGACCTGAAAC
 TGGAGAATGTGCTGCTGGACTCCGAGGGCCACATTGCTCCTCACGACTTCGGGCTGAGCAAGGAGTTCCT
 GACGGAGGAGAAAGAGCGGACCTTCTCCTTCTGTGGCACCATCGAGTACATGGCCCCGAAATCATCCGT
 AGCAAGACGGGGCATGGCAAGGCTGTGGACTGGTGGAGCCTGGGCATCTTGCTCTTCGAGCTGCTGACGG
 GGGCTCGCCCTTACCCTGGAGGGCGAGAGGAACACGACGGCTGAGGTGTCTCGACGGATCCTGAAGTG
 CTCCTCCCTTCCCCCTCGGATCGGGCCCGTGGCGCAGGACCTGCTGCAGCGGCTGCTTTGTAAGGAT
 CCTAAGAAGCGATTGGGCGCGGGGCCAGGGGGCACAAGAAGTCCGGAACCATCCCTTCTTCCAGGGCC
 TCGATTGGGTGGCTCGGCTGCCAGGAAGATCCAGCCCCATCCGGCCCCAAATCCGCTCAGAGCTGGA
 TGTGGCAACTTTGCGGAGGAATCACTCGGCTGGAGCCTGTCTACTACCCCTGGCAGCCCCCACT
 GGGGACCCCGAATCTTTCAGGGATACTCTTTGTGGCACCCTCATTCTCTTTGACCACAACAACGCGG
 TGATGACCGATGGGCTGGAAGCGCTGGTGTGGAGACCGCCAGGTGGGCAGCGGTGGCCAGGACGCG
 TATGATGCAGGACTCGCCCTTCTCCAGCAGTACGAGCTGGACCTGCGGGAGCCTGCGCTGGGCCAGGGC
 AGCTTTTCTGTGTGTCGCGCTGCCGCCAGCGCCAGAGCGGCCAGGAGTTCGAGTCAAGATCCTCAGTC
 GCAGGCTGGAGGCGAACACGACGCGCAAGTGGCTGCCCTGCGCCTGTGCCAGTACACCCCAACGTGGT
 GAATCTGCACGAGGTGCATCACGACCAGCTGCACACGTACCTGGTCTGGAGCTGCTGCGGGCGGGGAG
 CTGCTGGAGCACATCCGCAAGAAGCGGCACTTCAGCGAGTCGGAAGCAAGCCAGATCCTGCGCAGCCTCG
 TGTGCGCCGTGAGCTTATGCACGAGGAGGCGGGCGTGGTGCACCGGACCTCAAGCCGAGAACATCCT
 GTACGCCGACGACACGCCCGGGCCCGGTGAAAATCATCGACTTCGGGTTGCGCGGTTGCGGCCGAG
 AGTCCCGGGTGGCCATGCAGACGCCCTGCTTACGCTGCAGTACGCTGCCCCGAGCTGCTGGCGCAGC
 AGGGCTACGACGAGTCTGCGACCTCTGGAGCCTGGGCGTCATTCTGTACATGATGCTGTGCGGGCAGGT
 CCCCTTCCAGGGGCTCTGGCCAGGGCGGGCAGAGCCAGGCGGCCGAGATCATGTGCAAAATCCGCGAG
 GGGCGCTTCTCCCTTGACGGGGAGGCTGGCAGGGTGTATCCGAGGAAGCCAAGGAGCTGGTCCGAGGGC
 TCCTGACCGTGGACCCCGCAAGCGGCTGAAGCTCGAGGGACTGCGGGCAGCTCGTGGCTGCAGGACGG
 CAGCGCGCTCCTCGCCCCGCTCCGGACGCCGACGTGCTCGAGTCTCTGGGCCGAGTGCCTCG
 GGTCTCAACGCCACCTTATGGCATTCAACCGGGCAAGCGGGAGGGCTTCTTCTGAAGAGCGTGGAGA
 ATGCACCCCTGGCAAGCGGCGGAAGCAGAAGCTGCGGAGCGCCACCGCTCCCGCGGGGCTCCCTGC
 ACCAGCCAACCGGGCCGAGCCCCGCTCGCTCCAAGGGGCCCCCGCCGAGCCAACGGCCCCCTGCC
 CCCTCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG212570 representing NM_003942
Red=Cloning site Green=Tags(s)

```

MGDEDDDESCAVELRITEANLTGHEEKVSVENFELLKVLGTGAYGKVFLVRKAGGHDAGKLYAMKVLRKA
ALVQRAKTQEHTRTERSIVLELVRQAPFLVTLHYAFQTDAKLHLILDYVSGGEMFTHLYQRQYFKEAEVRV
YGGEIVLALEHLHKLGIYRDLKLENVLLDSEGHIVLTDGFLSKEFLTEEKERTSFSCGTIEYMAPEIIR
SKTGHGKAVDWWSLGILLFELLTGASPFLEGERNTQAEVSRRIKLCSPFPFPRIGPVAQDLLQRLLCKD
PKKRLGAGPQGAQEVNRHPFFQGLDWWALAARKIPAPFRPQIRSELDVGNFAEEFTRLEPVYSPPGSPPP
GDPRIFFQGYSFVAPSILFDHNNVMTDGL EAPGAGDRPGRAAVARSAMMQDSPFFQYELDLREPALGQG
SFSVCRRCRQRQSGQEFVAVKILSRRLEANTQREVAALRLCQSHPNVNLHEVHHDQLHTYLVLELLRGGE
LLEHIRKKRHFSESEASQILRSLVSAVSMHEEAGVVHRDLKPENILYADDTPGAPVKIIDFGFARLRPQ
SPGVPMQTPCFTLQYAAPELLAQQGYDESCDLWSLGVILYMMLSGQVPFQGASGQGGQSQAAEIMCKIRE
GRFSLDGEAWQGVSEEAKEKELVRGLLTVDPKRLKLEGLRGSSWLQDGSARSSPPLRTPDVLSSGPVARS
GLNATFMAFNRGKREGFFLKSVENAPLAKRRKQKLRSATASRRGSPAPANPGRAPVASKGAPRRANGPLP
PS
    
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_003942


ORF Size: 2316 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_003942.3
RefSeq Size:	3149 bp
RefSeq ORF:	2319 bp
Locus ID:	8986
UniProt ID:	O75676
Cytogenetics:	11q13.1
Domains:	pkinese, S_TK_X
Protein Families:	Druggable Genome, Protein Kinase, Transcription Factors
Protein Pathways:	MAPK signaling pathway, Neurotrophin signaling pathway
Gene Summary:	This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains 2 non-identical kinase catalytic domains and phosphorylates various substrates, including CREB1 and ATF1. The encoded protein can also phosphorylate histone H3 to regulate certain inflammatory genes. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2016]

Product images:



Circular map for RG212570