

## Product datasheet for **RG212743**

### **GH2 (NM\_022558) Human Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	GH2 (NM_022558) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	GH2
Synonyms:	GH-V; GHB2; GHL; GHV; hGH-V
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG212743 representing NM_022558 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCTGCAGGCTCCCGACGTCCTGCTCCTGGCTTTTGGCCTGCTCTGCCTGCTCCTGGCTTCAAGAGG  
GCAGTGCCTTCCCAACCATTCCTTATCCAGGCTTTTGGACAACGCTATGCTCCGCGCCCGTCGCCTGTA  
CCAGCTGGCATAATGACACCTATCAGGAGTTTGAAGAAGCCTATATCCTGAAGGAGCAGAAGTATTCATTC  
CTGCAGAACCCAGACCTCCCTCTGCTTCTCAGAGTCTATTCCAACACCTTCCAACAGGGTGAAAACGC  
AGCAGAAATCTAACCTAGAGCTGCTCCGCATCTCCCTGCTGCTCATCCAGTCATGGCTGGAGCCCGTGCA  
GCTCCTCAGGAGCGTCTTCGCCAACAGCCTGGTGTATGGCGCCTCGGACAGCAACGCTCTATCGCCACCTG  
AAGGACCTAGAGGAAGGCATCCAAACGCTGATAGGCTGGAAGATGGCAGCCCCGGACTGGGCAGATCTT  
CAATCAGTCTACAGCAAGTTTGCACAAAATCGCACAACGATGACGCACTGCTCAAGAACTACGGGCTG  
CTCTACTGCTTCAGGAAGGACATGGACAAGGTCGAGACATTCTGCGCATCGTGCAGTGCCGCTCTGTGG  
AGGGCAGCTGTGGCTTCTAGCTGCCCGGGTGGCATCCCTGTGACCCCTCCCAGTGCCCTCTCTGGTCTGT  
GGAAGTGCTACTCCAGTGCCACCAGCCTTGCC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG212743 representing NM\_022558  
Red=Cloning site Green=Tags(s)

MAAGSRTSLLAFGLLCLSWLQEGSAFPTIPLSRLFDNAMLRARLYQLAYDTYQEFEEAYILKEQKYSF  
 LQNPQTSLFCSESIPTPSNRVKTQQKSNLELLRISLLLQSWLEPVQLLRVAFANSLVYGASDSNVYRHL  
 KDLEEGIQTLIGWKMAAPGLGRSSISPTASLTQNRRTMTHCSRTTGCSTASGRWTWTRSRHSCASCSAALW  
 RAAVASSCPGGIPVTPPQCLSWSWKVLLQCPPALS

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**



**ACCN:** NM\_022558

**ORF Size:** 735 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_022558.4</a>
<b>RefSeq Size:</b>	917 bp
<b>RefSeq ORF:</b>	738 bp
<b>Locus ID:</b>	2689
<b>UniProt ID:</b>	<a href="#">P01242</a>
<b>Cytogenetics:</b>	17q23.3
<b>Protein Families:</b>	Druggable Genome, Secreted Protein
<b>Protein Pathways:</b>	Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway, Neuroactive ligand-receptor interaction
<b>Gene Summary:</b>	<p>The protein encoded by this gene is a member of the somatotropin/prolactin family of hormones which play an important role in growth control. The gene, along with four other related genes, is located at the growth hormone locus on chromosome 17 where they are interspersed in the same transcriptional orientation; an arrangement which is thought to have evolved by a series of gene duplications. The five genes share a remarkably high degree of sequence identity. Alternative splicing generates additional isoforms of each of the five growth hormones, leading to further diversity and potential for specialization. As in the case of its pituitary counterpart, growth hormone 1, the predominant isoform of this particular family member shows similar somatogenic activity, with reduced lactogenic activity. Mutations in this gene lead to placental growth hormone/lactogen deficiency. [provided by RefSeq, Jul 2008]</p>