

Product datasheet for **RG212808**

MTH1 (NUDT1) (NM_198953) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: MTH1 (NUDT1) (NM_198953) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: NUDT1
Synonyms: MTH1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG212808 representing NM_198953
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCGCCTCCAGGCTCTATACCCTGGTCTGGTCTGCAGCCTCAGCGAGTTCTCTGGGCATGAAAA
AGCGAGGCTTCGGGGCCGCGGTGGAATGGCTTTGGGGCAAAGTGCAAGAAGGAGAGACCATCGAGGA
TGGGGCTAGGAGGGAGCTGCAGGAGGAGCGGTCTGACAGTGGACGCCCTGCACAAGGTGGCCAGATC
GTGTTTGAGTTCGTGGGCGAGCCTGAGCTCATGGACGTGCATGTCTTCTGCACAGACAGCATCCAGGGGA
CCCCGTGGAGAGCGACGAAATGCGCCCATGCTGGTTCAGCTGGATCAGATCCCCTTCAAGGACATGTG
GCCCGACGACAGCTACTGGTTTCACTCCTGCTTCAGAAGAAGAAATCCACGGGTACTTCAAGTCCAG
GGTCAGGACACCATCCTGGACTACACTCCGCGAGGTGGACACGGTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG212808 representing NM_198953
Red=Cloning site Green=Tags(s)
MGASRLYTLVVLVLPQRVLLGMKKRFGAGRWNGFGGKVQEGETIEDGARRELQEESSLTVDALHKVQGI
VFEFVGEPELMDVHVFTDSIQGTPVESDEMPCWFQLDQIPFKDMWPDDSYWFPLLLQKKKFKHGYFKFQ
GQDTILDYTLREVDTV

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_198953

ORF Size: 468 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_198953.1, NP_945191.1</u>
RefSeq Size:	728 bp
RefSeq ORF:	471 bp
Locus ID:	4521
UniProt ID:	<u>P36639</u>
Cytogenetics:	7p22.3
Protein Families:	Stem cell - Pluripotency
Gene Summary:	Misincorporation of oxidized nucleoside triphosphates into DNA/RNA during replication and transcription can cause mutations that may result in carcinogenesis or neurodegeneration. The protein encoded by this gene is an enzyme that hydrolyzes oxidized purine nucleoside triphosphates, such as 8-oxo-dGTP, 8-oxo-dATP, 2-hydroxy-dATP, and 2-hydroxy rATP, to monophosphates, thereby preventing misincorporation. The encoded protein is localized mainly in the cytoplasm, with some in the mitochondria, suggesting that it is involved in the sanitization of nucleotide pools both for nuclear and mitochondrial genomes. Several alternatively spliced transcript variants, some of which encode distinct isoforms, have been identified. Additional variants have been observed, but their full-length natures have not been determined. A rare single-nucleotide polymorphism that results in the production of an additional, longer isoform (p26) has been described. [provided by RefSeq, Dec 2018]