

## Product datasheet for **RG213037**

### Macrophage inflammatory protein 5 (CCL15) (NM\_032965) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Macrophage inflammatory protein 5 (CCL15) (NM\_032965) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** CCL15  
**Synonyms:** HCC-2; HMRP-2B; LKN-1; LKN1; MIP-1 delta; MIP-1D; MIP-5; MRP-2B; NCC-3; NCC3; SCYA15; SCYL3; SY15  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG213037 representing NM\_032965  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCCGGATCGCC

ATGAAGGTCTCCGTGGCTGCCCTCCTGCCTCATGCTTGTTGCTGTCTTGGATCCCAGGCCAGTTCA  
TAAATGATGCAGAGACAGAGTTAATGATGTCAAAGCTTCCACTGGAAAATCCAGTAGTTCTGAACAGCTT  
TCACTTTGCTGCTGACTGCTGCACCTCCTACATCTCACAAAGCATCCCGTGTTCACTCATGAAAAGTTAT  
TTTGAACGAGCAGCGAGTGCTCCAAGCCAGGTGCATATTCCTACCAAGAAGGGCGGCAAGTCTGTG  
CCAAACCCAGTGGTCCGGGAGTTCAGGATTGCATGAAAAAGCTGAAGCCCTACTCAATA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG213037 representing NM\_032965  
Red=Cloning site Green=Tags(s)  
MKVSVAAALSCLMLVAVLGSQAQFINDAETELMMSKLPLENPVVLSNFHFAADCCTSYISQSIPLMSKSY  
FETSSECSKPGVIFLTKKGRQVCAKPSGPGVQDCMKKLPYSI

TRTRPLE - GFP Tag - V

**Restriction Sites:** Sgfl-Mlul




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_032965

ORF Size: 339 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_032965.4</a>
<b>RefSeq Size:</b>	955 bp
<b>RefSeq ORF:</b>	342 bp
<b>Locus ID:</b>	6359
<b>UniProt ID:</b>	<a href="#">Q16663</a>
<b>Cytogenetics:</b>	17q12
<b>Domains:</b>	IL8
<b>Protein Families:</b>	Druggable Genome, Secreted Protein, Transmembrane
<b>Protein Pathways:</b>	Chemokine signaling pathway, Cytokine-cytokine receptor interaction
<b>Gene Summary:</b>	This gene is located in a cluster of similar genes in the same region of chromosome 17. These genes encode CC cytokines, which are secreted proteins characterized by two adjacent cysteines. The product of this gene is chemotactic for T cells and monocytes, and acts through C-C chemokine receptor type 1 (CCR1). The proprotein is further processed into numerous smaller functional peptides. Naturally-occurring readthrough transcripts occur from this gene into the downstream gene, CCL14 (chemokine (C-C motif) ligand 14). [provided by RefSeq, Jan 2013]