

## Product datasheet for **RG213051**

### Tcl1 (TCL1A) (NM\_001098725) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tcl1 (TCL1A) (NM\_001098725) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Tcl1  
**Synonyms:** TCL1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG213051 representing NM\_001098725  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCCGAGTGCCCGACACTCGGGGAGGCAGTCACCGACCACCCGACCGCCTGTGGCCTGGGAGAAGT  
 TCGTGTATTTGGACGAGAAGCAGCACGCCTGGCTGCCCTTAACCATCGAGATAAAGGATAGGTTACAGTT  
 ACGGGTGCTCTTGCCTCGGGAAGACGTCGTCCTGGGGAGGCCTATGACCCCAACCCAGATAGGCCCAAGC  
 CTGCTGCCTATCATGTGCAGCTCTACCCTGATGGACGATACCGATCCTCAGACTCCAGTTTCTGGCGCT  
 TAGTGTACCACATCAAGATTGACGGCGTGGAGGACATGCTTCTCGAGCTGCTGCCAGATGAC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG213051 representing NM\_001098725  
 Red=Cloning site Green=Tags(s)  
 MAECPTLGEAVTDHPDRLWAWKEFVYLDEKQHWLPLTIEIKDRLQLRVLLRREDVVLGRPMTPTQIGPS  
 LLPIMWQLYPDGRYRSSDSSFWRLVYHIKIDGVEDMLELLPDD


**TRTRPLE** - GFP Tag - V

**Restriction Sites:** SgfI-MluI



[View online »](#)

## Cloning Scheme:



ACCN: NM\_001098725

ORF Size: 342 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM\\_001098725.2](#)

RefSeq Size: 1405 bp

RefSeq ORF: 345 bp

Locus ID:	8115
UniProt ID:	<a href="#">P56279</a>
Cytogenetics:	14q32.13
Protein Families:	Druggable Genome
Gene Summary:	Overexpression of the TCL1 gene in humans has been implicated in the development of mature T cell leukemia, in which chromosomal rearrangements bring the TCL1 gene in close proximity to the T-cell antigen receptor (TCR)-alpha (MIM 186880) or TCR-beta (MIM 186930) regulatory elements (summarized by Virgilio et al., 1998 [PubMed 9520462]). In normal T cells TCL1 is expressed in CD4-/CD8- cells, but not in cells at later stages of differentiation. TCL1 functions as a coactivator of the cell survival kinase AKT (MIM 164730) (Laine et al., 2000 [PubMed 10983986]).[supplied by OMIM, Jul 2010]

## Product images:



Circular map for RG213051