

## Product datasheet for **RG213082**

### Transglutaminase 7 (TGM7) (NM\_052955) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Transglutaminase 7 (TGM7) (NM_052955) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TGM7
Synonyms:	TGMZ
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG213082 representing NM\_052955  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGATCAGGTGGCAACCTTGC GGCTTGTAGTCTGTGACCTGCAGAGCTCCAGGAACAACAAGGAGCACC  
 ACACGCAGGAGATGGGCGTCAAGCGGCTCACTGTGCGCCGCGCCAGCCCTTCTACCTCCGGCTGAGCTT  
 CAGCCGACCCCTCCAGTCCCAGAACGACCACATCACCTTTGTGGCTGAGACCGGACCCAAAGCCGTACAG  
 CTGCTGGGGACCCGAGCCACATTCTTCTCACCCGGTCCAGCCCGGAATGTCTGGAGCGCTTCTGATT  
 TCACCATGACTCCAACCTCTCTCAAGTTTCCCTTTTACACCAGCCAATGCAGTTATTGGCCATTACAC  
 TCTGAAAATAGAGATCTCTCAGGGCCAAGGTACAGTGTGACTTACCCGCTGGGAACTTTCATCCTACTT  
 TTTAACCCCTTGGAGTCCAGAGGACGACGTCTACCTGCCAAGTAAATACTGCTGCAGGAGTATATCATGC  
 GAGATTATGGCTTTGTTACAAGGGTCATGAAAGATTCATCACCTCTGGCCCTGGAACTACGGGCAGTT  
 TGAAGAGGACATCATAGACATCTGCTTTGAGATCTGAACAAGAGCCTGTATCACTTAAAGAACCCGGCC  
 AAAGACTGTTCCAGCGGAACGACGTGGTGTATGTGTGCAGGGTGGTGTGATCAACAGCAACG  
 ATGACAATGGCGTGTGTCAGGGGAACTGGGGCAGGACTACTCCAAGGGGTGAGTCCCTCTGGAGTGAA  
 GGGCAGTGTGGCCATCTACAGCAGTGGTACGCCAGGGGCGGGCAGCCTGTGAAGTACGGACAGTGTGG  
 GTCTTCGCCTCTGTTATGTGCACCGTAATGAGATGCTTAGGTGTTCCAACCCGTGTTGTTTCCAATTTCC  
 GTTCCGCGCACACGTGGATAGGAACTTGACATCGATACGACTATGACCGAAATGCCGAGATGCTGTC  
 AACTCAGAAACGAGACAAAATA TGGAACTTCCACGTCTGGAATGAGTGTGGATGATCCGAAAGATCTC  
 CCACCAGGATACAACGGGTGGCAGGTTCTGGACCCACTCCCCAGCAGACCAGCAGTGGGCTGTTCTGCT  
 GTGGCCCTGCCTCTGTGAAGGCCATCAGGGAAGGGGATGTCCACCTGGCCTATGACACCCCTTTTGTGTA  
 TGCCGAGGTGAACGCCGATGAAGTCATTTGGCTCCTTGGGGATGGCCAGGCCAGGAAATCTGGCCAC  
 AACACCAGTTCATCGGGAAGGAGATCAGCACTAAGATGGTGGGGTACAGCCAGCGCCAGAGCATACCA  
 GCTCTACAAGTACCCAGAAGGATCCCCTGAGGAGAGAGCTGTCTTATGAAGGCTTCTCGGAAAATGCT  
 GGGCCCCAAAGAGCTTCTTTGCCCTTCTGGATCTCCTGGAGTCTGGGGTCTTAGGGATCAGCCAGCG  
 CAGCTGCAGTTCACCTGGCCAGGATACCCGAGTGGGGCCAGGACCTGCAGCTGCTGCTGCGTATCCAGA  
 GGGTGCCAGACAGCACCCACCCTCGGGGGCCATCGGACTGGTGGTGCCTTCTGTGCACAGGCCCTGCT  
 GCATGGGGTGGTACCCAGAAGCCCTTCTGGAGGCACACAGTGGGATGAACCTGGACTTTGGGAAGGAG  
 ACACAGTGGCCGCTCCTCCTGCCCTACAGCAATTACAGAAACAAGCTAACGGACGAAAAGCTCATCCGCG  
 TGTCTGGCATCGCGGAGTTGAAGAGACAGGGAGGTCCATGCTGGTCTAAAAGATATCTGTCTGGAGCC  
 TCCCCACTTGCTATTGAGGTGTCTGAGAGGGCTGAGGTGGGCAAGGCGCTGAGAGTCCATGTACCCCTC  
 ACCAACACCTTAATGGTGGCTCTGAGCAGCTGCACGATGGTGTGGAAGGAAGCGGCCTCATCAATGGGC  
 AGATAGCAAAGGACCTTGGGACTCTGGTGGCCGGACACACCCTCCAAATCAACTGGACCTTACCCGAC  
 CAAAGCTGGACCCCGCAGCTCCAGTTCTCATCAGCAGCAACGAGGTCAAGGAGATCAAAGGCTACAAG  
 GACATATTCGTCACTGTGGCTGGGGTCCC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG213082 representing NM\_052955  
 Red=Cloning site Green=Tags(s)

```

MDQVATLRLESVDLQSSRNKEHHTQEMGVKRLTVRRGQPFYLRLSFSRPFQSQNDHITFVAETGPKPSE
LLGTRATFFLTRVQPGNVWSASDFTIDSNSLQVSLFTPANAVIGHYTLKIEISQGQGHSVTYPLGTFILL
FNPWSPEDDVYLPSEILLQEYIMRDYGFVYKGERHERFITSWPWNYGQFEEDIIDICFEILNKSLYHLKNPA
KDCSQRNDVVYVCRVVSAMINSNDDNGVLQGNWGEDYSKGVSPLEWKGSVAILQQWSARGGQPVKYGQCW
VFASVMCTVMRCLGVPTRVVSNFRSAHNVDRNLTIIDTYDRNAEMLSTQKRDKIWNFHVWNECWMIKRD
PPGYNGWQVLDPTPQQTSSGLFCCGPASVKAIREGDVHLAYDTPFVYAEVNADEVIWLLGDGQAQEILAH
NTSSIGKEISTKMVGSDQRQSITSSYKYPEGSPEERAVFMKASRMLGPPQRASLPFLDLLESGLRDQPA
QLQLHLARIPWQDLQLLLRIQRVPDSTHPRGPIGLVVRFCQALLHGGGTQKPFWRHTVRMNLDFGKE
TQWPLLLPYSNYRNKLTDEKLIRVSGIAEVEETGRSMLVLKDICLEPPHLSIEVSERAEVGKALRVHVT
TNTLMVALSSCTMVLESGLINGQIAKDLGTLVAGHTLQIQDLDYPTKAGPRQLQVLISSNEVKEIKGYK
DIFVTVAGAP
  
```


TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**


**ACCN:** NM\_052955

**ORF Size:** 2130 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_052955.1](#), [NP\\_443187.1](#)

**RefSeq Size:** 2312 bp

**RefSeq ORF:** 2133 bp

**Locus ID:** 116179

**UniProt ID:** [Q96PF1](#)

**Cytogenetics:** 15q15.2-q15.3

**Protein Families:** Druggable Genome

**Gene Summary:** Transglutaminases (TGM; EC 2.3.2.13) are a family of structurally and functionally related enzymes that stabilize protein assemblies through the formation of gamma-glutamyl-epsilon lysine crosslinks. For additional background information on transglutaminases, see TGM1 (MIM 190195).[supplied by OMIM, Jul 2002]