

Product datasheet for **RG213252**

AGAP8 (AGAP4) (NM_133446) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AGAP8 (AGAP4) (NM_133446) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	AGAP4
Synonyms:	AGAP-4; AGAP-8; AGAP8; CTGLF1; CTGLF5; MRIP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG213252 representing NM_133446
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGGATCGCC

ATGGGGAACATACTGACCTGTCGTGTGCACCCTAGCGTCAGCCTCGAGTTTGACCAGCAGCAGGGGTCGG
TGTGTCCCTCTGAATCTGAGATCTATGAGGCAGGAGCTGGGGACAGGATGGCAGGAGCGCCCATGGCTGC
TGCTGTACAGCCTGCTGAGGTGACTGTTGAAGTTGGTGAGGACCTCCACATGCACCACGTTCTGTGACCGG
GAGATGCCTGAAGCTTTGGAGTTAACCCTTCTGCCAATCCAGAGGCAAGCACAATATCCAGAGGAACT
CTCAAACAGATGTTGTAGAAATAAGAAGAAGCAACTGTACAAACCATGTATCTACTGTGCGTTTCAGTCA
ACAATACAGCTTGTGTTGACAATATTCCTTGATGACAGCACAGCCATCCAGCATTATCTTACAATGACA
ATAATATCTGTGACCTTGGAGATACCTCATCATATCACACAAAGAGATGCAGATAGATCTTTGAGCATA
CTGATGAACAGTTACACTCATTTGCAGTTCCACCGTGACATTATGAAGAAAAGAAATGGAGGTGGGAG
TTTAAATAACTATTCCTCTCCATTCCATCGACTCCCAGCACCAGCCAGGAGACCCTCAGTTCAGTGT
CCTCCCACTGCCAACACACCACCCCGTTTGCAAGCGGTCCATGCGCTGGTCCAACCTGTTTACATCTG
AGAAAGGGAGTGACCCAGACAAAGAGAGAAAGCCCGGAGAATCATGCTGACACCATCGGGAGCGGCAG
AGCCATCCCCATTAACAGGGCATGCTCTTAAAGCGAAGTGGAAATGGCTGAAGACATGGAAAAAGAAA
TACGTCACCCTGTGTTCCAATGGCGTGCTCACCTATTATTCAAGCTTAGGTGATTATATGAAGAATATTC
ATAAAAAAGAGATTGACCTTCAGACATCTACCATCAAAGTCCCAGGAAAGTGGCCATCCCTAGCCACATC
GGCCTGCACACCCATCTCCAGCTCTAAAAGCAATGGCCTATCCAAGGACATGGACACCGGGCTGGGTGAC
TCCATATGCTTCAGCCCCAGTATCTCCAGCACCACCAGCCCCAAGCTCAACCCGCCCTCTCCTCATG
CTAATAAAAAAGAACACCTAAAGAAGAAAAGCACCAACAACCTTATGATTGTGTCTGCCACTGGCCAAAC
GTGGCACTTTGAAGCCACGACGTATGAGGAGCGGGATGCCTGGTCCAAGCCATCCAGAGCCAGATCCTG
GCCAGCCTGCAGTCATGCGAGAGCAGTAAAAGCAAGTCCCAGCTGACCAGCCAGGCAAGGCCATGGCCC
TGCAGTCGATCCAAAACATGCGTGGGAACGCCACTGTGTGGACTGTGAGACCCAGAATCTAAGTGGGC
CAGTTTGAACCTGGGAGTCTCATGTGTATTGAATGCTCAGGTATCCACCGCAGTCTTGGACCCGCCTT
TCCCGTGTGCGATCTCTGGAGCTGGATGACTGGCCAGTTGAGCTCAGGAAGGTTATGTCATCTATTGGCA
ATGACCTAGCCAACAGCATCTGGGAAGGGAGCAGCCAGGGGCAGACAAAACCTCAGAAAAGTCCACGAG
GGAAGAGAAGGAACGGTGGATCCGTTCCAATATGAGGAGAAGCTTTTCTGGCCCCACTACCCTGCACT
GAGCTGTCCCTGGGCCAGCAGCTGCTGCGGGCCACCGCTGATGAGGACCTGCAGACAGCCATCCTGTGC
TGGCAGATGGCTCCCGTGAGGAGGTGAACGAGACCTGTGGGGAGGGAGACGGCTGCACGGCGCTCCATCT
GGCCTGCCGCAAGGGGAATGTGGTCTGGCGCAGCTCCTGATCTGGTACGGGGTGGACGTATGGCCCGA
GATGCCACGGGAACACAGCGCTGACCTACGCCGGCAGGCCCTCCAGCCAGGAGTGCATCAACGTGCTTC
TGCAGTACGGCTGCCCGACGAGTGCCTG

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG213252 representing NM_133446
 Red=Cloning site Green=Tags(s)


MGNILTCRVHPSVSLFDQQQGSVCPSESEIYEAGAGDRMAGAPMAAAVQPAEVTVEVGEDLHMHHVDR
 EMPEALEFNPSANPEASTIFQRNSQTDVVEIRRSNCTNHVSTVRFSSQYSLCSTIFLDDSTAIQHYLTMT
 IISVTLIPHHITQRDADRSLIPDEQLHSFAVSTVHIMKKRNGGSLNNYSSSIPSTPSTSQEDPQFSV
 PPTANTPTPVCKRSMRWSNLFSEKGSDDPKERKAPENHADTIGSGRAIPIKQGMMLLKRSGKWLKTWKKK
 YVTLCSNGLVLYYSSLGDMYKNIHKKEIDLQTSTIKVPGKWPLATSACTPISSSKSNGLSKDMDTGLGD
 SICFSPSISSTTSPKLNPPSPHANKKKHLKKKSTNNFMIVSATGQTWHFEEATTYEERDAWVQAIQSQIL
 ASLQSCSSKSKSQLTSQSKAMALQSIQNMGRNAHCVDCEQNPKWASLNLGVLMCIECSIHRSLGTRL
 SRVRSLELDDWPVELRKVMSSIGNDLANSIWEGSSQGQTKPSEKSTREEKERWIRSKYEEKLFLAPLPCT
 ELSLQQLLRATADEDLQTAIILLAHGSREEVNETCGEGDGCTALHLACRKGNNVLAQLLIWYGVDMAR
 DAHGNTALTYARQASSQECINVLLQYGCPECV

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_133446

ORF Size: 1989 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_133446.1](#), [NP_597703.1](#)

RefSeq Size: 2420 bp
RefSeq ORF: 1992 bp
Locus ID: 119016
UniProt ID: [Q96P64](#)
Cytogenetics: 10q11.22
Domains: ArfGap, PH, ANK
Gene Summary: Putative GTPase-activating protein.[UniProtKB/Swiss-Prot Function]