

Product datasheet for **RG213347**

C9orf66 (NM_152569) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	C9orf66
Synonyms:	RP11-5906.1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG213347 representing NM_152569
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAGACACTCCGTGGCGCGCCCAACACGGTTACCCAGGCGCCTATCACCTTCTGGGACCCGGCCACCT
GTAAGAACTGGAGGGCGGAGCGGGCGAGGTTGTTTCGGGAAGGGATCCCGCGGGCTGAGGACATCCAG
GAGCACCGAGATTCTTGGGGAGGGCTCGCCGGGCCAGCGCAGGAGCCGCGCGTCCAGCGGCTCCA
CCGCCGAGCCTCGGGAACCGGGCGCCCAAGTCTCAGGCGCGCCCGCGCGGACGCGAATGGACTCGT
CCGGCCTCGGCCCTGCAGCGAAGCTCCGCTTACACCAGCGCTGGGCTGTCCGGCCGCAACTGCGCGC
CGCGGGGGCGTCTTACCTGTTGATCTTGAGCGCGAACGCGCGGCGCTCTGCGCTCGGCAGAGTGCCAT
GGCCCCCGGCGGTTGGTGGCTTCTAGGGTCGCGGGTGCGGAAAGCGGGGGCTGGCGCGCCGAGAG
TAGCCGCGGAACATGCCAGCCAAGCGCAAACCTCGTCTGTCGGTCGGCGCTGGAACCTCCGCATTTC
ACCTTCTAAGCCGAAATCACCCCGGGGAGGGTCCGAGCTCGGAGCAGCGACGCGCGGCTCCGCCACCCG
GCCTGGCGCGCAGGACGCGGGCGAGGAGGGCGTGGGCCAGCGCGGAAGTGGCCTCCCGATACTGGG
GCAGACGGCGCGCTCTGCCGGGGCCGCCACCTCCGGCCGAAGGAGCGCGGGCGGACGACAGGCGACC
ACTGCGCGCAGGCCGAAGCTCCATCTGCCGAAGTCCGAGGTCGCCAGGGAATGTAGGCAAGTCGGG
GAGCCCCACAAAGCTGGCGAGGTCGGAATCACCCAGGGATTCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTAA



Protein Sequence: >RG213347 representing NM_152569
 Red=Cloning site Green=Tags(s)

MRHSVARPTRLPRRLSPFWDPATCKNLEGGAGEVVRGRDPRLRTSRSTEILGEGLAGPSAGAAARPAAP
 PPQPREPGAPGLRRAPPRTRMDSSGLGPCSEAPLHTSAGLSGRNLRAAGVLPVDLERERAAALCARQSGH
 GPPAVRWLLGSRGAE SGLARRVAAEHAQPSANLVCRSALETSAFPSPKPSRGRVRRSSDGLRHP
 AWRAGSGRGRGPSAELASRYWGRRRALPGAADLRPKGARADRRRLRAGRKLHLPEAARLPGNVGKSG
 EPHKAGEVGNHPRDS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_152569

ORF Size: 885 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_152569.1](#), [NP_689782.1](#)


RefSeq Size: 2279 bp

RefSeq ORF: 888 bp

Locus ID: 157983

Cytogenetics: 9p24.3

Product images:



Circular map for RG213347