

Product datasheet for **RG213638**

UQCRQ (NM_014402) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: UQCRQ (NM_014402) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: UQCRQ
Synonyms: MC3DN4; QCR8; QP-C; QPC; UQCR7
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG213638 representing NM_014402
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCCGCGAGTTTGGGAATCTGACGCGGATGCGGCATGTGATCAGCTACAGCTTGTCACCGTTTCGAGC
AGCGCGCCTATCCGCACGTCTTCACTAAAGGAATCCCAATGTTCTGCGCCGATTTCGGGAGTCTTTCTT
TCGCGTGGTGCCGAGTTGTAGTGTATCTTATCTACACATGGGGACTGAAGAGTTCGAGAGATCC
AAGAGGAGGATCCAGCTGCCTATGAAAATGACAAATGAGCAACGCATCCGGATGACGGTTCCTGTCTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG213638 representing NM_014402
Red=Cloning site Green=Tags(s)

MGREFGNLTRMRHVISYSLSPFEQRAYPHVFTKGIPNVLRRIRESFFRVVPQFVVFYLIYTWGTEEFERS
KRRIQLPMKMTNEQRIRMTVPCL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_014402

ORF Size: 279 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_014402.2</u> , <u>NP_055217.1</u>
RefSeq Size:	388 bp
RefSeq ORF:	249 bp
Locus ID:	27089
UniProt ID:	<u>O14949</u>
Cytogenetics:	5q31.1
Domains:	UcrQ
Protein Pathways:	Alzheimer's disease, Cardiac muscle contraction, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease
Gene Summary:	This gene encodes a ubiquinone-binding protein of low molecular mass. This protein is a small core-associated protein and a subunit of ubiquinol-cytochrome c reductase complex III, which is part of the mitochondrial respiratory chain. [provided by RefSeq, Jul 2008]