


Product datasheet for **RG213640**

Teneurin 4 (TENM4) (NM_001098816) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Teneurin 4 (TENM4) (NM_001098816) Human Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | Teneurin 4 |
| Synonyms: | Doc4; ETM5; ODZ4; ten-4; Ten-M4; TEN4; TNM4 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| Restriction Sites: | Sgfl-RsrII |
| Cloning Scheme: | |

Cloning sites used for ORF Shutting:



Kozac
Consensus

EcoR I BamH I Kpn I RBS Sgf I Asc I

CTATAGGCGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCCGCGATCGCCGGCCGCCAGATCT

Hind III Nhe I Rsr II Mlu I Not I Xho I GFP Tag

CAAGCTTAAGCTAGCTAGCGGACCG ACG CGT ACG CGG CCG CTC GAG ATG GAG AGC GAC ---

T R T R P L E M E S D - - -

Pme I Fse I

--- GAA GAA AGA GTT TAA ACGGCCGGCCGGAGCT

- - - E E R V Stop

| | |
|-----------|--------------|
| ACCN: | NM_001098816 |
| ORF Size: | 8307 bp |



[View online >](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001098816.3](#)

RefSeq Size: 13566 bp

RefSeq ORF: 8310 bp

Locus ID: 26011


UniProt ID: [Q6N022](#)

Cytogenetics: 11q14.1

Protein Families: Druggable Genome, Transmembrane

Gene Summary: The protein encoded by this gene plays a role in establishing proper neuronal connectivity during development. Defects in this gene have been associated with hereditary essential tremor-5. [provided by RefSeq, Oct 2016]

Product images:



Circular map for RG213640