

Product datasheet for **RG213681**

Semaphorin 3A (SEMA3A) (NM_006080) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Semaphorin 3A (SEMA3A) (NM_006080) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Semaphorin 3A
Synonyms:	coll-1; COLL1; HH16; Hsema-I; Hsema-III; SEMA1; SEMAD; SEMAIII; SEMAL; SemD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG213681 representing NM_006080
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGCTGGTAACTAGGATTGTCTGTCTTTTCTGGGAGTATTACTTACAGCAAGACAACTATCAGA
 ATGGGAAGAACAATGTGCCAAGGCTGAAATTATCTTACAAAGAAATGTTGGAATCCAACAATGTGATCAC
 TTTCAATGGCTTGGCCAACAGCTCCAGTTATCATACCTTCCTTTTGGATGAGGAACGGAGTAGGCTGTAT
 GTTGGAGCAAAGGATCACATATTTTCATTCGACCTGGTTAATATCAAGGATTTTCAAAGATTGTGTGGC
 CAGTATCTTACACCAGAAGAGATGAATGCAAGTGGGCTGGAAAAGACATCCTGAAAGAATGTGCTAATTT
 CATCAAGTACTTAAGGCATATAATCAGACTCACTTGTACGCCTGTGGAACGGGGCTTTTCATCCAATT
 TGCACCTACATTGAAATGGACATCATCCTGAGGACAATATTTTTAAGCTGGAGAACTCACATTTTGAAA
 ACGGCCGTGGGAAGAGTCCATATGACCCTAAGCTGCTGACAGCATCCCTTTTAATAGATGGAGAATTATA
 CTCTGAACTGCAGCTGATTTTATGGGGCGAGACTTTGCTATCTCCGAACCTTTGGGCACCACCACCCA
 ATCAGGACAGAGCAGCATGATTCAGGTTGGCTCAATGATCCAAAGTTCATTAGTGCCACCTCATCTCAG
 AGAGTGACAATCCTGAAGATGACAAAGTATACTTTTTCTCCGTGAAAATGCAATAGATGGAGAACACTC
 TGGAAAAGCTACTCACGTAGAATAGGTCAGATATGCAAGAATGACTTTGGAGGGCACAGAAGTCTGGTG
 AATAAATGGACAACATTCCTCAAAGCTCGTCTGATTTGCTCAGTGCCAGGTCCAAATGGCATTGACACTC
 ATTTTGATGAACTGCAGGATGTATTCTAATGAACTTTAAAGATCCTAAAAATCCAGTTGTATATGGAGT
 GTTTACGACTTCCAGTAACATTTTCAAGGGATCAGCCGTGTGTATGTATAGCATGAGTGTGTGAGAAGG
 GTGTTCCCTGGTCCATATGCCACAGGGATGGACCAACTATCAATGGGTGCCTTATCAAGGAAGAGTCC
 CCTATCCACGGCCAGGAACCTTGTCCAGCAAAACATTTGGTGGTTTTGACTCTACAAAGGACCTTCTGGA
 TGATGTTATAACCTTTGCAAGAAGTCATCCAGCCATGTACAATCCAGTGTTTCTATGAACAATCGCCCA
 ATAGTGATCAAAAACGGATGTAATTTATCAATTTACACAAATTGTCGTAGACCGAGTGGATGCAGAAGATG
 GACAGTATGATGTTATGTTTATCGGAACAGATGTTGGGACCGTCTTAAAGTAGTTTCAATTCTAAGGA
 GACTTGGTATGATTTAGAAGAGGTTCTGCTGGAAGAAATGACAGTTTTTTCGGGAACCGACTGCTATTTCA
 GCAATGGAGCTTTCCACTAAGCAGCAACAATAATATTGGTTCAACGGCTGGGTTGCCAGCTCCCTT
 TACACCGGTGTGATATTTACGGGAAAGCGTGTGCTGAGTGTGCCTCGCCGAGACCCCTTACTGTGCTTG
 GGATGGTTCTGCATGTTCTCGCTATTTTCCACTGCAAAGAGACGCACAAGACGACAAGATATAAGAAAT
 GGAGACCCACTGACTCACTGTTTCAAGCTTACACCATGATAATCACCATGGCCACAGCCCTGAAGAGAGAA
 TCATCTATGGTGTAGAGAATAGTAGCACATTTTGGAAATGCAGTCCGAAGTCGCAGAGAGCGCTGGTCTA
 TTGGCAATTCAGAGGCGAAATGAAGAGCGAAAAGAAGAGATCAGAGTGGATGATCATATCATCAGGACA
 GATCAAGGCCCTTCTGCTACGTAGTCTACAACAGAAGGATTCAGGCAATTACCTCTGCCATGCGGTGGAAC
 ATGGGTTCCATACAACTCTTCTTAAGGTAACCTGGAAGTCATTGACACAGAGCATTTGGAAGAACTTCT
 TCATAAAGATGATGATGGAGATGGCTCTAAGACCAAGAAATGTCCAATAGCATGACACCTAGCCAGAAG
 GTCTGGTACAGAGACTTCATGCAGCTCATCAACCACCCCAATCTCAACACGATGGATGAGTTCTGTGAAC
 AAGTTTGGAAAAGGGACCGAAAACAACGTCGGCAAAGGCCAGGACATACCCAGGGAACAGTAACAAATG
 GAAGCACTTACAAGAAAATAAGAAAGGTAGAAAACAGGAGGACCCACGAATTTGAGAGGGCACCCAGGAGT
 GTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG213681 representing NM_006080
 Red=Cloning site Green=Tags(s)


```

MGWLTRIVCLFWGVLLTARANYQNGKNNVPRCLKLSYKEMLESNNVITFNGLANSSSYHTFLLDEERSRLY
VGAKDHIFSDLVNIKDFQKIVWPVSYTRRDECKWAGKDILKECANFIKVLKAYNQTHLYACGTGAFHPI
CTYIEIGHHPEDNIFKLENSHFENGRGKSPYDPKLLTASLLIDGELYSGTAAADFMRDFAIFRTLGHHP
IRTEQHDSRWLNDPKFISAHLISESDNPEDDKVYFFFRENAIDGEGHSGKATHARIGQICKNDFGGHRSLV
NKWTTFLKARLICSVPGPNGIDTHFDELQDVFLMNFKDPKNPVVYGVFTTSSNIFKGSAVCMYSMSDVRR
VFLGPYHRDGPYQWVPYQGRVYPRPGTGPCSKTFGGFDSTKDLRDDVITFARSHPAMYNPVFPMMNRP
IVIKTDVNYQFTQIVVDRVDAEDGQYDVMFIGTDVGTVLKVVSIKETYDLEEVLLLEEMTVFREPTAIS
AMELSTKQQQLYIGSTAGVAQLPLHRCDIYGKACAECCLARDPYCAWDGSACSRYPYPTAKRRTRRQDIRN
GDPLTHCSDLHHDNHGHSPEERIIYGVENSSTFLECSPKSQRALVYWQFQRRNEERKEEIRVDDHIIRT
DQGLLLRSLQKQDSGNYLCHAVEHGF IQTLLKVTLEVIDTEHLEELLHKDDDGDSKTKEMSNSMTPSQK
VWYRDFMQLINHPLNTMDEFCEQVWKRDRKQRRQRPGHTPGNSNKWKHLQENKGRNRRTHEFERAPRS
V
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_006080


ORF Size: 2313 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_006080.1</u> , <u>NP_006071.1</u>
RefSeq Size:	2530 bp
RefSeq ORF:	2316 bp
Locus ID:	10371
UniProt ID:	<u>Q14563</u>
Cytogenetics:	7q21.11
Domains:	Sema, ig, IG, PSI
Protein Families:	Secreted Protein
Protein Pathways:	Axon guidance
Gene Summary:	<p>This gene is a member of the semaphorin family and encodes a protein with an Ig-like C2-type (immunoglobulin-like) domain, a PSI domain and a Sema domain. This secreted protein can function as either a chemorepulsive agent, inhibiting axonal outgrowth, or as a chemoattractive agent, stimulating the growth of apical dendrites. In both cases, the protein is vital for normal neuronal pattern development. Increased expression of this protein is associated with schizophrenia and is seen in a variety of human tumor cell lines. Also, aberrant release of this protein is associated with the progression of Alzheimer's disease. [provided by RefSeq, Jul 2008]</p>

Product images:



Circular map for RG213681