

Product datasheet for **RG213692**

SAR1B (NM_001033503) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SAR1B (NM_001033503) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: SAR1B
Synonyms: ANDD; CMRD; GTBPB; SARA2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG213692 representing NM_001033503
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCCTTCATATTTGATTGGATTTACAGTGGTTTCAGCAGTGTGCTACAGTTTTAGGATTATATAAGA
AAACTGGTAACTGGTATTTCTGGATTGGATAATGCAGGAAAAACAACATTGCTACACATGCTAAAAGA
TGACAGACTTGGACAACATGTCCAACATTACATCCCCTCCGAAGAACTGACCATTGCTGGCATGACG
TTTACAACTTTTGATCTGGTGGACATGTTCAAGCTCGAAGAGTGTGGAAAACTACCTTCTGCTATCA
ATGGCATTGTATTTCTGGTGGATTGTGCAGACCACGAAAGGCTGTAGAGTCAAAGAAGAACTTGATTC
ACTAATGACAGATGAAACCATTGCTAATGTGCCTATACTGATTCTTGGGAATAAGATCGACAGACCTGAA
GCCATCAGTGAAGAGAGGTTGCGAGAGATGTTTGGTTTATATGGTCAGACAACAGGAAAGGGGAGTATAT
CTCTGAAAGAAGTGAATGCCCGACCCTTAGAAGTTTTTCATGTGTAGTGTCTCAAAGACAAGGTTACGG
AGAAGGCTTCGCTGGATGGCACAGTACATTGAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG213692 representing NM_001033503
Red=Cloning site Green=Tags(s)

MSFIFDWIYSGFSSVLQFLGLYKKTGKLVFLGLDNAGKTTLLHMLKDDRLGQHVPTLHPTSEELTIAGMT
FTTFDLGGHVQARRVWKNYLPAINGI VFLVDCADHERLLESKEELDSLMTDETIANVPILILGNKIDRPE
AISEERLREMFGLYGQTTGKGSISLKELNARPLEVFMCSVLKRQGYGEGFRWMAQYID

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001033503

ORF Size: 594 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033503.3](#)

RefSeq Size: 1377 bp

RefSeq ORF: 597 bp

Locus ID: 51128


UniProt ID: [Q9Y6B6](#)

Cytogenetics: 5q31.1

Gene Summary:

The protein encoded by this gene is a small GTPase that acts as a homodimer. The encoded protein is activated by the guanine nucleotide exchange factor PREB and is involved in protein transport from the endoplasmic reticulum to the Golgi. This protein is part of the COPII coat complex. Defects in this gene are a cause of chylomicron retention disease (CMRD), also known as Anderson disease (ANDD). Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Mar 2010]

Product images:



Circular map for RG213692