

Product datasheet for **RG213742**

CTRP3 (C1QTNF3) (NM_181435) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CTRP3 (C1QTNF3) (NM_181435) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CTRP3
Synonyms:	C1ATNF3; CORCS; CORS; CORS-26; CORS26; CTRP3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG213742 representing NM_181435 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTTTGGAGGCAGCTCATCTATTGGCAACTGCTGGCTTTGTTTTCTCCCTTTTGCCTGTGTCAAG
ATGAATACATGGAGGTGAGCGGAAGAATAATAAGTGGTGGCAAGAATAGTGCAAAGCCACCAGCAGAC
TGGCCGTAGCGGCTCCAGGAGGGAGAAAGTGAAGAGCGGAGCCATCCTAAAAGTGGACTGTGGATAAT
AACACTTCTACAGACCTAAAATCCCTGAGACCAGATGAGCTACCGCACCCGAGGTAGATGACCTAGCCC
AGATCACCATTTCTGGGGCCAGTCTCCACAAACCGGAGGACTACCCCAAGACTGCAGTAAGTGTGTCA
TGGAGACTACAGCTTTGAGGCTACCAAGGCCCCCTGGGCCACCGGGCCCTCCTGGCATTCCAGGAAAC
CATGGAACAATGGCAACAATGGAGCCACTGGTCATGAAGGAGCCAAAGGTGAGAAGGGCGACAAAGGTG
ACCTGGGGCTCGAGGGGAGCGGGGCAGCATGGCCCCAAAGGAGAGAAGGGCTACCCGGGGATTCCACC
AGAAGTTCAGATTGCATTCATGGCTTCTGGCAACCACTTCAGCAATCAGAACAGTGGGATTATCTTC
AGCAGTGTGAGACCAACATTGGAACTTCTTTGATGTCATGACTGGTAGATTTGGGGCCCCAGTATCAG
GTGTGATTTCTCACCTTCAGCATGATGAAGCATGAGGATGTTGAGGAAGTGTATGTGTACCTTATGCA
CAATGGCAACACAGTCTTCAGCATGTACAGCTATGAAATGAAGGGCAAATCAGATACATCCAGCAATCAT
GCTGTGCTGAAGCTAGCCAAAGGGGATGAGGTTGGCTGCGAATGGGCAATGGCGCTCTCCATGGGGACC
ACCAACGCTTCTCCACCTTTCAGGATTCTGCTCTTTGAACTAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG213742 representing NM_181435
Red=Cloning site Green=Tags(s)

MLWRQLIYWQLLALFFLPFLCQDEYMEVSGRTNKVVARIVQSHQQTGRSGSRREKVRERSHPKTGTVDN
 NTSTDCLKSLRPDELPHPEVDDLAQITTFWQSPQTGGLPPDCSKCCHGDYSFRGYQPPGPPGPPGIPGN
 HGNNGNNGATGHEGAKGEKGDGDLGPRGERGQHGPKEKGYPGIPPELQIAFMASLATHFSNQNSGIIF
 SSVETNIGNFVDVMTGRFGAPVSGVYFFTF SMMKHEDVEEVYVYL MHNNGTVF SMYSYEMK GKSDTSSNH
 AVLKLAGDEVWLRMGNGALHGDHQRFFSTFAGFLLFETK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_181435

ORF Size: 957 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_181435.6](#)

RefSeq Size: 3789 bp

RefSeq ORF: 960 bp


Locus ID: 114899

UniProt ID: [Q9BXJ4](#)

Cytogenetics: 5p13.2

Protein Families: Secreted Protein

Product images:



Circular map for RG213742