

Product datasheet for **RG213912**

Clorf141 (NM_001013674) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	Clorf141
Synonyms:	MGC48454
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG213912 representing NM_001013674
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAGAAAAATCCTAGAGAAGTTGGATGTCCTTGATAAGCAAGCAGAGATAATCTTGCCAGAAGAA
CAAAGATAAACAGGCTTCAGAGTGAAGGAAGAAAAACAATATGGCTATACCCCTGACATTTGATTTTCA
GTTGGAATTTGAAGAAGCTCTTGCTACATCCGCGTCTAAGGCAATATCAAAGATCAAAGAAGACAAGTCA
TGCAGCATTACAAAATCAAAAATGCATGTCTTTCAAATGTGAGCCTGAACCTAGAAAGAGTAATTTTG
AAAAGTCAAAATTTAAGACCATTCTTTATTCAAACAAATGTAAAAATAAAGAAAGTGAGTCAACAGCTCA
AATTGAAAAAACCTAGGAAACCATTGGATTCTGTTGGTCTCTTAGAAGGTGATAGAAATAAAGAAAA
AAATCTCCACAGATGAACGATTTTAAATATAAAGAAAAACAATCGGTGAGAAATTATCAATTAAGTAAGT
ATAGGTCAGTAAGAAAGAAAAGCTTGCTCCGTTGTGCTTTGAGGATGAATTGAAAAATCCACATGCCAA
GATAATCAACGTTAGTCCAACAAAGACAGTAACCTCTCACATGGAACAAAAGGACACAAATCCCATTAATT
TTCCATGACACAGAATATGTACGAATGTTACTTTTGACAAAAAATAGATTTTCTCTCATCCTTTGGAAA
ATGAAACATTTACCCACATAAAGAACAAATTTCAATTTAGAAAGAAATTGTGAAATCCTCAAAATCTAT
AATTGGCAATCAATCTATTTCTCTTTCAAACCCCAAAAACTATGCCTACAGTACAGAGAAAAGATATA
CAGATCCCTATGTCTTTTAAAGCGGGCCACAACTGTAGATGATAAACTAAAGAAAGAAAACATAAAGC
AGACACTAGAAAACAGATCTTGAATACACTCTATAATTTCTCACAGAATTTTCTAGCCTAACAAAACA
ATTTGTGGGTTACCTTGATAAAGCTGTTATTATGAAATGAGTGCCCAAACTGGAAAATTTGAAAGAATG
TTTTCTGCAGGAAAACACGAGCATACCCACATCCAGTGCCTTACCTGTCAAATGTTACTCAAAGCCTT
TTAAATATATATGAATAAATAATGTAACACCACTGGATAATTTGTTAACTTATCAAATGAAATTTT
AAATGCCTCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG213912 representing NM_001013674
 Red=Cloning site Green=Tags(s)

MAEKILEKLDVLDKQAEIILARRTKINRLQSEGRKTTMAIPLTFDFQLEFEEALATSASKAISKIKEDKS
 CSITKSKMHVSFKCEPEPRKSNFEKSNLRPFQITNVKNKESESTAQIEKKPRKPLDSVGLLEGDRNKRK
 KSPQMNDNFNIKENKSVRNYQLSKYRSVRKKSLLPLCFEDELKNPHAKIINVSPKTKVTSHMEQKDTNP
 FHDTEYVRMLLLTKNRFSSHPLENENIYPHKRTNFILERNCEILKSIIGNQSISLFPQKTMPTVQRKDI
 QIPMSFKAGHTTVDDKLKKTKNQTLENRSWNTLYNFSQNFSSLTKQFVGYLDKAVIHEMSAQTGKFERM
 FSAGKPTSIPTSSALPVKCYSKPFKYIYELNNVTPLDNLLNLSNEILNAS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001013674

ORF Size: 1200 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001013674.1](#), [NP_001013696.1](#)


RefSeq Size: 2143 bp

RefSeq ORF: 1202 bp

Locus ID: 400757

Cytogenetics: 1p31.3

Product images:



Circular map for RG213912