

Product datasheet for **RG213947**

DHPS (NM_013406) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DHPS (NM_013406) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	DHPS
Synonyms:	DHS; DS; MIG13; NEDSSWI
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG213947 representing NM_013406 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAGGTTCCCTGGAACGGGAGGCGCCAGCGGGGGCGCTGGCCGCCGTGCTAAAGCACAGCTCGACGT
TGCCGCCCGAAAGCACCCAGGTCCGGGGCTACGACTTCAACCGCGGTGTGAATTACCGCGCACTGCTGGA
GGCCTTCGGCACCACCGCTTCCAAGCAACCAACTTCGGGCGCGCTGTACAGCAAGTCAATGCCATGATC
GAGAAGAAGCTGGAACCACTGTACAGGATGAAGACCAGCACGGGACCTGACCCAGAGCCGCCGCCAC
TTACCAGCTGCACATTTTCTGGGATATACATCCAACCTCATCAGTTCAGGCATCCGTGAGACCATTTCG
CTACCTTGTGCAGCACAACATGGTGGACGTATTGGTGACCACAGCTGGCGCGTGGAGGAAGACCTCATC
AAGTGCCTGGCGCCACATACTTGGGCGAGTTTAGCCTCAGGGGAAGGAGCTCCGGGAGAACGGGATCA
ATAGGATCGAAAACCTGCTGGTGCCCAATGAGAATTACTGCAAGTTTGGAGACTGGCTGATGCCATTCT
GGACCAGATGGTGTGGAGCAGAACACAGAGGGTGTAAAGTGGACGCCTTCTAAGATGATCGCCCGGCTG
GGCAAGGAGATCAACAACCCAGAGTCCGTGTATTACTGGGCCAGAAGAACCACATCCCTGTGTTTATGC
CCGCACTTACAGACGGCTCGCTGGGCGACATGATCTTCTCCATTCTACAAGAACCAGGCGCTGGTCTCT
GGACATCGTTGAGGCCAGGAGTTTGTGGCTCTGACTCAGGTGCCCGACCAGACGAGGCTGTCTCTGG
GGCAAGATCCGGTGGATGCACAGCCGTCAAGGTCTATGCTGACGCCTCCCTGGTCTTCCCCTGCTTG
TGCTGAAACCTTGGCCAGAAGATGGATGCCTTCATGCATGAGAAGAACGAGGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG213947 representing NM_013406
 Red=Cloning site Green=Tags(s)

MEGLSLEAPAGALAAVLKHSSTLPPESTQVRGYDFNRGVNRYRALLEAFGTTGFQATNFGRAVQQVNAMI
 EKKLEPLSQDEDEDHADLTQSRRLTSCITFLGYTSNLISSGIRETIRYL VQHNMDVLTVTAGGVEEDLI
 KCLAPTYLGEFSLRGKELRENGINRIGNLLVPNENYCKFEDWLMPIILDQMVMQNTGKWTSPKMIARL
 GKEINNPESEVYYWAQKNHIPPVFPALTDGSLGDMIFFHSYKNPGLVLDIVAEQEFDGSDSGARPDEAVSW
 GKIRVDAQPVKYYADASLVFPLLVAETFAQKMDAFMHEKNED

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_013406

ORF Size: 966 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_013406.2](#)

RefSeq Size: 1210 bp

RefSeq ORF: 969 bp

Locus ID: 1725


UniProt ID: [P49366](#)

Cytogenetics: 19p13.13

Domains: DS

Gene Summary: This gene encodes a protein that is required for the formation of hypusine, a unique amino acid formed by the posttranslational modification of only one protein, eukaryotic translation initiation factor 5A. The encoded protein catalyzes the first step in hypusine formation by transferring the butylamine moiety of spermidine to a specific lysine residue of the eukaryotic translation initiation factor 5A precursor, forming an intermediate deoxyhypusine residue. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, May 2011]

Product images:



Circular map for RG213947