

Product datasheet for **RG213953**

C12orf71 (NM_001080406) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: TurboGFP
Symbol: C12orf71
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG213953 representing NM_001080406
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCATATTCATCCTCTAACAGCGACATAGAGGACGACAGCTCCAATCCAATCCAACCTGAGCCTCT
CTGTGGGCTATTTCCCTGTGAGGACACCCCTGTGAGGACACAACCTCTGGGAAGATGCACCTTCCAA
GGGTCTTCCATCCACTTTCTGCCTCCCGTCCAAGGGGCATGGGGGACTGAAAGGATAGGGAGACGCATG
AAGAGACAAGACCAATTCAGGATGAACCAGAGCAGTTTGTCAAATAAGCATCTTCTGGCCTGGGACG
TGGACATTGGCTCCGATAACACAGACTCAAGAGCTAATAGGCTTCTAAATGGAGACAACCTGTGGATAGA
CAAGTTACCAAAAGAGAGAACAACTGTCTGTTGGCAAATGAATAATCTTGCAAGAGTTTCAGATA
TTTCTAGAAAATCTGAAAGATGATGACGTGTATTTCCTGAACTGCTCAGCAAGATTTCAGCTATCCA
GCGGCTCCCTCCGAAATGGTTCAGATGATAAGCCAGGCAACCGCCAGCCAAAGGACAAGCGCTCCAGA
GATCTCCTCAATCCTGTCAGAGCAGCCAGAGAAGGATGACACTCCTTCACACACAGGCCAGTGCTGC
CTGAACCTTGGTGGGCTTCAGCTGGCTGAGGCAGCGTATCTCCCTCTCTGCTGAGGAGGGATCACC
CTGTGAATGCCACCAAAAGTCCCATCGGTGAGCACCACGAAAGACTCTTTCACAGAGGCAAGAGAAT
TCAACCTCAAGAAACCCTTGAATTAGGACATCCCAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG213953 representing NM_001080406
Red=Cloning site Green=Tags(s)


MAYSSNSDIEDSSKSNLSLVGYFPCEDTPCEDTTSWEDAPSKGPSIHFLPPVQGAWTERIGRRM
KRQDQIQDEPEQFCKLSIFLAWDVDIGSDNTDSRANRLNGDNLWIDKLPKERTKLSVGKLNLLVQEFQI
FLENLKDDDAVPETAQQDFQLSSGSPPMVQMISQATASQRTSAPEISSILSEQPEKDDTPSHTQAQCC
LNFGWAFSWLRQILPSLLRRDHPVNATKSPHRSAPTKRLFHRGKRIQPQETLELGHPI

TRTRPLE - GFP Tag - V



Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001080406

ORF Size: 807 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001080406.1](#), [NP_001073875.1](#)

RefSeq Size: 966 bp


RefSeq ORF: 810 bp

Locus ID: 728858

UniProt ID: A8MTZ7

Cytogenetics: 12p11.23

Product images:



Circular map for RG213953