

Product datasheet for **RG213958**

Stathmin 1 (STMN1) (NM_005563) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Stathmin 1 (STMN1) (NM_005563) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: Stathmin 1
Synonyms: C1orf215; Lag; LAP18; OP18; PP17; PP19; PR22; SMN
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG213958 representing NM_005563
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTTCTTCTGATATCCAGGTGAAAGAACTGGAGAAGCGTGCCTCAGGCCAGGCTTTTGAGCTGATTC
TCAGCCCTCGGTCAAAGAATCTGTTCCAGAATCCCCCTTCCCTCCAAGAAGAAGGATCTTCCCT
GGAGGAAATTCAGAAGAAATTAGAAGCTGCAGAAGAAAGCGCAAGTCCCATGAAGCTGAGGTCTTGAAG
CAGCTGGCTGAGAAACGAGAGCACGAGAAAGAAGTGCTTCAGAAGGCAATAGAAGAGAACAACAATTCA
GTAAAATGGCAGAAGAGAACTGACCCACAAAATGGAAGCTAATAAAGAGAACCAGAGGCACAAATGGC
TGCCAAACTGGAACGTTTGCAGAGAAGGATAAGCACATTGAAGAAGTGCAGGAAGAACAAGAATCCAAA
GACCCTGCTGACGAGACTGAAGCTGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG213958 representing NM_005563
Red=Cloning site Green=Tags(s)

MASSDIQVKELEKRASGQAFELILSPRSKESVPEFPLSPPKKKDLLEEIQKKLEAAEERRKSHEAEVLK
QLAEKREHEKEVLQKAIENNNSKMAEEKLTHKMEANKENREAQMAAKLERLREKDKHIEEVRKNKESK
DPADETEAD

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_005563

ORF Size: 447 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_005563.4](#)

RefSeq Size: 1542 bp


RefSeq ORF: 450 bp

Locus ID: 3925

UniProt ID: [P16949](#)

Cytogenetics:	1p36.11
Domains:	Stathmin
Protein Pathways:	MAPK signaling pathway
Gene Summary:	This gene belongs to the stathmin family of genes. It encodes a ubiquitous cytosolic phosphoprotein proposed to function as an intracellular relay integrating regulatory signals of the cellular environment. The encoded protein is involved in the regulation of the microtubule filament system by destabilizing microtubules. It prevents assembly and promotes disassembly of microtubules. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2009]

Product images:



Circular map for RG213958