

## Product datasheet for **RG214016**

### PRDM13 (NM\_021620) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PRDM13 (NM_021620) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PRDM13
Synonyms:	MU-MB-20.220; PFM10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>RG214016 representing NM\_021620  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCACGGAGCCGCCAGAGCGCCAGCCACCAGCGTGAAGTCCGACTGCTGCATCCCGGCCGGCTTGCGCC  
TCGGACCGGTGCCTGGTACCTTCAAGCTGGCAAGTACCTGTGACACCGCAGGGAGCCCGGGCCTAAGAA  
AAAGGTATTGACCATTACAGCTCTGCCACCAGGTGCGCATGGTGAGAGGGGAGCTGGTGGACGAGTCG  
GGGGGCTCCCCTCTGGAGTGGATAGGGTTAATCCGGGCAGCCAGAACTCCAGGAACAGACTCTGGAAG  
CTATTGCAGACTTACCCGGAGGACAGATCTTCTACCGAGCATTGCGAGACGTCCAGCCAGGGGAGGAGCT  
GACAGTGTGGTATTCTAACTCCTTGCTCAGTGGTTCGACATCCCACCACAGCGACTCCGACTCACGAC  
GAGAAAGGGGAGGAGCGCTACATCTGCTGGTACTGCTGGAGGACGTTTAGATACCCCAACAGCCTTAAGG  
CACACCTGCGTTTCCACTGCGTGTTCAGCGCGGTGGAGGCGGCCCTTCTGCACCACGAACACGCGGC  
TCGCCAAGGCGCGTCCCAGCGGTGATGGCTCGGTCTCTCCCAAACCCCGCGCCGATTTTCGCC  
GGCCTTCCCAGGCAGGAACCTTTCGACCCACCCCTGGGCCCGCCACCAAGTTCAGGCCTGCGGTGCGC  
GGGAGGGCATCAAGCGCGAGGCCTTTCGCGCCCTCGGCCACCTCGCCGACCCAGGCAAGTGGGGCA  
GCCAAGAAGGGCAAGGAGCAGCTGGACCGTGCCTGGACATGAGCGGAGCCGCCGAGGACAAGGGCAC  
TTCTCGGCATCGTGGCGGGTCTCGCGGGGGTTCGGCAGCCTGGCTTCTACCCGGCGTGCCTCAG  
CTTTAAGCCCGCCGCTAGCGAGGGCGCGCGCGCTCACGGCGACCCCTACCGGGAGGAGAGCAG  
CAGCAAGCAAGGAGCGCGCTCGCTTGGGCGAGCTGCTGGGCGGGGCGGGGCTGCGGGCGCCCGGG  
AGCGGGGAGAAGTCCGGCGGGCGCGCGGGTACCACCATACCACCACGCGCACACCACCACCATC  
CCAAGTGCCTGCTCGCTGGGACCCCGCCCGCGCCGCTGGCTGCCCTGCTCTGGGGCCCTGCG  
CGGCTTCCCTCTGCTCTCCGTCCTCCCGAAGAGGCGTCCGCCTTCAAGCACGTGGAGCGCGCCCGCC  
GCAGCCGCGCGCTGCCAGGAGCGCTTATGCGCAGCTGCCCTGCGCCGGGGTTCGCCCTCGAGCGCT  
GCGCGCTGCCGCCCTCGACCCGGCGGTCTCAAAGCCTATCCGGGTGGTGAAGTGCAGCCACCTGCCCGC  
CGTCATGCCGGCTTTACAGTCTACAACGGGAGCTGCTTACGGCTCACCGGCCACCACCGCTTATTAC  
CCGCTCAAATTGCACTTCGGCGGGTCTGAAGTATCCGGAGTCCATCTCTACTTCAGCGGGCTGCAG  
CGGCCGCCCTAAGCCCGCCGAGCTGGGGTCTGCGCCAGCATCGACCGAGAGATCGCCATGCACAATCA  
GCAGCTGTCCGAGATGGTGCCTGGGAAGGGTCCGGACGCTGGACTCGGGGACGTTGCCACCGGCCGTC  
GCGCGCGGGGAGGCACCGGGGCGCGCAGCGAGGCAGCGGCGCAGGTAAGCCCAAGACCGGCCACC  
TGTGCCTTACTGTGGCAAGCTGTACTCGCGCAAGTATGGGCTCAAGATCCACATGCGGACGCACACGGG  
CTACAAGCCACTCAAGTCAAAAGTCTGTCTGCGGCCCTTCGGCGACCCAGCAATCTCAACAAGCACATC  
CGGCTGCACGCCGAGGGCAATACGCCCTACCGCTGCGAGTTCTGCGGCAAGGTAAGTGTGCGCCCGGGG  
ACCTGGAGCGACATGTCAAGTCCCGCCACCCTGGCCAGAGTCTGCTCGCCAAAGCGGGGACGCGCCGGG  
TGCCGAGCCCGGCTATCCCCGGAGCCTGGGGATCCCAAGAGCGACGACGTGACGTGGACGTCTGCTTC  
ACAGACGACCAGAGCGACCCCGAGTTGGGGCGCGGGGAGCGCGACTTG

**ACGCGT**ACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG214016 representing NM\_021620  
 Red=Cloning site Green=Tags(s)


MHGAARAPATSVSADCCIPAGLRLGPVPGTFKLGKYLSDRREPGPKKKVLTIQSAHQVRMVRGELVDES  
 GGSPLIEWIGLIRAARNSQEQTLEAIADLPGGQIFYRALRDVQPGEELTVWYSNSLAQWFDIPTTATPTHD  
 EKGEERYICWYCWRTFRYPNSLKAHLRFHCVFSGGGGGAFLHHEHAARQGAVPAADGLGLSPKPPAPDFA  
 APSQAGTLRPHPLGPPPQACGAREGIKREASSAPSATSPTPGKWGQPKKGEQLDRALDMSGARGQGH  
 FLGIVGGSSAGVGLAFYPGVRSFAFKPAGLARAAAAAHGDPYREESSKQAGLALGRLLGGGRACGRPG  
 SGENSAAGGAGHHHHHHHHHPKCLLAGDPPPPPPGLPCSGALRGFPLLSVPPEEASAFKHVERAPP  
 AAAALPGARYAQLPPAPGLPLERCALPPLDPGGLKAYPGGECSHLPVMPAFTVYNGELLYGSPATTAYY  
 PLKLFHFGLLKYPESISYFSGPAAAALSPAELGSLASIDREIAMHNQQLSEMAAGKGRGRLDSTLPPAV  
 AAAGGTGGGGSGGSGAGPKTGHLCLYCGKLYSRKYGLKIHMRTHTGYKPLKCKVCLRPFGDPSNLNKHI  
 RLHAEGNTPYRCEFCGKVLVRRRDLERHVKSRRHPGQSLAKAGDGPGEAEPGYPPEPDPKSDSDVDVCF  
 TDDQSDPEVGGGGERDL

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



## Plasmid Map:



ACCN: NM\_021620

ORF Size: 3201 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_021620.2](#), [NP\\_067633.1](#)

**RefSeq Size:** 3234 bp

**RefSeq ORF:** 2124 bp

**Locus ID:** 59336

**UniProt ID:** [Q9H4Q3](#)

**Cytogenetics:** 6q16.2

**Gene Summary:** May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]