

Product datasheet for **RG214266**

PI 3 Kinase p85 alpha (PIK3R1) (NM_181504) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PI 3 Kinase p85 alpha (PIK3R1) (NM_181504) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PI 3 Kinase p85 alpha
Synonyms:	AGM7; GRB1; IMD36; p85; p85-ALPHA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG214266 representing NM_181504
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTACAATACTGTTTGGAAATATGGAAGACCTGGATTTAGAATATGCCAAGACAGATATAAATTGTGGCA
 CAGACTTGATGTTTTATAGAAATGGACCCACCAGCACTGCCTCCTAAACCACCAAAACCTACTACTGT
 AGCCAACAACGGTATGAATAACAATATGTCCTTACAAGATGCTGAATGGTACTGGGGAGATATCTCGAGG
 GAAGAAGTGAATGAAAACTTCGAGATACAGCAGACGGGACCTTTTTGGTACGAGATGCGTCTACTAAAA
 TGCATGGTGATTACTCTTACACTAAGGAAAGGGGAAATAACAAATTAATCAAATATTTTCATCGAGA
 TGGGAAATATGGCTTCTGACCCATTAACTTCAGTTCTGTGGTTGAATTAATAAACCCTACCGGAAT
 GAATCTCTAGCTCAGTATAATCCCAAATGGATGTGAAATTAATTTATCCAGTATCCAAATACCAACAGG
 ATCAAGTTGTCAAAGAAGATAATATTGAAGCTGTAGGAAAAAATACATGAATATAACTCAGTTTCA
 AGAAAAAGTCGAGAATATGATAGATTATGAAGAATATACCCGCACATCCAGGAAATCCAAATGAAA
 AGGACAGCTATTGAAGCATTAAATGAAACCATAAAAAATTTTGAAGAACAGTGCCAGACCAAGAGCGGT
 ACAGCAAAGAATACATAGAAAAGTTTAAACGTGAAGGCAATGAGAAAGAAATACAAGGATTATGCATAA
 TTATGATAAGTTGAAGTCTCGAATCAGTGAAATATTGACAGTAGAAGAAGATTGGAAGAAGACTTGAAG
 AAGCAGGCAGCTGAGTATCGAGAAATGACAAACGTATGAACAGCATTAAACCAGACCTTATCCAGCTGA
 GAAAGACGAGAGACCAATACTTGATGTGGTTGACTCAAAAAGGTGTTCCGGCAAAGAAGTTGAACGAGTG
 GTTGGGCAATGAAACACTGAAGACCAATTTCACTGGTGAAGATGATGAAGATTTGCCCATCATGAT
 GAGAAGACATGGAATGTTGGAAGCAGCAACCGAAACAAAGCTGAAAACCTGTTGCGAGGGAAGCGAGATG
 GCACTTTTCTTCCGGGAGAGCAGTAAACAGGGCTGCTATGCCTGCTCTGTAGTGGTGGACGGCGAAGT
 AAAGCATTGTGTATAAAACAAAACAGCAACTGGCTATGGCTTTGCCGAGCCCTATAAATTGTACAGCTCT
 CTGAAAGAAGTGGTGTACATTACCAACACACCTCCCTTGTGCAGCACAACGACTCCCTCAATGTACAC
 TAGCTACCCAGTATATGCACAGCAGAGGCGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG214266 representing NM_181504
 Red=Cloning site Green=Tags(s)

MYNTVWNMEDLDLEYAKTDINCGTDLMFYIEMDPPALPPKPPKPTTVANNGMNNMSLQDAEWYWGDISR
 EEVNEKLRDADGTFLVRDASTKMHGDYTLRLKGGNNKLIKIFHRDGKYGFSDPLTFSSVVELINHYRN
 ESLAQYNPKLDVKLLYPVSKYQQDQVVKEDNIEAVGKKLHEYNTQFQEKREYDRLYEEYTRTSQEIQMK
 RTAIEAFNETIKIFEEQCQTQERYKEYIEKFKREGNEKEIQRIMHNYDKLKSRISEIIDSRRRLEEDLK
 KQAAEYREIDKRMNSIKPDLIQLRKRTRDQYLMWLTQKGVROKKLNEWLGNENTEDQYSLVEDDELPHHD
 EKTWNVGSNRNKAENLLRGRDGTFLVRESSKQGCYACSVVVDGEVKHCVINKTATGYGFAEPYNLYSS
 LKELVLHYQHTSLVQHNSLNVTLAYPVYAQRR

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_181504

ORF Size: 1362 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_181504.4](#)

RefSeq Size: 5663 bp

RefSeq ORF: 1365 bp

Locus ID: 5295

UniProt ID: [P27986](#)


Cytogenetics: 5q13.1

Protein Families: Druggable Genome

Protein Pathways: Acute myeloid leukemia, Apoptosis, B cell receptor signaling pathway, Chemokine signaling pathway, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, ErbB signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, Glioma, Insulin signaling pathway, Jak-STAT signaling pathway, Leukocyte transendothelial migration, Melanoma, mTOR signaling pathway, Natural killer cell mediated cytotoxicity, Neurotrophin signaling pathway, Non-small cell lung cancer, Pancreatic cancer, Pathways in cancer, Phosphatidylinositol signaling system, Progesterone-mediated oocyte maturation, Prostate cancer, Regulation of actin cytoskeleton, Renal cell carcinoma, Small cell lung cancer, T cell receptor signaling pathway, Toll-like receptor signaling pathway, Type II diabetes mellitus, VEGF signaling pathway

Gene Summary: Phosphatidylinositol 3-kinase phosphorylates the inositol ring of phosphatidylinositol at the 3-prime position. The enzyme comprises a 110 kD catalytic subunit and a regulatory subunit of either 85, 55, or 50 kD. This gene encodes the 85 kD regulatory subunit. Phosphatidylinositol 3-kinase plays an important role in the metabolic actions of insulin, and a mutation in this gene has been associated with insulin resistance. Alternative splicing of this gene results in four transcript variants encoding different isoforms. [provided by RefSeq, Jun 2011]

Product images:



Circular map for RG214266