

Product datasheet for **RG214312**

SARG (C1orf116) (NM_001083924) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	SARG
Synonyms:	SARG
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG214312 representing NM_001083924
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCCAAAAAGCCAAGGAAACAGTCTCAACCAGGTACACACAACCCAGCCTCCTCTGCAGGGTTGC
CTCAGAATGCAAGAGCTGAAGATGCTCCCCTCTCATCAGGGGAGGACCCAAACAGCCGACTAGTCCCCT
CACAACCCCTAAGCCCCGGAAGCTGCCACCTAATATTGTTCTGAAGAGCAGCCGAAGCAGTTCCACAGT
GACCCCAGCACTGGCTGTCCCGCCACACTGAGGCTGCCCCGAGATTCTGGCCTGATCTCCTGTTTAC
TGCAAGAGCAGAGAAAAGCACGTAAGAAGCTCTAGAGAAGCTGGGGCTACCCAGGATCAAGATGAGCC
TGGACTCCACTTAAGTAAGCCACAGCTCCATCAGACCCAAGGAGACACGGGCCAGCATCTGTCCCA
GCTCCAGGTCTGGCTCAGCCTGCAGCTCCAGCCAGGCCTCAGCAGCTATTCTGCTGCTGGGAAGGCTC
TGGCTCAAGCTCCGGCTCCAGCTCCAGGTCCAGCTCAGGGACCTTTGCCAATGAAGTCTCCAGCTCCAGG
CAATGTTGAGCTAGCAAACTATGCAATTTCTATCCCTAAGGCCCAAGGGCAAACAGTCCCCTGACT
CCACCGAAGCCAGAGTCAGGGCTGACTCTCCAGGAGAGCAACACCCCTGGCCTGAGACAGATGAACCTCA
AGTCCAACACTCTGGAGCGCTCAGGCGTGGGACTGAGCAGCTACCTTTCAACTGAGAAAGATGCCAGCCC
CAAAACCAGCACTTCTCTGGGAAAGGGCTCCTTCTTGACAAGATCTCGCCAGTGTCTTACGTAATTCT
CGGCCCCCGCCGCTCCTTGGGCACGGGGAAGATTTGCAAGTATCCAGGTAGGCAAGCTGGCTGACC
TGGAGCAGGAGCAGAGCTCCAAGCGCTGTCTACCAAGGACAGAGCCGTGACAAGCTTCTCGCCCCC
CTGTGTGCTCAAGATCTCCCAAGGGTGTCCCAATGAACACAGAAGGGAGGCCCTGAAGAAGCTG
GACTGTTGAAGGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG214312 representing NM_001083924
Red=Cloning site Green=Tags(s)

MSQKAKETVSTRYTQPQPPAGLPQNARAEDAPLSSGEDPNSRLAPLTTPKPRKLPPNIVLKSSRSSFHSDPQHWLSRHTAAPGDSGLISCSLQEQRKARKEALEKLGLPQDQDEPGLHLSKPTSSIRPKETRAQHLSPAPGLAQPAAPAQASAAIPAAGKALAQAPAPAGPAQGGLPMKSPAPGNVAASKMPIISIPKAPRANSALT PPKPESGLTLQESNTPLRQMNFKSNTLERSGVGLSSYLSTEKDASPKTSTSLGKGSFLDKISPSVLRNS RPRPASLGTGKDFAGIQVGLADLEQEQSSKRLSYQQQSRDKLPRPPCVSVKISPKGVPNEHRREALKKL GLLKE

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001083924

ORF Size: 1065 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001083924.2](#)

RefSeq Size: 5316 bp

RefSeq ORF: 1068 bp


Locus ID: 79098

UniProt ID: [Q9BW04](#)

Cytogenetics: 1q32.1

Gene Summary: Putative androgen-specific receptor.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG214312