

Product datasheet for **RG214326**

FAM55D (NXPE4) (NM_017678) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM55D (NXPE4) (NM_017678) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	NXPE4
Synonyms:	C11orf33; FAM55D
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG214326 representing NM_017678 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAAAATTCAATACAATTAGTGTCTCCAAATGCAACAAAGAAACAGTTGCAATGAAAGAGAAATGCA
AGTTTGAATGACATCCACAATCCCCAGTGGGCATGTCTGGAGAAACACATGGAATCCTGTCTCCTGTAG
TTTGGCTACAGTCAAAATGAAGGAATGCCTGAGAGGAAAACATACCTAATGGGAGATTCCACGATC
CGCCAGTGGATGGAATACTTCAAAGCCAGTATCAACACACTGAAGTCAGTGGATCTGCATGAATCTGGAA
AATTGCAACACCAGCTTGTGTGGATTTGGATAGGAACATCAACATCCAGTGGCAAAAATATTGTTATCC
CTTGATAGGATCAATGACCTATTCAGTCAAAGAGATGGAGTACCTCACCCGGGCCATTGACAGAACTGGA
GGAGAAAAAATACTGTCATTGTTATTTCCCTGGGCCAGCATTTTCAGACCCTTTCCATTGATGTTTTTA
TCCGAAGGGCCCTCAATGTCCACAAAGCCATTCAGCATCTTCTTCTGAGAAGCCCAGACACTATGTTAT
CATCAAAACAGAAAACATCAGGGAGATGTACAATGATGCAGAAAGATTTAGTGACTTTTATGTTTACATT
CAATATCTCATATAAAGGACATTTCCAGGATCTCAGTGTGAGTATCATTGATGCCTGGGATATAACAA
TTGCATATGGCACAAATAATGTACCCACCTCAACATGTAGTCGGAAATCAGATTAATATATTATTTAA
CTATATTTGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG214326 representing NM_017678
Red=Cloning site Green=Tags(s)

MEKFNTISVSKCNKETVAMKEKCKFGMTSTIPSGHVWRNTWNPVSCSLATVKMKECLRGKLIYLMGDSTI
 RQWMEYFKASINTLKSVDLHESGKLQHLAVDLDRNINIQQWKYCYPLIGSMTYSVKEMEYLTRAIDRTG
 GEKNTVIVISLQGHFRFPIDVFIRRALNVHKAIQHLLLRSPDTMVIK TENIREMYNDAERFSDFHGYI
 QYLI IKDIFQDLSVSIIDAWDIT IAYGTNNVHPPQHVVGNQINILLNYIC

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_017678

ORF Size: 780 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_017678.3
RefSeq Size:	1715 bp
RefSeq ORF:	783 bp
Locus ID:	54827
UniProt ID:	Q6UWF7
Cytogenetics:	11q23.2
Protein Families:	Secreted Protein, Transmembrane