

Product datasheet for **RG214446**

ZNF358 (NM_018083) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ZNF358
Synonyms:	ZFEND
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >RG214446 representing NM_018083
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGCGGCGCTCAGTCTGGTCAGGAACCCAGGCCACAAGGCCTGAGACCCGTTTATGAAGAGCTCGACT
 CTGACTCCGAGGACCTAGACCCCAATCCTGAAGATCTGGACCCGGTTTCTGAAGACCCAGAGCCTGATCC
 TGAAGACCTCAACTGTCCCGGAAGACGTGGACCCAGCTATGAAGATCTGGAGCCCGTCTCGGAGGAT
 CTGGACCCCGACGCCGAAGCTCCGGGCTCGGAACCCCAAGATCCCGACCCCATGTCTTCGAGTTTCGACC
 TCGATCCAGATGTGATTGGCCCCGTACCCCTGATTCTCGATCCTAACAGCGACACCCCTCAGCCCCGGGGA
 TCCAAAAGTGGACCCCATCTCCTCTGGCCTCACTGCCACCCCCAGGTCTTGGCCACAGCCCCGGGCTG
 CTCCCCGCCCCCGCAGCCCGCCCGGCCCTTCTCCTGCCGGATTGCGGGCGAGCCTTCCGCCGAGCT
 CCGGGCTGAGCCAGCATCGCCGCACGCACAGCGCGAGAAGCCGTACCGCTGCCCGACTGCGGGAAGTC
 CTTAGCCACGGTGCCACCCTGGCTCAGCACCGTGGCATCCACACTGGGGCGCGGCCGTACCACTGCGCG
 GGCTGCGGCAAGGCCTTCGGCTGGCGCTCCACGCTGCTGAAACATCGCAGCAGCCACAGCGGGGAGAAGC
 CGCACCACTGCCGGTGTGTGGCAAGGCCTTCGGGCACGGCTCGCTCCTGGCACAGCACCTGCGCACGCA
 CGGGCGGCCCGGCCCAAGTGCCCGGTGTGCGCCAAGGGCTTCGGCCAGGGCTCTGCGCTGCTCAAA
 CACCTGCGCACGCACACGGGCGAGCGGCCCTACCCGTGTCCGACGTGCGCAAGGCCTTCGGGCAGAGCT
 CGGCGCTGCTGCAGCACAGCGCACACACAGGCGGAGCGGCCCTACCGCTGCCCCCACTGCGGCAAGC
 CTTGCGGCAGAGCTCCTCACTTGAACACACCTGCGCATCCACACGGGCGAGCGGCCCTACGCTGCCCG
 CACTGCTCAAGGCCTTCGGCCAGAGCTCAGCGCTGCTCCAGCACCTGCACGTGCATTCGGGCGAGCGTC
 CCTATCGTGTGAGCTCTGCGGAAGGCCTTCGGCCAGGCCTCCAGCCTACCAAGCACAAACGGGTGCA
 TGAGGGTGCAGCCGCTGCTGCAGCTGCCGCGGCCGCTGCAGCTGCAGCAGCGGCCCGCCGCTGGGCTC
 GGGCCTGGCCTAAGCCCTGCATCCATGATGAGGCCGGGGCAGGTCTCCCTCCTGGGTCTGATGCTGTTT
 CTGTGCTCGGCTCTGGCTTGGGCTCAGCCCTGGCACCAGCTCTGGCCGAACCCCTGACCCCTGGCTCTGG
 GCCGGGCACTCTGCCGATCCAGCTCCAAACCCCTCCCGGCTCCAGATCCACCCAGCCCTACTCCT
 GTGGAATCTTCTGACCCAAAGGCTGGGCACGACGCTGGTCCCGACCTTGTGCCAGCCAGACCTTGATC
 CTGTGCCCAGCCAGACCTGATCCTGTGCCCAGCCCTGATCCCAACCTGTGTCTGCCCTGACCCCTG
 TTCTCCCACTCGTGGCACTGTGAGCCAGCCCTCCCTACCGCGAGAGTCCAGAGTGGGTACAGGAGCAA
 GGGCACTGCTGGGCTGATGGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG214446 representing NM_018083
 Red=Cloning site Green=Tags(s)

MRRSVLVRNPGHKGLRPVYEELDSSELDPNPEDLDPVSEDPEPDPELNTVPEDVDPSEYEDLEPVSED
 LDPDAEAPGSEPQDPDPMSFFDLDPDVIIGPVPLILDPNSDTLSPGDPKVDPISSGLTATPQVLATSPAV
 LPAPASPPRPFSCPDGCAFRSSGLSQHRRTHSGEKPYPYRCPDCGKSFSGATLAQHRGIHTGARPYQCA
 ACGKAFGWRSTLLKRRSSHSGEKPHHCPVCGKAFGHGSLLAQHLRTHGGPRPHKCPVCAKFGQGSALLK
 HLRTHTGERPYPCPQCGKAFGQSSALLQHQRTHTAERPYPYRCPHCGKAFGQSSNLQHHLRIHTGERPYACP
 HCSKAFGQSSALLQHLHVHSGERPYPYRQQLCGKAFGQASSLTKHKRVHEGAAAAAAAAAAAAAAAAAGLGL
 GPGLSPASMMRPQVSLGPDVSVLGSGLGLSPGTSSGRNPDGSGPGTLDPSSKPLPGSRSTPSPTP
 VESSDPKAGHDAGPDLVPSPLDPVSPDPDPVSPDPNPVSCPDPCSPTRGTVPALPTGESPEWVQE
 GALLGPDG

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_018083

ORF Size: 1704 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_018083.5](#)

RefSeq Size: 2003 bp

RefSeq ORF:	1707 bp
Locus ID:	140467
UniProt ID:	Q9NW07
Cytogenetics:	19p13.2
Domains:	zf-C2H2
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:


Circular map for RG214446